

<p align="center">2 DETERMINATION OF SPECIES ORIGIN</p>	<p align="center">Page 1 of 9</p>
<p align="center">PRESUMPTIVE AND CONFIRMATORY TESTS FOR BIOLOGICAL SUBSTANCES – FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION II</p>	<p>Amendment Designator: 2A</p>
	<p>Effective Date: 17-June-2005</p>
<p>2 DETERMINATION OF SPECIES ORIGIN</p> <p>2.1 Quality Control of Antiserums</p> <p>2.1.1 Before using any new lot number of precipitating antiserum for testing casework samples, the specificity must be tested and appropriately documented in the laboratory's quality control (QC) records. Routine quality control testing will be performed at the time the antiserum is reconstituted, or if the antiserum is received in liquid form, within one week of receipt.</p> <p>2.1.2 Anti-human serum, not anti-human hemoglobin, will be used with the procedure in this section for determining whether a sample is of human origin. Anti-human serum, as well as all animal antisera in the "Routine Species Collection" specified below, must be tested against all available species (normal or whole sera and known bloods) in the "Routine Species Collection".</p> <p>2.1.3 A positive control, a host control (typically normal rabbit serum or normal goat serum), and a negative control (distilled water) must be included in the specificity testing. The host control (representing the animal in which the antiserum was prepared) is used to demonstrate that the antiserum is not reacting with any proteins in the animal in which it was made.</p> <p>2.1.4 The quality control documentation will include:</p> <p>2.1.4.1 A diagrammatic representation of the placement of samples in the wells of the rosette on the Ouchterlony plate with the subsequent observations (i.e., precipitin lines).</p> <p>2.1.4.2 Date of the testing.</p> <p>2.1.4.3 Initials of the person conducting the testing.</p> <p>2.1.4.4 Lot number, date of receipt, and manufacturer of the antiserum being tested.</p> <p>2.1.4.5 Lot number, date of receipt, and manufacturer of the normal sera being used for the testing.</p> <p>2.1.4.6 Results of the testing.</p> <p>2.1.5 Once the appropriate testing has been performed on a particular lot number of antiserum, it need not be repeated for each case. If another vial of the same lot number is received on a different date, the QC testing described above must be repeated.</p> <p>2.2 Quality Control of Normal (Whole) Sera</p> <p>2.2.1 Before using any new lot number of normal serum for testing casework samples, the specificity must be tested and appropriately documented in the laboratory's quality control (QC) records. Routine quality control testing will be performed at the time the normal serum is reconstituted, or if the normal serum is received in liquid form, within one week of receipt.</p>	

2 DETERMINATION OF SPECIES ORIGIN	Page 2 of 9
PRESUMPTIVE AND CONFIRMATORY TESTS FOR BIOLOGICAL SUBSTANCES – FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION II	Amendment Designator: 2A
	Effective Date: 17-June-2005
<p>2.2.2 All normal serums in the “Routine Species Collection” must be tested against all available antisera in the “Routine Species Collection”.</p> <p>2.2.3 A positive control must be included in the specificity testing. Although normal human serum may be purchased, <u>a straw colored dilution</u> of known human blood may be used instead. Similarly, the use of known blood from other species may replace the purchase and use of normal serums from the species.</p> <p>2.2.4 It is not necessary to conduct quality control testing on the known bloods. Label known bloods with the species name and date of preparation/initials of person preparing the sample. Store known bloods in the freezer.</p> <p>2.2.5 The quality control documentation will include:</p> <p>2.2.5.1 A diagrammatic representation of the placement of samples in the wells of the rosette on the Ouchterlony plate with the subsequent observations (i.e., precipitin lines).</p> <p>2.2.5.2 Date of the testing.</p> <p>2.2.5.3 Initials of the person conducting the testing.</p> <p>2.2.5.4 Lot number, date of receipt, and manufacturer of the normal serum being tested.</p> <p>2.2.5.5 Lot number, date of receipt, and manufacturer of the antisera being used for the testing.</p> <p>2.2.5.6 Results of the testing.</p> <p>2.2.6 Once the appropriate testing has been performed on a particular lot number of normal serum, it need not be repeated for each case. If another vial of the same lot number is received on a different date, the QC testing described above must be repeated.</p> <p>2.3 Routine Species Collection</p> <p>2.3.1 The following will be maintained in each laboratory as the “Routine Species Collection” and will undergo QC testing for specificity as outlined above:</p> <p>2.3.1.1 Bovine antiserum and normal bovine serum or known blood</p> <p>2.3.1.2 Swine antiserum and normal swine serum or known blood</p> <p>2.3.1.3 Cat antiserum and normal cat serum or known blood</p> <p>2.3.1.4 Dog antiserum and normal dog serum or known blood</p> <p>2.3.1.5 Rabbit antiserum and normal rabbit serum or known blood</p> <p>2.3.1.6 Sheep antiserum and normal sheep serum or known blood</p> <p>2.3.1.7 Deer antiserum and normal deer serum or known blood</p>	

2 DETERMINATION OF SPECIES ORIGIN	Page 3 of 9
PRESUMPTIVE AND CONFIRMATORY TESTS FOR BIOLOGICAL SUBSTANCES – FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION II	Amendment Designator: 2A
	Effective Date: 17-June-2005
<p>2.3.1.8 Human antiserum and normal human serum or known blood</p> <p>2.3.1.9 Normal goat serum or known blood (goat antiserum is unavailable)</p> <p>NOTE: It is recommended that the antisera and normal sera listed above be reconstituted and QC tested at the time of receipt to ensure ready availability.</p> <p>2.3.2 Other antisera/normal sera, such as bear, rodent, fowl, horse, etc. should also be maintained for use in special cases and must also undergo QC testing as specified above. However, since these are used only in special cases, it is recommended that they not be reconstituted, aliquoted, and QC tested until it is determined that there is a specific need to do so.</p> <p>2.4 Storage of Antiserum/Normal Serum</p> <p>2.4.1 Small aliquots of the antiserum/normal serum will be prepared for routine use and frozen within one week of reconstitution (when antiserum is lyophilized) or upon receipt (when antiserum is liquid). All frozen aliquots have an indefinite expiration date.</p> <p>2.4.2 A thawed aliquot may be stored refrigerated for up to 1 month. If the aliquot is to be maintained in this manner, the expiration date must be clearly marked on the vial. Otherwise, the aliquot will be immediately discarded following its use.</p> <p>2.5 Labeling of Antiserum/Normal Serum</p> <p>2.5.1 Labels on each aliquot will include:</p> <p>2.5.1.1 The manufacturer</p> <p>2.5.1.2 Type of antiserum or normal serum</p> <p>2.5.1.3 Lot number</p> <p>2.5.1.4 Date reconstituted/date frozen</p> <p>2.5.1.5 Initials of the person preparing the aliquot</p> <p>2.6 OUCHTERLONY (DOUBLE DIFFUSION) TEST (Reference 5, pp. 221-241, Appendix A)</p> <p>2.6.1 Equipment</p> <p>2.6.1.1 Punch</p> <p>2.6.1.2 Aspirator</p> <p>2.6.1.3 100 ml and 500 ml graduated cylinders</p> <p>2.6.1.4 Balance</p>	

2 DETERMINATION OF SPECIES ORIGIN	Page 4 of 9
PRESUMPTIVE AND CONFIRMATORY TESTS FOR BIOLOGICAL SUBSTANCES – FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION II	Amendment Designator: 2A
	Effective Date: 17-June-2005
<div>2.6.1.5 Spatula</div> <div>2.6.1.6 Scissors</div> <div>2.6.1.7 Tweezers</div> <div>2.6.1.8 Hot plat or oven (37° C)</div> <div>2.6.1.9 Incubator (optional)</div> <div>2.6.1.10 Magnetic stir plate</div> <div>2.6.1.11 Refrigerator (optional)</div> <div>2.6.2 Materials</div> <div>2.6.2.1 Petri dishes, slides, or comparable containers</div> <div>2.6.2.2 Test tubes</div> <div>2.6.2.3 Weigh boat or weigh paper</div> <div>2.6.2.4 Moisture chamber</div> <div>2.6.2.5 Disposable pipets</div> <div>2.6.2.6 Capillary tubes</div> <div>2.6.3 Reagents</div> <div>2.6.3.1 Normal saline (0.9%)</div> <div>2.6.3.2 Agarose gel (1%)</div> <div>2.6.3.3 Distilled water</div> <div>2.6.3.4 Antiserum</div> <div>2.6.3.5 Normal serum or known blood (positive control)</div> <div>2.6.4 Agarose Gel Preparation</div> <div>2.6.4.1 Normal saline (0.9% NaCl):</div> <div> <ul style="list-style-type: none"> • 9 g Sodium chloride • 1000 ml Distilled water </div>	

2 DETERMINATION OF SPECIES ORIGIN	Page 5 of 9
PRESUMPTIVE AND CONFIRMATORY TESTS FOR BIOLOGICAL SUBSTANCES – FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION II	Amendment Designator: 2A
	Effective Date: 17-June-2005
<div data-bbox="462 317 867 348" data-label="List-Group"> <ul style="list-style-type: none"> • Mix thoroughly until dissolved. </div> <div data-bbox="347 384 683 415" data-label="Section-Header"> <p>2.6.4.2 Agarose gel (1%):</p> </div> <div data-bbox="462 453 1136 556" data-label="List-Group"> <ul style="list-style-type: none"> • 1 g Type I agarose • 100 ml Normal saline (0.9% NaCl) • Heat until agarose is dissolved. Allow to cool slightly. </div> <div data-bbox="347 590 1541 653" data-label="Text"> <p>2.6.4.3 Pour the agarose into a petri dish, onto a slide, or into a comparable container to a thickness of 2-3 mm and allow to cool.</p> </div> <div data-bbox="347 688 1482 756" data-label="Text"> <p>2.6.4.4 Cut wells in a rosette pattern (refer to the diagram below) into the gel using a punch or disposable pipet connected to an aspirator.</p> </div> <div data-bbox="724 772 984 1052" data-label="Image"> </div> <div data-bbox="250 1092 594 1125" data-label="Section-Header"> <p>2.6.5 Storage and Labeling</p> </div> <div data-bbox="347 1159 1546 1260" data-label="Text"> <p>2.6.5.1 When a batch of plates is prepared, the plates should be numbered consecutively and placed in a moisture chamber in the refrigerator. Label the moisture chamber with the lot number of the batch (date of preparation/initials of person preparing the plates).</p> </div> <div data-bbox="347 1291 1330 1327" data-label="Text"> <p>2.6.5.2 There is no expiration date (see 2.6.6 Minimum Standards and Controls).</p> </div> <div data-bbox="250 1358 742 1392" data-label="Section-Header"> <p>2.6.6 Minimum Standards and Controls</p> </div> <div data-bbox="347 1428 1524 1528" data-label="Text"> <p>2.6.6.1 A positive control (known sample against which the antiserum is directed) and a substrate control (or if not available, distilled water) must be tested on each plate, unless the stain is on a cotton swab. It is not necessary to test submitted control swabs.</p> </div> <div data-bbox="250 1560 1031 1593" data-label="Section-Header"> <p>2.6.7 OUCHTERLONY DOUBLE DIFFUSION PROCEDURE</p> </div> <div data-bbox="347 1627 1546 1761" data-label="Text"> <p>2.6.7.1 To prepare an extract of the stain, place a small cutting of the stain in distilled water until a <u>straw color</u> is obtained. A small piece of stained material, which is moistened with distilled water, can be used in lieu of an extract. Treat the substrate control in the same manner as the stain.</p> </div> <div data-bbox="347 1795 1451 1862" data-label="Text"> <p>2.6.7.2 Add antiserum in the center well of the Ouchterlony plate with a disposable pipet or capillary tube.</p> </div>	

2 DETERMINATION OF SPECIES ORIGIN	Page 6 of 9
PRESUMPTIVE AND CONFIRMATORY TESTS FOR BIOLOGICAL SUBSTANCES – FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION II	Amendment Designator: 2A
	Effective Date: 17-June-2005
<p>2.6.7.3 Add appropriate extracts/pieces of stained material, the positive control, and negative control(s) to the surrounding wells. Do not overfill the wells. Avoid getting bubbles in the wells. Document the placement of samples in the wells of the rosette on the Ouchterlony work sheet found in Appendix B.</p> <p>2.6.7.4 Record the lot numbers of antiserums and normal serums used for the testing procedure.</p> <p>NOTE: Alternatively, known normal serum/known blood extract may be placed in the center well with appropriate antiserums in the surrounding wells or the stain extract/piece of stained material may be placed in the center well with appropriate antiserums in the surrounding wells.</p> <p>2.6.7.5 Incubate the plate in a moisture chamber at 37° C for 3-4 hours. Alternatively, it may be left overnight at room temperature or 4° C.</p> <p>2.6.7.6 Record observations (precipitin lines) on the diagram, and interpret and record the results.</p> <p>2.6.7.7 All controls must give the expected results before a conclusion can be reached on an unknown sample, i.e., white precipitin lines must be observed between the antiserum and positive control (known serum) and no precipitin lines should be observed between the antiserum and the negative control.</p> <p>2.6.7.8 Interpretation</p> <p>2.6.7.8.1 Positive Result = White precipitin lines between the antiserum well and the sample well</p> <p>2.6.7.8.2 Negative Result = No precipitin lines between the antiserum well and the sample well</p> <p>2.6.7.8.3 Inconclusive Result = Questionable precipitin lines between the antiserum well and the sample well, OR precipitin lines between the antiserum well and the negative control well(s), OR no precipitin line between the antiserum well and the positive control well. If sufficient sample remains, an inconclusive result should be repeated.</p> <p>NOTE: The prozone phenomenon can result in a soluble antigen-antibody complex due to too many antibodies present to form a complete lattice (Reference 6, Appendix A). Because of this phenomenon, weak precipitin lines may be observed initially, but disappear upon staining with Coomassie Brilliant Blue R250. This is considered an inconclusive result. Therefore, it is recommended that the testing results be recorded PRIOR TO staining as well as after staining.</p> <p>2.6.7.9 Staining the plate with Coomassie Brilliant Blue R250 may be necessary to visualize weak reactions. Refer to 2.7 for the Coomassie Brilliant Blue R250 staining procedure.</p> <p>2.6.7.10 Reporting Results</p> <p>2.6.7.10.1 Report positive test results as “(species tested according to the label on the <u>antiserum</u>) protein was detected...”</p>	

2 DETERMINATION OF SPECIES ORIGIN	Page 7 of 9
PRESUMPTIVE AND CONFIRMATORY TESTS FOR BIOLOGICAL SUBSTANCES – FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION II	Amendment Designator: 2A
	Effective Date: 17-June-2005
<div data-bbox="467 346 1528 520"> <p>2.6.7.10.2 Report negative test results as “no (<u>species tested according to the label on the antiserum</u>) protein was detected...”</p> <p>2.6.7.10.3 Report inconclusive test results as “the test for (<u>species tested according to the label on the bottle</u>) protein was inconclusive...”</p> </div> <div data-bbox="196 548 1393 583"> <p>2.7 COOMASSIE BRILLIANT BLUE R250 STAINING PROCEDURE (Reference 6, Appendix A)</p> </div> <div data-bbox="253 615 602 651"> <p>2.7.1 Safety Considerations</p> </div> <div data-bbox="347 680 1471 884"> <p>2.7.1.1 Coomassie Brilliant Blue R250 - Caution! Avoid contact and inhalation! Emits toxic fumes under fire conditions! Container explosion may occur under fire conditions!</p> <p>2.7.1.2 Methanol - Caution! Irritant! Dangerous when exposed to heat or flame!</p> <p>2.7.1.3 Glacial acetic acid - Caution! Corrosive! Flammable!</p> </div> <div data-bbox="253 915 474 951"> <p>2.7.2 Equipment</p> </div> <div data-bbox="347 982 1002 1423"> <p>2.7.2.1 Weight</p> <p>2.7.2.2 Oven or incubator (40-60°C)</p> <p>2.7.2.3 Rotator (optional)</p> <p>2.7.2.4 10 ml, 50 ml, and 500 ml graduated cylinders</p> <p>2.7.2.5 Balance</p> <p>2.7.2.6 Spatula</p> <p>2.7.2.7 Trays for staining and destaining</p> </div> <div data-bbox="253 1453 454 1488"> <p>2.7.3 Materials</p> </div> <div data-bbox="347 1520 1019 1759"> <p>2.7.3.1 Gel bond, glass plate, or other support medium</p> <p>2.7.3.2 Weigh boats or weigh paper</p> <p>2.7.3.3 Whatman #1 filter paper</p> <p>2.7.3.4 Paper towels</p> </div> <div data-bbox="253 1789 451 1824"> <p>2.7.4 Reagents</p> </div> <div data-bbox="347 1856 673 1892"> <p>2.7.4.1 Staining Solution</p> </div>	

2 DETERMINATION OF SPECIES ORIGIN	Page 8 of 9
PRESUMPTIVE AND CONFIRMATORY TESTS FOR BIOLOGICAL SUBSTANCES – FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION II	Amendment Designator: 2A
	Effective Date: 17-June-2005
<div data-bbox="349 310 703 346">2.7.4.2 Destaining Solution</div> <div data-bbox="349 380 989 415">2.7.4.3 Normal saline (0.9% NaCl) – refer to 2.7.8.1</div> <div data-bbox="349 449 643 485">2.7.4.4 Distilled water</div> <div data-bbox="256 518 915 554">2.7.5 Preparation of Staining and Destaining Solutions</div> <div data-bbox="349 585 673 621">2.7.5.1 Staining Solution</div> <div data-bbox="462 655 1122 831"> <ul style="list-style-type: none"> • 0.1 g Coomassie Brilliant Blue R250 • 45.0 ml Methanol • 10.0 ml Glacial acetic acid • 45.0 ml Distilled water • Mix the above ingredients until thoroughly dissolved. </div> <div data-bbox="349 863 703 898">2.7.5.2 Destaining Solution</div> <div data-bbox="462 932 1122 1073"> <ul style="list-style-type: none"> • 45.0 ml Methanol • 10.0 ml Glacial acetic acid • 45.0 ml Distilled water • Mix the above ingredients until thoroughly dissolved. </div> <div data-bbox="256 1104 433 1140">2.7.6 Storage</div> <div data-bbox="349 1173 1287 1209">2.7.6.1 The Staining and Destaining Solutions are stable at room temperature.</div> <div data-bbox="256 1241 448 1276">2.7.7 Labeling</div> <div data-bbox="349 1310 1528 1377">2.7.7.1 Label as Staining or Destaining Solution with the lot number (date of preparation followed by the initials of the person preparing the solution).</div> <div data-bbox="349 1409 789 1444">2.7.7.2 There is no expiration date.</div> <div data-bbox="256 1476 1154 1512">2.7.8 COOMASSIE BRILLIANT BLUE R250 STAINING PROCEDURE</div> <div data-bbox="349 1545 1468 1612">2.7.8.1 Wash the plate overnight in normal saline solution (0.9% NaCl - 9 g NaCl in 1000 ml distilled water) to remove unprecipitated proteins.</div> <div data-bbox="462 1644 1536 1711">2.7.8.1.1 If pieces of stained material were used in lieu of extracts, remove these prior to washing the plate.</div> <div data-bbox="462 1743 1471 1848">2.7.8.1.2 The gel may detach from the plate during the washing process. Mark the orientation of the gel to ensure that it can be re-oriented properly after the washing has been completed.</div>	

2 DETERMINATION OF SPECIES ORIGIN	Page 9 of 9
PRESUMPTIVE AND CONFIRMATORY TESTS FOR BIOLOGICAL SUBSTANCES – FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION II	Amendment Designator: 2A
	Effective Date: 17-June-2005
<div data-bbox="467 310 1555 1791"> <div>2.7.8.1.3</div> <p>Alternatively, the gel may be removed from the plate before beginning the washing process. If this is done, mark the gel to ensure that it can be re-oriented properly after the washing has been completed.</p> <div>2.7.8.2</div> <p>The next day wash the gel for approximately fifteen minutes in distilled water. Rinse the gel and repeat the wash. Ensuring proper orientation, place the gel (face up) on the hydrophilic side of a piece of gel bond or on some other support medium such as a glass plate.</p> <div>2.7.8.3</div> <p>Cover the gel with a piece of Whatman #1 filter paper moistened with distilled water. Add a layer of paper towels on top of the filter paper and press with a weight for approximately 30 minutes. Remove paper towels and filter paper and dry the gel in a 40-60° C oven.</p> <div>2.7.8.4</div> <p>Place the gel in the staining solution and allow it to soak for less than 1 minute up to 10 minutes. This may be done on a rotator. Intermittently check staining progress to prevent over staining.</p> <div>2.7.8.5</div> <p>Place the gel in the destaining solution until the background is clear or until no more dye leaches from the gel. This may be done on a rotator. Change the destaining solution and destain further if desired.</p> <div>2.7.8.6</div> <p>Record observations (precipitin lines) on the diagram, and interpret and record the results. All controls must give the expected results before a conclusion can be reached on an unknown sample, i.e., blue precipitin lines must be observed between the antiserum and positive control (known serum) and no precipitin lines should be observed between the antiserum and the negative control.</p> <div>2.7.8.7</div> <p>Interpretation</p> <div>2.7.8.7.1</div> <p>Positive Result = Blue precipitin lines between the antiserum well and the sample well.</p> <div>2.7.8.7.2</div> <p>Negative Result = No precipitin lines between the antiserum well and the sample well.</p> <div>2.7.8.7.3</div> <p>Inconclusive Result = Questionable precipitin lines between the antiserum well and the sample well, OR precipitin lines between the antiserum well and the negative control well(s), OR no precipitin line observed between the antiserum well and the positive control well. If sufficient sample remains, an inconclusive result should be repeated.</p> <div>2.7.8.8</div> <p>Reporting Results</p> <div>2.7.8.8.1</div> <p>Refer to 2.6.7.10.</p> </div> <div data-bbox="1393 1892 1468 1923">◆END</div>	