

8 FLUORESCENT DETECTION OF THE ELECTROPHORESIS GEL	Page 1 of 1
<b>FLUORESCENT DETECTION PCR-BASED STR DNA PROTOCOL:POWERPLEX®16 BIO SYSTEM - FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION III</b>	Issue No. 3
	Effective Date: 6-March-2006
<p><b>8 FLUORESCENT DETECTION OF THE ELECTROPHORESIS GEL</b></p> <p><b>NOTE:</b> Review the FMBIO manual before proceeding. Refer to this manual for further details on the following procedures.</p> <p>8.1 EQUIPMENT</p> <p>8.1.1 FMBIO II or III Plus Fluorescent Imaging Analysis System (Hitachi/MiraiBio)</p> <p>8.1.2 Power Mac or Mac Clone Computer</p> <p>8.1.3 Laser Printer</p> <p>8.2 MATERIALS</p> <p>8.2.1 Kimwipes</p> <p>8.3 REAGENTS</p> <p>8.3.1 Type 1 Water</p> <p>8.3.2 70 % Ethanol</p> <p>8.4 FLUORESCENT DETECTION FOR POWERPLEX® SYSTEM GELS</p> <p>8.4.1 The FMBIO and Macintosh/PC computer should be turned on at least 30 minutes before scanning a typing gel.</p> <p>8.4.2 Clean the glass plate setup by placing it under running water and removing all the acrylamide and urea from the outside.</p> <p>8.4.3 Rinse the glass plate setup with Type I Water and dry with a Kimwipe. Remove any smuggles using 70% Ethanol and a Kimwipe.</p> <p>8.4.4 Refer to Appendix F (FMBIO II) or Appendix G (FMBIO III Plus) for procedure for operating the instrument.</p> <p><b>NOTE:</b> If the polyacrylamide gel will be reused after it has been scanned, refer to Appendix N, Reuse of A Polyacrylamide Electrophoresis Gels, for the procedure.</p> <p style="text-align: right;"><b>◆END</b></p>	