

APPENDIX A: CONVERSION FACTORS AND OTHER USEFUL INFORMATION	Page 1 of 2
FLUORESCENT DETECTION PCR-BASED STR DNA PROTOCOL:POWERPLEX®16 BIO SYSTEM - FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION III	Issue No. 3
	Effective Date: 6-March-2006
<p>APPENDIX A: CONVERSION FACTORS AND OTHER USEFUL INFORMATION</p> <p>I. CONVERSION FACTORS</p> <p>One microgram = 1 µg = 1 x 10⁻⁶ g = 1000 ng</p> <p>One nanogram = 1 ng = 1 x 10⁻⁹ g = 1000 pg</p> <p>One picogram = 1 pg = 1 x 10⁻¹² g = 1 x 10⁻³ ng</p> <p>One liter = 1 L = 1000 mL</p> <p>One milliliter = 1 mL = 1 x 10⁻³ (0.001) L = 1000 µL</p> <p>One microliter = 1 µL = 1 x 10⁻³ mL = 0.001 mL</p> <p>II. CONCENTRATIONS</p> <p>A. Concentration can be expressed several ways:</p> <ol style="list-style-type: none"> Weight percent = (mass A/total mass of solution) x 100 or simplified: Wt %_A = (g_A/100 mL of solution) x 100 Volume percent = (volume A/total volume of solution) x 100 or simplified: Volume %_A = (volume_A/100 mL of solution) x 100 Molarity (<i>M</i>): <i>M</i> = no. moles solute A/no. liters solution = molecular weight of solute A in 1000 mL solution where 1 mole of A = 1 gram formula weight of A Normality (<i>N</i>) = no. MW/no. liters solution where in acid-base reactions: MW acid = weight of acid which reacts with 1 mole of OH⁻ MW base = weight of base which reacts with 1 mole of H⁺ The normality of a given reagent depends on the reaction in which it participates. (Example: 1 L of 1M H₃PO₄ which can have <i>N</i> = 1, 2 or 3 depending upon the reaction in which it is involved.) Because of ambiguities, the concept of normality is to be used carefully. 	

APPENDIX A: CONVERSION FACTORS AND OTHER USEFUL INFORMATION	Page 2 of 2
FLUORESCENT DETECTION PCR-BASED STR DNA PROTOCOL:POWERPLEX®16 BIO SYSTEM - FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION III	Issue No. 3
	Effective Date: 6-March-2006
<p>B. Examples:</p> <p>1. Prepare 50% solution of polyethylene glycol (PEG):</p> <p>Place 50 g of PEG in a flask and bring to a volume of 100 mL with H₂O.</p> <p>Wt % PEG = (50 g/100 mL) x 100 = 50%</p> <p>2. Prepare 100 mL of 2<i>M</i> HCl from a stock solution of 12<i>M</i> HCl. (The question is how much stock solution of 12<i>M</i> HCl is needed?)</p> <p>Conc of Stock Soln = C_s = 12<i>M</i> Conc of Final Soln = C_f = 2<i>M</i></p> <p>Volume of Stock Soln = V_s = ? Volume of Final Soln = V_f = 0.1 L</p> <p>$C_s \cdot V_s = C_f \cdot V_f$</p> <p>$V_s = (C_f \cdot V_f)/C_s = 2(0.1)/12 = 0.2/12 = 0.0167 \text{ L}$</p> <p>$V_s = 0.0167 \text{ L} = 16.7 \text{ mL}$</p> <p>Answer: Take 16.7 mL of 12<i>M</i> HCl and dilute to final volume of 100 mL to give 2<i>M</i> HCl.</p> <p style="text-align: right;">◆END</p>	