

<p align="center">1 ISOLATING DNA USING THE DNA IQ™ ISOLATION SYSTEM</p>	<p align="center">Page 1 of 8</p>
<p align="center">BIOMEK®2000 AUTOMATION WORKSTATION PROCEDURES MANUAL- FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION IV</p>	<p align="center">Issue No. : 4</p>
	<p align="center">Effective Date: 8-August-2005</p>
<p>1 ISOLATING DNA USING THE DNA IQ™ ISOLATION SYSTEM</p> <p>When the DNA examiner has completed his/her analysis, all remaining evidence samples, including extracted DNA which has been dried down using a DNA concentrator/evaporator (Refer to Appendix K, Procedure For Drying Down And Resolubilizing Extracted DNA, of Section III of the Forensic Biology Section Manual) or other suitable methods, will be returned to the submitting agency. However, amplified DNA samples WILL NOT be returned to the submitting agency, but instead will be discarded in the designated area to prevent possible transfer of amplified DNA to the remaining evidence samples. If the biological sample deposited on the evidence is consumed during the analysis then the extracted DNA sample and the cutting (in separate tubes) will be returned to the submitting agency with the evidence. Proper case file documentation and chain of custody documentation for these samples must be maintained.</p> <p>DNA for PCR amplification and analysis may be extracted from bloodstains, sperm cells, buccal cells, hair, tissue, bone, and other samples. Slightly different extraction procedures are required for each type of specimen and therefore are outlined in this chapter.</p> <p>It is important to handle all samples aseptically to prevent contamination by extraneous DNA. It is also important to prepare evidence samples at a separate time and/or space from reference samples to prevent possible cross-contamination.</p> <p>NOTE: For tracking purposes the samples will be listed on the worksheet in the order in which they were processed/handled and will be processed in accordance with the procedures outlined in the <u>Commonwealth of Virginia Department of Forensic Science Forensic Biology Section Procedure Manual, Section VI - Quality Assurance Program DNA Typing of Biological Materials</u>, Chapter 6, Dedicated PCR Facilities.</p> <p>Special Precautions:</p> <ul style="list-style-type: none"> • The DNA extraction and PCR setup of evidence samples will be performed at a separate time or space from the DNA extraction and PCR setup of reference samples. This helps to prevent potential cross-contamination between evidence samples and reference samples. • DNA extraction from samples containing high levels of DNA (for example, tissue) will be performed separately from samples expected to contain low levels of DNA (single hairs, small bloodstains, etc.) to minimize the potential for sample-to-sample contamination. • Disposable gloves will be used at all times. Gloves will be changed frequently to avoid sample-to-sample contamination with DNA and whenever moving between work areas. Gloves will be changed if suspected direct contamination has occurred from the sample DNA. • Scissors will be thoroughly cleaned with a 10% solution of bleach or a solution that will remove/degrade the DNA after cutting each item/stain. Subsequently use Isopropyl Alcohol to remove the residue left by the chemicals, using special care to remove all residue left on surfaces. A fresh scalpel blade may also be used to cut each item/stain. • A clean cutting surface will be used for each piece of evidence. • Disposable plugged pipette tips and microcentrifuge tubes will be used. 	

<p align="center">1 ISOLATING DNA USING THE DNA IQ™ ISOLATION SYSTEM</p>	<p align="center">Page 2 of 8</p>
<p align="center">BIOMEK®2000 AUTOMATION WORKSTATION PROCEDURES MANUAL- FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION IV</p>	<p align="center">Issue No. : 4</p>
	<p align="center">Effective Date: 8-August-2005</p>
<ul style="list-style-type: none"> • Pipette tips will be changed between samples. • Reagents will be stored in small quantities to reduce the risk of possible contamination to the stock solution. • To avoid splashing and minimize aerosols, all liquid will be centrifuged to the bottom of the closed tube before the tube is opened. • Reagent blank controls will be included with each set of DNA extractions to check for the presence of contaminating DNA in the reagents. • All work surfaces will be thoroughly cleaned with a 10% bleach solution or a solution that will remove/degrade the DNA. Subsequently use Isopropyl Alcohol to remove the residue left by the chemicals, using special care to remove all residue left on surfaces before setting up the DNA Extraction Work Area. Disposable bench paper will be used to prevent the accumulation of human DNA on permanent work surfaces. • The quantity of samples handled during a single analysis will be limited to a manageable number. This precaution reduces the risk of sample mix-up and the potential for sample-to-sample contamination. • A dedicated lab coat will be worn for pre-amplification sample handling when working in the DNA Extraction Work Area. • A dedicated lab coat will be worn when working with amplified DNA in the PCR Post Amplification Work Area. <p>1.1 TECHNICAL NOTES</p> <p>1.1.1 The DNA IQ™ Isolation System is designed to rapidly purify small quantities of DNA, approximately 100 ng or less, and becomes more efficient with samples containing less than 50 ng of DNA.</p> <p>1.1.2 The DNA IQ™ Isolation System will isolate genomic DNA in general and is not human specific.</p> <p>1.1.3 A random sample will be run with each set of convicted offender and arrestee sample extractions to serve as a verification that the samples are successfully being entered into Combined DNA Index System (CODIS) and the search algorithm is working properly. If a sample must be re-extracted a new random sample must be extracted along with the sample.</p> <p>1.1.4 Routinely the DNA IQ™ extraction method is used to isolate DNA from blood/buccal samples obtained in criminal cases and the QIAGEN® BioRobot™ 9604 is used to isolate DNA from convicted offenders blood/buccal or arrestee buccal samples in accordance with the <u>Commonwealth of Virginia Department of Forensic Science Forensic Biology Section Procedure Manual, Section V - QIAGEN® BioRobot™ 9604 Procedure Manual</u>. However, when necessary convicted offender blood/buccal and arrestee buccal samples may also be extracted using the procedures outlined in this manual.</p>	

1 ISOLATING DNA USING THE DNA IQ™ ISOLATION SYSTEM	Page 3 of 8
BIOMEK®2000 AUTOMATION WORKSTATION PROCEDURES MANUAL- FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION IV	Issue No. : 4
	Effective Date: 8-August-2005
<div data-bbox="248 321 513 352">1.2 EQUIPMENT</div> <div data-bbox="345 390 1053 1157"> <div data-bbox="345 390 787 422">1.2.1 Heat block or incubator, 37°C</div> <div data-bbox="345 459 787 491">1.2.2 Heat block or incubator, 56°C</div> <div data-bbox="345 529 626 560">1.2.3 Microcentrifuge</div> <div data-bbox="345 598 594 630">1.2.4 Vortex mixer</div> <div data-bbox="345 667 532 699">1.2.5 Scissors</div> <div data-bbox="345 737 548 768">1.2.6 Tweezers</div> <div data-bbox="345 806 743 837">1.2.7 Microcentrifuge tube rack</div> <div data-bbox="345 875 992 907">1.2.8 Pipettes - 10 µL, 100 µL, 200 µL, and 1000 µL</div> <div data-bbox="345 945 678 976">1.2.9 Refrigerator/Freezer</div> <div data-bbox="345 1014 1053 1045">1.2.10 Scalpel and blades (for DNA extraction from swabs)</div> <div data-bbox="345 1083 1026 1115">1.2.11 Stereo microscope (for DNA extraction from hair)</div> <div data-bbox="345 1152 984 1184">1.2.12 Beaker, 50 mL (for DNA extraction from hair)</div> </div> <div data-bbox="248 1222 509 1253">1.3 MATERIALS</div> <div data-bbox="345 1291 1229 1724"> <div data-bbox="345 1291 797 1323">1.3.1 Microcentrifuge tubes, 1.5 mL</div> <div data-bbox="345 1360 634 1392">1.3.2 Transfer pipettes</div> <div data-bbox="345 1430 1229 1461">1.3.3 Sterile ART tips for pipettes - 10 µL, 100 µL, 200 µL, and 1000 µL</div> <div data-bbox="345 1499 735 1530">1.3.4 Microcentrifuge tube lids</div> <div data-bbox="345 1568 558 1600">1.3.5 Kimwipes</div> <div data-bbox="345 1638 768 1669">1.3.6 Spin-Ease basket (Optional)</div> <div data-bbox="345 1707 521 1738">1.3.7 Gloves</div> </div> <div data-bbox="248 1761 495 1793">1.4 REAGENTS</div> <div data-bbox="345 1831 909 1927"> <div data-bbox="345 1831 909 1862">1.4.1 Proteinase K - 20 mg/mL (Keep on ice.)</div> <div data-bbox="345 1900 495 1932">1.4.2 TNE</div> </div>	

1 ISOLATING DNA USING THE DNA IQ™ ISOLATION SYSTEM	Page 4 of 8
BIOMEK®2000 AUTOMATION WORKSTATION PROCEDURES MANUAL- FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION IV	Issue No. : 4
	Effective Date: 8-August-2005
<div data-bbox="342 321 821 966"> <ul style="list-style-type: none"> 1.4.3 20% Sarkosyl 1.4.4 0.39M Dithiothreitol (DTT) 1.4.5 Sterile Type I Water 1.4.6 PCR Digestion Buffer 1.4.7 95% Ethanol 1.4.8 DNA IQ™ System Reaction Kit 1.4.9 DNA IQ™ Lysis Buffer 1.4.10 DNA IQ™ Wash Buffer 1.4.11 DNA IQ™ Elution Buffer 1.4.12 DNA IQ™ Resin </div> <div data-bbox="248 1035 1356 1098"> <p>1.5 DNA IQ™ EXTRACTION METHOD FOR BUCCAL CELL TYPE SAMPLES AND BLOODSTAINS</p> </div> <div data-bbox="342 1136 1549 1304"> <p>The portion size of the swab or bloodstain removed for DNA extraction should be judged based on a number of criteria, such as whether the stain appears dilute (for bloodstains) or if the sample may be heavily soiled or possibly degraded. Examine the bloodstain or buccal swab and remove a reasonable portion for DNA extraction. For example, only a 3 mm² section may be all that is necessary to remove from a heavily bloodstained item.</p> </div> <div data-bbox="342 1337 1084 1371"> <p>NOTE: Process a reagent blank along with each set of samples.</p> </div> <div data-bbox="342 1402 1549 1570"> <p>1.5.1 Cut an approximate 3 mm² blood stain or other biological stain and place into a labeled 1.5 mL microcentrifuge tube with a depression in the lid. If a stain is smeared over a large area of fabric, more than one microcentrifuge tube may be needed to extract the stain. If the sample is a buccal swab, remove a small portion of the swab and place into a labeled microcentrifuge tube.</p> </div> <div data-bbox="342 1606 1549 1675"> <p>NOTE: Prior to the addition of the DNA IQ™ Lysis buffer, DTT at the concentration described in the reagent preparation section (Appendix B) MUST be added.</p> </div> <div data-bbox="342 1709 1549 1778"> <p>1.5.2 Add 250 uL of DNA IQ™ Lysis buffer. If the biological sample is dispersed over a large area of the substrate, add up to 325 uL of DNA IQ™ Lysis buffer to the microcentrifuge tube.</p> </div> <div data-bbox="342 1812 1422 1845"> <p>1.5.3 Vortex vigorously for 20-30 seconds, then pulse spin to force the cutting into liquid.</p> </div>	

1 ISOLATING DNA USING THE DNA IQ™ ISOLATION SYSTEM	Page 5 of 8
BIOMEK® 2000 AUTOMATION WORKSTATION PROCEDURES MANUAL- FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION IV	Issue No. : 4
	Effective Date: 8-August-2005
<div data-bbox="342 289 1549 386"> <p>1.5.4 Place the microcentrifuge tube into a 56°C heat block for a minimum of 30 minutes. If the sample is deposited on FTA paper, place the microcentrifuge tube into a 95° C heat block for a minimum of 30 minutes.</p> </div> <div data-bbox="342 424 1549 621"> <p>1.5.5 Vortex vigorously for 20-30 seconds, then pulse spin the tube, punch 2-3 holes in the lid of the tube, remove the cutting from the liquid and place in the lid. Spin the tube for 5 minutes in a microcentrifuge at a minimum of 10,000 rpm to remove the excess liquid from the cutting. Alternatively, a Spin-Ease basket may be used instead of placing the cutting in the lid of the tube. If the entire biological sample deposited on the evidence was consumed the cutting will be returned with the evidence.</p> </div> <div data-bbox="342 659 1549 831"> <p>NOTE: Samples can be left at room temperature in DNA IQ™ Lysis buffer for up to 24 hours after heating and centrifugation before proceeding to the BioMek® 2000 Automation Workstation. If the samples are stored in a refrigerator before proceeding to Section 2, Preparation of the BioMek® 2000 Automation Workstation, place the samples in a 56°C heat block for 5 minutes to resolubilize the samples.</p> </div> <div data-bbox="342 869 1393 900"> <p>1.5.6 Proceed to Section 2, Preparation of the BioMek® 2000 Automation Workstation.</p> </div>	

1 ISOLATING DNA USING THE DNA IQ™ ISOLATION SYSTEM		Page 6 of 8
BIOMEK® 2000 AUTOMATION WORKSTATION PROCEDURES		Issue No. : 4
MANUAL- FORENSIC BIOLOGY SECTION PROCEDURE		Effective Date: 8-August-2005
MANUAL, SECTION IV		
1.6	DNA IQ™ EXTRACTION METHOD FOR MIXED BODY FLUID STAINS (DIFFERENTIAL PROCEDURE)	
1.6.1	Cut an approximate 3 mm ² stain or a portion of a cotton swab proportional to the number of sperm identified. Place the stain into a labeled 1.5 mL microcentrifuge tube with a depression in the lid. If the stain is smeared over a large area on the fabric or the body fluids on a swab are weak and more than one swab must be used, more than one microcentrifuge tube may be needed to extract the stain.	
1.6.1.1	If a sample contains a weak smear over a large surface area the sample should be placed in several 1.5 mL microcentrifuge tubes and the entire sample condensed into one sample tube during the purification step.	
NOTE: Process a reagent blank along with each set of samples.		
1.6.2	Add:	
	400 µL TNE 25 µL 20% Sarkosyl 75 µL Sterile Type I Water 5 µL Proteinase K	
	in proportional amounts to saturate the cutting.	
1.6.3	Mix by hand or light vortexing then pulse spin to force the cutting into the liquid.	
1.6.4	Place the tube into a 37°C incubator or heat block for a minimum of 2 hours.	
1.6.5	Pulse spin the tube, punch 2-3 holes in the lid of the tube, remove the cutting from the liquid and place in the lid. Spin the tube for 5 minutes in a microcentrifuge at a minimum of 10,000 rpm to remove the excess liquid from the cutting. Alternatively, a Spin-Ease basket may be used instead of placing the cutting in the lid of the tube.	
1.6.6	Using a pipette, carefully transfer all but approximately 50 µL of the supernatant into a new 1.5 mL labeled tube with a lid. Be careful not to dislodge or disturb the pellet on the bottom of the tube. The supernatant removed from the pellet is the NON-SPERM FRACTION.	
1.6.7	At this stage set the non-sperm fraction tube aside and wait until the sperm fraction is ready, and then proceed to Section 2, Preparation of the BioMek® 2000 Automation Workstation, with both the sperm and non-sperm fractions.	
1.6.8	Remove and discard the old lid and cutting from the original tube containing the pellet. (If the entire biological sample deposited on the evidence was consumed the cutting will be returned with the evidence). Place a new colored lid on the tube. This tube contains the SPERM FRACTION.	
1.6.9	Wash the pellet as follows: Resuspend the pellet in 500 µL of PCR digestion buffer by vortexing briefly. Spin the tube for 5 minutes in a microcentrifuge at a minimum of 10,000	

1 ISOLATING DNA USING THE DNA IQ™ ISOLATION SYSTEM	Page 7 of 8
BIOMEK® 2000 AUTOMATION WORKSTATION PROCEDURES MANUAL- FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION IV	Issue No. : 4
	Effective Date: 8-August-2005
<p>rpm. Using a pipette with a sterile 1 mL pipette tip, remove all but 50 µL of the supernatant and discard.</p> <p>1.6.10 Repeat the wash in step 1.5.16 an additional 2 times. If a low sperm count has been determined, the sperm pellet may be washed up to 5 times. After the final wash, remove all but 50 µL of the wash buffer and discard. The sperm cells are not lysed at this point. The sperm cells are lysed during the isolation step using the BioMek® 2000 Automation Workstation</p> <p>NOTE: If necessary the sample(s) once capped can remain at room temperature overnight or be stored in a refrigerator before proceeding to Section 2. If stored in a refrigerator place the samples in a 56°C heat block for 5 minutes to resolubilize the DNA before proceeding.</p> <p>1.6.11 Proceed to Section 2, Preparation of the BioMek® 2000 Automation Workstation, with both the sperm and non-sperm fractions.</p> <p>OPTION: If a sperm search has not previously been conducted, remove 3 µL of the supernatant and spot the sample onto a glass microscope slide for examination.</p>	

1 ISOLATING DNA USING THE DNA IQ™ ISOLATION SYSTEM	Page 8 of 8
BIOMEK®2000 AUTOMATION WORKSTATION PROCEDURES MANUAL- FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION IV	Issue No. : 4
	Effective Date: 8-August-2005
<p>1.7 DNA IQ™ EXTRACTION METHOD FOR HAIR AND LOW LEVEL SAMPLES</p> <p>NOTE: This procedure may also be used for samples that are believed to have a low concentration of DNA (i.e., envelopes, stamps, cigarette butts) or highly concentrated bloodstain (i.e., dried blood flakes).</p> <p>1.7.1 Follow the steps outlined below for the extraction of hairs. Proceed to step 1.6.2 for all other sample types:</p> <p>1.7.1.1 Using clean tweezers, place the hair on a clean piece of white or black paper, whichever is appropriate, and examine under a stereo microscope for the presence of sheath material. Note the presence of any body fluids on the hair.</p> <p>1.7.1.2 Wash the hair to reduce surface dirt and contaminants by immersing the hair in sterile Type I Water in a clean 50 mL beaker. If the hair contains a biological fluid that is important to the investigation DO NOT wash the hair.</p> <p>1.7.1.3 Return the hair to the stereo microscope. Use a clean scalpel blade to cut a 0.5 to 1 cm portion from the root end of the hair and then place the hair root into a 1.5 mL microcentrifuge tube and proceed to step 1.6.3.</p> <p>NOTE: Process a reagent blank along with each set of samples.</p> <p>1.7.2 Cut an approximate 3 mm² blood stain or other biological stain and place into a labeled 1.5 mL microcentrifuge tube with a depression in the lid</p> <p>1.7.3 Add to the 1.5 mL microcentrifuge tube:</p> <p>37.5 µL TNE 12.5 µL 20% Sarkosyl 10.0 µL 0.39M DTT 32.5 µL Sterile Type I Water 10.0 µL Proteinase K</p> <p>in proportional amounts to saturate the cutting. NOTE: the maximum volume of sample loaded into the deep well plate is 100 µL.</p> <p>1.7.4 Mix by hand or lightly vortex, then pulse spin the microcentrifuge tube to force the sample into the liquid.</p> <p>1.7.5 Place the tube into a 56°C incubator or heat block for minimum of 1 hour.</p> <p>1.7.6 Pulse spin the tube in a microcentrifuge for 10 seconds to force the condensate to the bottom of the tube.</p> <p>1.7.7 Proceed to Section 2, Preparation of the BioMek® 2000 Automation Workstation.</p> <p style="text-align: right;">◆END</p>	