

APPENDIX E: CLEANING THE LUMINOMETER LENS	Page 1 of 1
BIOMEK®2000 AUTOMATION WORKSTATION PROCEDURES MANUAL-FORENSIC BIOLOGY SECTION PROCEDURE MANUAL, SECTION IV	Issue No.: 4
	Effective Date: 8-August-2005
<p>APPENDIX E: CLEANING THE LUMINOMETER LENS</p> <ol style="list-style-type: none"> 1. Turn off the Luminometer and unplug the unit. 2. Unscrew the four (4) Allen screws in the back of the Luminometer and take off the top cover. 3. Unscrew the four (4) finger screws and remove the cover from the optical unit. 4. Unscrew the three (3) long finger screws that hold the optical unit to the base. 5. Unplug the multi-colored ribbon cable from the optical unit. 6. Remove the optical unit from the base. The lens is on the underside of the optical unit. 7. Clean the lens first with lens paper. If lens is still dirty use 70% isopropyl alcohol on a swab to clean the lens. <u>DO NOT</u> use 70% EtOH or 95% EtOH to clean the lens. 8. Remove the light shield from the base (it sits in the hole directly underneath the lens when the lens is attached to the base). 9. Clean the light shield with 70% Isopropyl alcohol. 10. Make sure that the lens and light shield are dry before replacing them. 11. Replacement is the reverse of the removal. <p style="text-align: right;">◆END</p>	