

6 DETERMINATION OF SPECIES ORIGIN	Page 1 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<div data-bbox="151 296 716 331"> <p>6 DETERMINATION OF SPECIES ORIGIN</p> </div> <div data-bbox="199 363 354 399"> <p>6.1 GOALS</p> </div> <div data-bbox="256 430 1485 667"> <p>6.1.1 To acquire a basic understanding of immunology, including the theory and procedures for species origin determination.</p> <p>6.1.2 To acquire a thorough understanding of the use of controls.</p> <p>6.1.3 To become acquainted with the specificity, sensitivity, and limitations of the Ouchterlony double diffusion method.</p> </div> <div data-bbox="199 699 349 735"> <p>6.2 TASKS</p> </div> <div data-bbox="256 766 1539 1717"> <p>6.2.1 Determine the sensitivity of anti-human precipitin serum by testing various dilutions in Normal Saline (neat, 1:5, 1:10, 1:20, 1:50, 1:100) of known human blood and bloodstains using Ouchterlony double diffusion.</p> <div data-bbox="354 898 873 1077"> <p>6.2.1.1 Normal Saline (0.9%):</p> <ul style="list-style-type: none"> • 9 g Sodium chloride • 1000 ml Distilled water • Mix thoroughly until dissolved. </div> <p>6.2.2 Test at least 15 bloodstains subjected to various environmental conditions (heat, moisture, heat and moisture combined), and decomposition using Ouchterlony double diffusion.</p> <p>6.2.3 Test at least 15 bloodstains exposed to various contaminants (including, but not limited to superglue, fingerprint powder, ninhydrin, redwop power - rhodamine base, bleach, soap, motor oil, luminol, and mold) using Ouchterlony double diffusion.</p> <p>6.2.4 Test the precipitating antisera in the routine species collection addressed in 6.4.6 for cross-reactivity using Ouchterlony double diffusion. Observe specificity.</p> <p>6.2.5 Stain Ouchterlony plates using Coomassie Brilliant Blue R250. Note enhancement of weak reactions.</p> <p>6.2.6 Instruction by and observation of qualified examiners performing routine examinations of case material.</p> <p>6.2.7 Test at least 10 unknown stains provided by training coordinator or designee.</p> <p>6.2.8 Read applicable literature. Refer to Appendix A and Appendix B.</p> </div> <div data-bbox="199 1749 592 1785"> <p>6.3 TRAINING EVALUATION</p> </div> <div data-bbox="256 1816 1185 1917"> <p>6.3.1 Knowledge</p> <div data-bbox="354 1879 1185 1917"> <p>6.3.1.1 Review of notes in training notebook by training coordinator.</p> </div> </div>	

6 DETERMINATION OF SPECIES ORIGIN	Page 2 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>6.3.1.2 Mini-mock trials/oral and practical examinations.</p> <p>6.3.1.3 Completion of checklist by training coordinator.</p> <p>6.3.2 Skills</p> <p>6.3.2.1 Observation by training coordinator or designee.</p> <p>6.3.2.2 Review of notes in training notebook by training coordinator.</p> <p>6.3.2.3 Mini-mock trials/oral and practical examinations.</p> <p>6.3.2.4 Completion of checklist by training coordinator.</p> <p>6.4 SPECIES IDENTIFICATION – TECHNICAL NOTES</p> <p>6.4.1 The procedure described in this section requires the use of precipitating antiserum to determine species origin. This procedure is an immunological procedure for the identification of animal or human protein. The identification is made by comparing the reaction of the unknown protein (in the case sample) against a known antiserum with the reaction of a known protein against a known antiserum.</p> <p>6.4.2 A positive control (a known sample of normal serum or blood against which the antiserum is directed) and a substrate control (when available) must be run on each plate when the testing is performed on a case. If a substrate control is not available, distilled water will be used. The substrate control monitors for contaminating protein activity (which could cause a false positive reaction) in the unstained portion of the substrate as well as for contaminating protein activity in the reagent (distilled water) used for the extraction.</p> <p>6.4.3 Alternatively, a human DNA quantitation method may be used to determine that a sample is of human origin.</p> <p>6.4.4 Quality Control of Antiserums</p> <p>6.4.4.1 Before using any new lot number of precipitating antiserum for testing casework samples, the specificity must be tested and appropriately documented in the laboratory's quality control (QC) records. Routine quality control testing will be performed at the time the antiserum is reconstituted, or if the antiserum is received in liquid form, within one week of receipt.</p> <p>6.4.4.2 Anti-human serum, not anti-human hemoglobin, will be used with the procedure in this section for determining whether a sample is of human origin. Anti-human serum, as well as all animal antisera in the "Routine Species Collection" specified below, must be tested against all available species (normal or whole serums and known bloods) in the "Routine Species Collection".</p> <p>6.4.4.3 A positive control, a host control (typically normal rabbit serum or normal goat serum), and a negative control (distilled water) must be included in the specificity testing. The host control (representing the animal in which the antiserum was prepared) is used to demonstrate that the antiserum is not reacting to any proteins in the animal in which it was made.</p>	

6 DETERMINATION OF SPECIES ORIGIN	Page 3 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>6.4.4.4 The quality control documentation will include:</p> <p>6.4.4.4.1 A diagrammatic representation of the placement of samples in the wells of the rosette on the Ouchterlony plate with the subsequent observations (i.e., precipitin lines).</p> <p>6.4.4.4.2 Date of the testing.</p> <p>6.4.4.4.3 Initials of the person conducting the testing.</p> <p>6.4.4.4.4 Lot number, date of receipt, and manufacturer of the antiserum being tested.</p> <p>6.4.4.4.5 Lot number, date of receipt, and manufacturer of the normal serums being used for the testing.</p> <p>6.4.4.4.6 Results of the testing.</p> <p>6.4.4.5 Once the appropriate testing has been performed on a particular lot number of antiserum, it need not be repeated for each case. If another vial of the same lot number is received on a different date, the QC testing described above must be repeated.</p> <p>6.4.5 Quality Control of Normal (Whole) Serums</p> <p>6.4.5.1 Before using any new lot number of normal serum for testing casework samples, the specificity must be tested and appropriately documented in the laboratory's quality control (QC) records. Routine quality control testing will be performed at the time the normal serum is reconstituted, or if the normal serum is received in liquid form, within one week of receipt.</p> <p>6.4.5.2 All normal serums in the "Routine Species Collection" must be tested against all available antisera in the "Routine Species Collection".</p> <p>6.4.5.3 A positive control must be included in the specificity testing. Although normal human serum may be purchased, <u>a straw colored dilution</u> of known human blood may be used instead. Similarly, the use of known blood from other species may replace the purchase and use of normal serums from the species.</p> <p>6.4.5.4 It is not necessary to conduct quality control testing on the known bloods. Label known bloods with the species name and date of preparation/initials of person preparing the sample. Store known bloods in the freezer.</p> <p>6.4.5.5 The quality control documentation will include:</p> <p>6.4.5.5.1 A diagrammatic representation of the placement of samples in the wells of the rosette on the Ouchterlony plate with the subsequent observations (i.e., precipitin lines).</p> <p>6.4.5.5.2 Date of the testing.</p>	

6 DETERMINATION OF SPECIES ORIGIN	Page 4 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>6.4.5.5.3 Initials of the person conducting the testing.</p> <p>6.4.5.5.4 Lot number, date of receipt, and manufacturer of the normal serum being tested.</p> <p>6.4.5.5.5 Lot number, date of receipt, and manufacturer of the antiserums being used for the testing.</p> <p>6.4.5.5.6 Results of the testing</p> <p>6.4.5.6 Once the appropriate testing has been performed on a particular lot number of normal serum, it need not be repeated for each case. If another vial of the same lot number is received on a different date, the QC testing described above must be repeated.</p> <p>6.4.6 Routine Species Collection</p> <p>6.4.6.1 The following will be maintained in each laboratory as the “Routine Species Collection” and will undergo QC testing for specificity as outlined above:</p> <p>6.4.6.1.1 Bovine antiserum and normal bovine serum or known blood</p> <p>6.4.6.1.2 Swine antiserum and normal swine serum or known blood</p> <p>6.4.6.1.3 Cat antiserum and normal cat serum or known blood</p> <p>6.4.6.1.4 Dog antiserum and normal dog serum or known blood</p> <p>6.4.6.1.5 Rabbit antiserum and normal rabbit serum or known blood</p> <p>6.4.6.1.6 Sheep antiserum and normal sheep serum or known blood</p> <p>6.4.6.1.7 Deer antiserum and normal deer serum or known blood</p> <p>6.4.6.1.8 Human antiserum and normal human serum or known blood</p> <p>6.4.6.1.9 Normal goat serum or known blood (goat antiserum is unavailable)</p> <p>NOTE: It is recommended that the antiserums and normal serums listed above be reconstituted and QC tested at the time of receipt to ensure ready availability.</p> <p>6.4.6.2 Other antiserums/normal serums, such as bear, rodent, fowl, horse, etc. should also be maintained for use in special cases and must also undergo QC testing as specified above. However, since these are used only in special cases, it is recommended that they not be reconstituted, aliquoted, and QC tested until it is determined that there is a specific need to do so.</p> <p>6.4.7 Storage of Antiserum/Normal Serum</p> <p>6.4.7.1 Small aliquots of the antiserum/normal serum will be prepared for routine use and frozen within one week of reconstitution (when antiserum is lyophilized) or upon receipt (when antiserum is liquid). All frozen aliquots have an indefinite expiration date.</p>	

6 DETERMINATION OF SPECIES ORIGIN	Page 5 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>6.4.7.2 A thawed aliquot may be stored refrigerated for up to 1 month. If the aliquot is to be maintained in this manner, the expiration date must be clearly marked on the vial. Otherwise, the aliquot will be immediately discarded following its use.</p> <p>6.4.8 Labeling of Antiserum/Normal Serum</p> <p>6.4.8.1 Labels on each aliquot will include:</p> <p>6.4.8.1.1 The manufacturer</p> <p>6.4.8.1.2 Type of antiserum or normal serum</p> <p>6.4.8.1.3 Lot number</p> <p>6.4.8.1.4 Date reconstituted/date frozen</p> <p>6.4.8.1.5 Initials of the person preparing the aliquot</p> <p>6.4.9 OUCHTERLONY (DOUBLE DIFFUSION) TEST (Reference 6, pp. 221-241, Appendix B)</p> <p>6.4.9.1 Equipment</p> <p>6.4.9.1.1 Punch</p> <p>6.4.9.1.2 Aspirator</p> <p>6.4.9.1.3 100 ml and 500 ml graduated cylinders</p> <p>6.4.9.1.4 Balance</p> <p>6.4.9.1.5 Spatula</p> <p>6.4.9.1.6 Scissors</p> <p>6.4.9.1.7 Tweezers</p> <p>6.4.9.1.8 Hot plate or oven (37° C)</p> <p>6.4.9.1.9 Incubator (optional)</p> <p>6.4.9.1.10 Magnetic stir plate</p> <p>6.4.9.1.11 Refrigerator (optional)</p> <p>6.4.9.2 Materials</p> <p>6.4.9.2.1 Petri dishes, slides, or comparable containers</p> <p>6.4.9.2.2 Test tubes</p>	

6 DETERMINATION OF SPECIES ORIGIN	Page 6 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<div data-bbox="467 296 930 331">6.4.9.2.3 Weigh boat or weigh paper</div> <div data-bbox="467 363 824 399">6.4.9.2.4 Moisture chamber</div> <div data-bbox="467 430 818 466">6.4.9.2.5 Disposable pipets</div> <div data-bbox="467 497 789 533">6.4.9.2.6 Capillary tubes</div> <div data-bbox="354 564 574 600">6.4.9.3 Reagents</div> <div data-bbox="467 632 862 667">6.4.9.3.1 Normal saline (0.9%)</div> <div data-bbox="467 699 818 735">6.4.9.3.2 Agarose gel (1%)</div> <div data-bbox="467 766 784 802">6.4.9.3.3 Distilled water</div> <div data-bbox="467 833 734 869">6.4.9.3.4 Antiserum</div> <div data-bbox="467 900 1179 936">6.4.9.3.5 Normal serum or known blood (positive control)</div> <div data-bbox="354 968 756 1003">6.4.9.4 Agarose Gel Preparation</div> <div data-bbox="467 1035 940 1071">6.4.9.4.1 Normal saline (0.9% NaCl):</div> <div data-bbox="610 1102 1008 1209"> <ul style="list-style-type: none"> • 9 g Sodium chloride • 1000 ml Distilled water • Mix thoroughly until dissolved. </div> <div data-bbox="467 1241 824 1276">6.4.9.4.2 Agarose gel (1%):</div> <div data-bbox="610 1308 1279 1415"> <ul style="list-style-type: none"> • 1 g Type I agarose • 100 ml Normal saline (0.9% NaCl) • Heat until agarose is dissolved. Allow to cool slightly. </div> <div data-bbox="467 1446 1544 1516">6.4.9.4.3 Pore the agarose into a petri dish, onto a slide, or into a comparable container to a thickness of 2-3 mm and allow to cool.</div> <div data-bbox="467 1547 1511 1617">6.4.9.4.4 Cut wells in a rosette pattern (refer to the diagram below) into the gel using a punch or disposable pipet connected to an aspirator.</div> <div data-bbox="724 1631 990 1913"> </div>	

6 DETERMINATION OF SPECIES ORIGIN	Page 7 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>6.4.9.5 Storage and Labeling</p> <p>6.4.9.5.1 When a batch of plates is prepared, the plates should be numbered consecutively and placed in a moisture chamber in the refrigerator. Label the moisture chamber with the lot number of the batch (date of preparation/initials of person preparing the plates).</p> <p>6.4.9.5.2 There is no expiration date (see 6.4.9.6 Minimum Standards and Controls).</p> <p>6.4.9.6 Minimum Standards and Controls</p> <p>6.4.9.6.1 A positive control (known sample against which the antiserum is directed) and a substrate control (or if not available, distilled water) must be tested on each plate, unless the stain is on a cotton swab. It is not necessary to test submitted control swabs.</p> <p>6.4.9.7 OUCHTERLONY DOUBLE DIFFUSION PROCEDURE</p> <p>6.4.9.7.1 To prepare an extract of the stain, place a small cutting of the stain in distilled water until a <u>straw color</u> is obtained. A small piece of stained material, which is moistened with distilled water, can be used in lieu of an extract. Treat the substrate control in the same manner as the stain.</p> <p>6.4.9.7.2 Add antiserum in the center well of the Ouchterlony plate with a disposable pipet or capillary tube.</p> <p>6.4.9.7.3 Add appropriate extracts/pieces of stained material, the positive control, and negative control(s) to the surrounding wells. Do not overfill the wells. Avoid getting bubbles in the wells. Document the placement of samples in the wells of the rosettes on the Ouchterlony work sheet found near the end of this chapter.</p> <p>6.4.9.7.4 Record the lot numbers of antisera and normal sera used for the testing procedure.</p> <p>NOTE: Alternatively, known normal serum/known blood extract may be placed in the center well with appropriate antisera in the surrounding wells or the stain extract/piece of stained material may be placed in the center well with appropriate antisera in the surrounding wells.</p> <p>6.4.9.7.5 Incubate the plate in a moisture chamber at 37° C for 3-4 hours. Alternatively, it may be left overnight at room temperature or 4° C.</p> <p>6.4.9.7.6 Record observations (precipitin lines) on the diagram, and interpret and record the results.</p> <p>6.4.9.7.7 All controls must give the expected results before a conclusion can be reached on an unknown sample, i.e., white precipitin lines must be observed between the antiserum and positive control (known serum) and no precipitin lines should be observed between the antiserum and the negative control.</p>	

6 DETERMINATION OF SPECIES ORIGIN	Page 8 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<div data-bbox="467 260 769 296">6.4.9.7.8 Interpretation</div> <div data-bbox="610 327 1549 737"> <div data-bbox="610 327 1549 401">6.4.9.7.8.1 Positive Result = White precipitin lines between the antiserum well and the sample well</div> <div data-bbox="610 432 1393 506">6.4.9.7.8.2 Negative Result = No precipitin lines between the antiserum well and the sample well</div> <div data-bbox="610 537 1549 737">6.4.9.7.8.3 Inconclusive Result = Questionable precipitin lines between the antiserum well and the sample well, OR precipitin lines between the antiserum well and the negative control well(s), OR no precipitin line between the antiserum well and the positive control well. If sufficient sample remains, an inconclusive result should be repeated.</div> </div> <div data-bbox="293 768 1511 936"> <p>NOTE: The prozone phenomenon can result in a soluble antigen-antibody complex due to too many antibodies present to form a complete lattice (Reference 8, Appendix B). Because of this phenomenon, weak precipitin lines may be observed initially, but disappear upon staining with Coomassie Brilliant Blue R250. This is considered an inconclusive result. Therefore, it is recommended that the testing results be recorded PRIOR TO staining as well as after staining.</p> </div> <div data-bbox="467 968 1490 1073"> <div data-bbox="467 968 1490 1073">6.4.9.7.9 Staining the plate with Coomassie Brilliant Blue R250 may be necessary to visualize weak reactions. Refer to 6.4.10 for the Coomassie Brilliant Blue R250 staining procedure.</div> </div> <div data-bbox="467 1104 821 1140">6.4.9.7.10 Reporting Results</div> <div data-bbox="610 1171 1549 1472"> <div data-bbox="610 1171 1549 1245">6.4.9.7.10.1 Report positive test results as “(species tested according to the label on the antiserum) protein was detected...”</div> <div data-bbox="610 1276 1549 1350">6.4.9.7.10.2 Report negative test results as “no (species tested according to the label on the antiserum) protein was detected...”</div> <div data-bbox="610 1381 1549 1472">6.4.9.7.10.3 Report inconclusive test results as “the test for (species tested according to the label on the bottle) protein was inconclusive...”</div> </div> <div data-bbox="256 1539 1495 1575">6.4.10 COOMASSIE BRILLIANT BLUE R250 STAINING PROCEDURE (Reference 13, Appendix B)</div> <div data-bbox="350 1606 724 1642">6.4.10.1 Safety Considerations</div> <div data-bbox="467 1673 1549 1913"> <div data-bbox="467 1673 1549 1778">6.4.10.1.1 Coomassie Brilliant Blue R250 - Caution! Avoid contact and inhalation! Emits toxic fumes under fire conditions! Container explosion may occur under fire conditions!</div> <div data-bbox="467 1810 1471 1845">6.4.10.1.2 Methanol - Caution! Irritant! Dangerous when exposed to heat or flame!</div> <div data-bbox="467 1877 1247 1913">6.4.10.1.3 Glacial acetic acid - Caution! Corrosive! Flammable!</div> </div>	

6 DETERMINATION OF SPECIES ORIGIN	Page 9 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>6.4.10.2 Equipment</p> <p>6.4.10.2.1 Weight</p> <p>6.4.10.2.2 Oven or Incubator (40-60°C)</p> <p>6.4.10.2.3 Rotator (optional)</p> <p>6.4.10.2.4 10 ml, 50 ml, and 500 ml graduated cylinders</p> <p>6.4.10.2.5 Balance</p> <p>6.4.10.2.6 Spatula</p> <p>6.4.10.2.7 Trays for staining and destaining</p> <p>6.4.10.3 Materials</p> <p>6.4.10.3.1 Gel bond, glass plate, or other support medium</p> <p>6.4.10.3.2 Weigh boats or weigh paper</p> <p>6.4.10.3.3 Whatman #1 filter paper</p> <p>6.4.10.3.4 Paper towels</p> <p>6.4.10.4 Reagents</p> <p>6.4.10.4.1 Staining Solution</p> <p>6.4.10.4.2 Destaining Solution</p> <p>6.4.10.4.3 Normal saline (0.9% NaCl) – refer to 6.4.10.8.1</p> <p>6.4.10.4.4 Distilled water</p> <p>6.4.10.5 Preparation of Stain and Destain Solutions</p> <p>6.4.10.5.1 Staining Solution</p> <ul style="list-style-type: none"> • 0.1 g Coomassie Brilliant Blue R250 • 45.0 ml Methanol • 10.0 ml Glacial acetic acid • 45.0 ml Distilled water • Mix the above ingredients until thoroughly dissolved. <p>6.4.10.5.2 Destaining Solution</p> <ul style="list-style-type: none"> • 45.0 ml Methanol • 10.0 ml Glacial acetic acid 	

6 DETERMINATION OF SPECIES ORIGIN	Page 10 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<div data-bbox="609 268 1266 336" data-label="List-Group"> <ul style="list-style-type: none"> • 45.0 ml Distilled water • Mix the above ingredients until thoroughly dissolved. </div> <div data-bbox="350 369 557 401" data-label="Section-Header"> <h4>6.4.10.6 Storage</h4> </div> <div data-bbox="462 434 1430 468" data-label="Text"> <p>6.4.10.6.1 The Staining and Destaining Solutions are stable at room temperature.</p> </div> <div data-bbox="350 501 573 533" data-label="Section-Header"> <h4>6.4.10.7 Labeling</h4> </div> <div data-bbox="462 567 1461 636" data-label="Text"> <p>6.4.10.7.1 Label as Staining or Destaining Solution with the lot number (date of preparation followed by the initials of the person preparing the solution).</p> </div> <div data-bbox="462 667 932 701" data-label="Text"> <p>6.4.10.7.2 There is no expiration date.</p> </div> <div data-bbox="350 735 1278 766" data-label="Section-Header"> <h4>6.4.10.8 COOMASSIE BRILLIANT BLUE R250 STAINING PROCEDURE</h4> </div> <div data-bbox="462 802 1510 869" data-label="Text"> <p>6.4.10.8.1 Wash the plate overnight in normal saline solution (0.9% NaCl - 9 g NaCl in 1000 ml distilled water) to remove unprecipitated proteins.</p> </div> <div data-bbox="605 903 1481 970" data-label="Text"> <p>6.4.10.8.1.1 If pieces of stained material were used in lieu of extracts, remove these prior to washing the plate.</p> </div> <div data-bbox="605 1003 1546 1104" data-label="Text"> <p>6.4.10.8.1.2 The gel may detach from the plate during the washing process. Mark the orientation of the gel to ensure that it can be re-oriented properly after the washing has been completed.</p> </div> <div data-bbox="605 1138 1546 1272" data-label="Text"> <p>6.4.10.8.1.3 Alternatively, the gel may be removed from the plate before beginning the washing process. If this is done, mark the gel to ensure that it can be re-oriented properly after the washing has been completed.</p> </div> <div data-bbox="462 1306 1531 1440" data-label="Text"> <p>6.4.10.8.2 The next day wash the gel for approximately fifteen minutes in distilled water. Rinse the gel and repeat the wash. Ensuring proper orientation, place the gel (face up) on the hydrophilic side of a piece of gel bond or on some other support medium such as a glass plate.</p> </div> <div data-bbox="462 1474 1534 1608" data-label="Text"> <p>6.4.10.8.3 Cover the gel with a piece of Whatman #1 filter paper moistened with distilled water. Add a layer of paper towels on top of the filter paper and press with a weight for approximately 30 minutes. Remove paper towels and filter paper and dry the gel in a 40-60° C oven.</p> </div> <div data-bbox="462 1642 1510 1740" data-label="Text"> <p>6.4.10.8.4 Place the gel in the staining solution and allow it to soak for 1 to 10 minutes. This may be done on a rotator. Intermittently check staining progress to prevent over staining.</p> </div> <div data-bbox="462 1774 1534 1875" data-label="Text"> <p>6.4.10.8.5 Place the gel in the destaining solution until the background is clear or until no more dye leaches from the gel. This may be done on a rotator. Change the destaining solution and destain further if desired.</p> </div>	

6 DETERMINATION OF SPECIES ORIGIN	Page 11 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p data-bbox="467 266 1528 432">6.4.10.8.6 Record observations (precipitin lines) on the diagram, and interpret and record the results. All controls must give the expected results before a conclusion can be reached on an unknown sample, i.e., blue precipitin lines must be observed between the antiserum and positive control (known serum) and no precipitin lines should be observed between the antiserum and the negative control.</p> <p data-bbox="467 468 769 499">6.4.10.8.7 Interpretation</p> <p data-bbox="610 535 1528 600">6.4.10.8.7.1 Positive Result = Blue precipitin lines between the antiserum well and the sample well.</p> <p data-bbox="610 636 1520 701">6.4.10.8.7.2 Negative Result = No precipitin lines between the antiserum well and the sample well.</p> <p data-bbox="610 737 1536 934">6.4.10.8.7.3 Inconclusive Result = Questionable precipitin lines between the antiserum well and the sample well, OR precipitin lines between the antiserum well and the negative control well(s), OR no precipitin line observed between the antiserum well and the positive control well. If sufficient sample remains, an inconclusive result should be repeated.</p> <p data-bbox="467 970 821 1001">6.4.10.8.8 Reporting Results</p> <p data-bbox="610 1037 1036 1068">6.4.10.8.8.1 Refer to 6.4.9.7.10.</p>	

6 DETERMINATION OF SPECIES ORIGIN

Page 12 of 15

**TRAINING MANUAL: CASE APPROACH AND
IDENTIFICATION OF BIOLOGICAL SUBSTANCES**

Amendment Designator:

Effective Date: 14-March-2006


**VIRGINIA DEPARTMENT OF FORENSIC SCIENCE
OUCHTERLONY FOR SPECIES DETERMINATION**

ANALYST: _____


FS LAB#: _____

DATE: _____

PLATE#



PLATE#



	SAMPLE	LOT# / SOURCE	RESULTS
1			
2			
3			
4			
5			
6			
7			

	SAMPLE	LOT# / SOURCE	RESULTS
1			
2			
3			
4			
5			
6			
7			

REAGENT	LOT# / SOURCE
Saline solution	
Coomassie blue staining solution	
Destaining solution	

COMMENTS:

6 DETERMINATION OF SPECIES ORIGIN	Page 13 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006

DETERMINATION OF SPECIES ORIGIN STUDY QUESTIONS

1. What is an antigen?

2. What are the conditions of antigenicity?

3. What is an antibody?

4. How is antiserum made? How is normal serum made?

5. What is meant by cross-reactivity?

6. What are the advantages of Ouchterlony double diffusion? Are there disadvantages?

7. Explain Ouchterlony double diffusion and the purpose of each component/control.

8. What does a positive result look like and what does it tell you?

9. When using the Coomassie Blue stain, what is the dye staining? Why is this used?

10. What is Prozone and Postzone?

11. What would you do if you had a very small possible bloodstain and the investigator wanted to know if it was human blood, but species testing would consume the stain?

12. Why is specificity testing of antisera and normal sera required prior to using these for testing case material?

6 DETERMINATION OF SPECIES ORIGIN	Page 14 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006

TRAINEE CHECKLIST FOR THE DETERMINATION OF SPECIES ORIGIN

Name of Trainee: _____

1. Completion of the following tasks using Ouchterlony double diffusion:

Tested known human blood and bloodstains of varying dilutions (neat, 1:5, 1:10, 1:20, 1:50, 1:100).

Date: _____ Training Coordinator: _____

Comments: _____

Tested bloodstains subjected to various environmental conditions (15 minimum).

Date: _____ Training Coordinator: _____

Comments: _____

Tested bloodstains exposed to various contaminants (15 minimum).

Date: _____ Training Coordinator: _____

Comments: _____

Tested the precipitating antisera in the routine species collection addressed in 6.4.6 for cross- reactivity.

Date: _____ Training Coordinator: _____

Comments: _____

Stained Ouchterlony plates using Coomassie Brilliant Blue R250.

Date: _____ Training Coordinator: _____

Comments: _____

Accurately tested at least 10 unknown stains (provided by training coordinator or designee).

Date: _____ Training Coordinator: _____

Comments: _____

6 DETERMINATION OF SPECIES ORIGIN	Page 15 of 15
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006

2. Trainee has developed a thorough understanding of the theory behind species origin determination, including basic immunology, specificity, sensitivity, and limitations of all methods, as well as the purpose and use of controls. An unquestionably sound technique has been developed for the use of Ouchterlony double diffusion.

Date: _____ Training Coordinator: _____

Comments: _____
3. Notebook is organized and complete.

Date: _____ Training Coordinator: _____

Comments: _____
4. Trainee has participated in a mock trial and/or practical or oral examinations. Performance was satisfactory.

Date: _____ Training Coordinator: _____

Comments: _____
5. Trainee has read and understands all applicable literature.

Date: _____ Training Coordinator: _____

Comments: _____

◆END