

7 SEMEN IDENTIFICATION	Page 1 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>7 SEMEN IDENTIFICATION</p> <p>7.1 GOALS</p> <p>7.1.1 To become proficient in the use of alternate light sources for locating semen stains.</p> <p>7.1.2 To learn the physical and chemical characteristics of semen (animal and human).</p> <p>7.1.3 To become proficient in extraction techniques, staining techniques, and microscopic examination for spermatozoa.</p> <p>7.1.4 To learn the theory behind the use of chemical (color) tests and immunological tests for semen.</p> <p>7.1.5 To become proficient in the use of the Acid Phosphatase Test and the p30 test by OneStep ABACard®, including the use of controls and possible sources of error.</p> <p>7.1.6 To develop an understanding of the sensitivity, specificity and limitations of the Acid Phosphatase Test (qualitative and quantitative) and p30 test by OneStep ABACard®.</p> <p>7.1.7 To be able to locate and evaluate stains on evidentiary material.</p> <p>7.1.8 To become proficient in techniques used to prevent cross-contamination of seminal fluid/spermatozoa between samples.</p> <p>7.2 TASKS</p> <p>7.2.1 Examine and compare at least 20 samples of known physiological fluids (including, but not limited to, different semen dilutions prepared in distilled water, blood, saliva, perspiration, and mixtures) and substances known to react to an alternate light source (including, but not limited to, milk, yogurt, lotion, and “bleach alternative” detergent) on different substrates with the aid of all alternate light sources available in the section.</p> <p>7.2.2 Examine several stained and unstained smears for spermatozoa using phase contrast microscopy and compare results.</p> <p>7.2.3 Perform presumptive and confirmatory tests, as appropriate, on a minimum of 50 known semen samples of varying ages, on various substrates, including mixtures and dilutions (neat to 1:100), and stains subjected to various contaminants and environmental conditions.</p> <p>7.2.4 Examine and compare 20 different prepared slides of animal spermatozoa in the reference collection.</p> <p>7.2.5 Perform presumptive and confirmatory tests, as appropriate, on at least 20 samples of various known physiological fluids, including different semen dilutions, mixtures, and aspermic semen samples.</p> <p>7.2.6 Test a minimum of 12 samples of varying dilutions of semen using the OneStep ABACard® p30 Test to determine the sensitivity of the p30 test. Compare results.</p> <p>7.2.7 Using presumptive and confirmatory tests, as appropriate, examine a series of unknown samples (25 minimum) for spermatozoa identification as provided by the training coordinator or designee. These</p>	

7 SEMEN IDENTIFICATION	Page 2 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>samples should consist of samples of varying dilutions of spermatozoa as well as samples with no spermatozoa. Use appropriate cleaning techniques between samples to ensure that no cross-contamination has occurred.</p> <p>7.2.8 Observe and obtain instruction from qualified examiners performing routine examinations of case material.</p> <p>7.2.9 Read applicable literature. Refer to Appendix A and Appendix B.</p> <p>7.3 TRAINING EVALUATION</p> <p>7.3.1 Knowledge</p> <p>7.3.1.1 Review of notes in training notebook by training coordinator.</p> <p>7.3.1.2 Mini-mock trials/oral and practical examinations.</p> <p>7.3.1.3 Completion of checklist by training coordinator.</p> <p>7.3.2 Skills</p> <p>7.3.2.1 Observation by training coordinator or designee.</p> <p>7.3.2.2 Review of notes in training notebook by training coordinator.</p> <p>7.3.2.3 Mini-mock trials/oral and practical examinations.</p> <p>7.3.2.4 Completion of checklist by training coordinator.</p> <p>7.4 TECHNICAL NOTES</p> <p>7.4.1 Screening items such as clothing or bedding for the presence of semen stains may be facilitated by the use of an alternate light source (ALS). Alternate light sources include a UV light (sometimes referred to as a “Wood’s Lamp” by Forensic Nurses), the Omnichrome FLS 5000, LumaLite™ 2000A, and Mini Crime Scope MCS400, to name a few. Users must read the directions accompanying each ALS in order to learn the best combination of wavelengths and filters, to avoid damaging the instrument during start up and shut down, and to protect their eyes from the powerful light. The use of appropriate goggles (dependent on the ALS) helps to make the reaction detectable to the eye, while simultaneously protecting the eyes from the light source. If proper eye protection is not worn, permanent damage to the eye may occur. The principle behind the light sources is that semen contains a component(s) which reacts to light between 450 and 455 nm wavelengths. While some sources cite flavins, other sources cite acid phosphatase as being the reactive component in semen. The reaction may either appear as a light stain against a dark background, or in some circumstances, the stain appears darker against a light background. The reaction must be interpreted with caution since other substances (such as, but not limited to, urine, saliva, makeup, yogurt, cleaners, bleach alternatives such as UV dyes) may also react to an ALS. Samples exhibiting a reaction to an ALS require further examination to detect and/or confirm the presence of semen.</p>	

7 SEMEN IDENTIFICATION	Page 3 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>7.4.2 When the presence of semen is suspected in a stain, the Acid Phosphatase Test, a preliminary chemical test used to screen stains for the presence of semen, is conducted initially. This test is based on the detection of acid phosphatase, a major component of semen. In the presence of acid phosphatase, the sodium α-naphthyl acid phosphate is hydrolyzed to α-naphthol, which diazotizes with the dye to yield a colored azo-dye. Samples giving a positive reaction to the screening test require further examination to confirm the presence of semen.</p> <p>7.4.3 Although a positive result with the Acid Phosphatase Test is strongly indicative of semen, confirmation of its presence must be established by the identification of spermatozoa or, in the absence of spermatozoa, the detection of p30, a human seminal plasma protein. In general, the presence of semen on swabs from a Physical Evidence Recovery Kit is confirmed by the finding of spermatozoa on the correspondingly labeled smears. Acid Phosphatase testing is optional when the correspondingly labeled smears are positive for spermatozoa. The presence of semen in stains is confirmed by the finding of spermatozoa in an extract of the stain. If the acid phosphatase test suggests the presence of semen, but no spermatozoa are identified on the correspondingly labeled smears or in an extract of the stain, semen may be confirmed by the identification of p30.</p> <p>7.5 ACID PHOSPHATASE TEST (Reference 6, pp. 162-163, Appendix B)</p> <p>7.5.1 Safety Considerations</p> <p>7.5.1.1 Glacial acetic acid - Caution! Corrosive! Flammable!</p> <p>7.5.1.2 Sodium acetate - Caution! Irritant!</p> <p>7.5.1.3 Sodium α-naphthyl acid phosphate - Caution! Irritant! Emits toxic fumes under fire conditions!</p> <p>7.5.1.4 o-Dianisidine (Naphthanil diazo blue B) - Caution! Highly toxic! Emits toxic fumes under fire conditions!</p> <p>7.5.1.5 Naphthanil diazo red - Caution! Avoid contact and inhalation! Emits toxic fumes under fire conditions!</p> <p>7.5.2 Equipment</p> <p>7.5.2.1 5 ml and 500 ml Graduated cylinders</p> <p>7.5.2.2 Balance</p> <p>7.5.2.3 Spatula</p> <p>7.5.2.4 Scissors</p> <p>7.5.2.5 Tweezers</p> <p>7.5.3 Materials</p> <p>7.5.3.1 Filter paper or microtiter plate (optional)</p>	

7 SEMEN IDENTIFICATION	Page 4 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<div data-bbox="354 296 797 331">7.5.3.2 Weigh boats or weigh paper</div> <div data-bbox="354 363 623 399">7.5.3.3 Cotton swabs</div> <div data-bbox="354 430 704 466">7.5.3.4 Test tubes or bottles</div> <div data-bbox="354 497 911 533">7.5.3.5 Disposable transfer pipets or droppers</div> <div data-bbox="256 564 565 600">7.5.4 Working Solutions</div> <div data-bbox="354 632 824 667">7.5.4.1 Acid Phosphatase (AP) Buffer</div> <div data-bbox="467 699 1122 842"> <ul style="list-style-type: none"> • 2.5 ml Glacial acetic acid • 10.0 g Sodium acetate (anhydrous) • 450.0 ml Distilled water • Mix the above ingredients until thoroughly dissolved. </div> <div data-bbox="467 873 699 909">7.5.4.1.1 Storage</div> <div data-bbox="610 940 1289 976">7.5.4.1.1.1 The AP Buffer is stable at room temperature.</div> <div data-bbox="467 1008 716 1043">7.5.4.1.2 Labeling</div> <div data-bbox="610 1075 1544 1241"> <div data-bbox="610 1075 1544 1178">7.5.4.1.2.1 Label the bottle as AP Buffer with a lot number (the date of preparation followed by the initials of the person preparing the stock solution). Example: AP Buffer Lot Number 100899JD was prepared by Jane Doe on October 8, 1999.</div> <div data-bbox="610 1272 1544 1344">7.5.4.1.2.2 There is no expiration date (see 7.5.5 Minimum Standards and Controls).</div> </div> <div data-bbox="354 1377 997 1413">7.5.4.2 Sodium α-Naphthyl Acid Phosphate Solution</div> <div data-bbox="467 1444 1425 1547"> <div data-bbox="467 1444 1425 1547">7.5.4.2.1 Add a small amount (approximately 4 mg) of sodium α-naphthyl acid phosphate to approximately 3 ml of Acid Phosphatase buffer in an appropriately labeled 10 X 75 mm test tube or bottle.</div> <div data-bbox="467 1579 1094 1614">7.5.4.2.2 Discard the solution at the end of the day.</div> </div> <div data-bbox="354 1646 623 1682">7.5.4.3 Dye Solution</div> <div data-bbox="467 1713 1533 1879"> <div data-bbox="467 1713 1533 1816">7.5.4.3.1 Add a small amount (approximately 4 mg) of o-dianisidine or naphthanil diazo red to approximately 3 ml of Acid Phosphatase buffer in an appropriately labeled 10 X 75 mm test tube or bottle.</div> <div data-bbox="467 1848 1094 1883">7.5.4.3.2 Discard the solution at the end of the day.</div> </div> <div data-bbox="354 1915 639 1950">7.5.4.4 Distilled water</div>	

7 SEMEN IDENTIFICATION	Page 5 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>7.5.5 Minimum Standards and Controls</p> <p>7.5.5.1 Test a positive reagent control (known semen stain) and a negative reagent control (distilled water) to ensure that the reagents are working properly. The results of this testing must be documented in the case file.</p> <p>7.5.5.2 If either control does not give the expected result, do not proceed with testing evidence samples until the problem has been resolved as demonstrated by testing another set of positive and negative reagent controls and achieving the expected results with both controls.</p> <p>7.5.5.3 If the results of the test are positive, a substrate control (if available) must also be tested, unless the stain is on a cotton swab, and the results of the testing documented in the case file. It is not necessary to test submitted control swabs.</p> <p>7.5.6 Acid Phosphatase (AP) Test Procedure</p> <p>7.5.6.1 Moisten filter paper/swab with distilled water. (Do not use AP buffer solution, as this will contaminate the stained area.) Press the filter paper against the suspected stain or gently rub the stained area with the moistened swab. Alternatively, a small piece of the stain/swab can be placed on filter paper, in a small test tube, or in a microtiter plate. Treat the substrate control in the same manner.</p> <p>7.5.6.2 Add 1-2 drops of sodium α-naphthyl acid phosphate solution.</p> <p>7.5.6.3 Add 1-2 drops of dye solution.</p> <p>7.5.6.4 The development of a blue/purple color with o-dianisidine or an orange/red color with naphthanil diazo red within 10 to 15 seconds is indicative of acid phosphatase levels in the semen range. Although the development of a pink/peach color may be observed with o-dianisidine, this is not indicative of seminal acid phosphatase and therefore, is not considered a positive reaction.</p> <p>7.5.6.5 The presence of semen in all samples exhibiting an inconclusive result or a positive result must be confirmed by identifying spermatozoa or, in the absence of spermatozoa, p30.</p> <p>7.5.6.6 Interpretation</p> <p>7.5.6.6.1 Positive Reaction = Blue/purple color with o-dianisidine within 10 to 15 seconds</p> <p style="text-align: center;">OR</p> <p>Orange/red color with naphthanil diazo red within 10 to 15 seconds</p> <p>7.5.6.6.2 Negative Reaction = No color development, slight/slow color development</p> <p>7.5.6.6.3 Inconclusive Reaction = Slow moderate to strong color development</p> <p>7.5.6.7 Refer to section 7.11 for reporting results.</p>	

7 SEMEN IDENTIFICATION	Page 6 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>7.6 EXTRACTION OF SPERMATOOZOA FROM A SUBSTRATE</p> <p>7.6.1 Equipment</p> <p>7.6.1.1 Rotator, vortex, sonicator, or centrifuge (depending on extraction method used)</p> <p>7.6.1.2 Scissors</p> <p>7.6.1.3 Tweezers</p> <p>7.6.1.4 Dissecting needle (optional)</p> <p>7.6.2 Materials</p> <p>7.6.2.1 Microscope slides</p> <p>7.6.2.2 Test tubes</p> <p>7.6.3 Reagents</p> <p>7.6.3.1 Distilled water</p> <p>7.6.4 Extraction Methods</p> <p>7.6.4.1 Cut a small portion of a stain and soak in a test tube overnight in distilled water.</p> <p>7.6.4.2 Soak a small portion of a stain in distilled water and rotate overnight.</p> <p>7.6.4.3 Soak a small portion of a stain in distilled water and sonicate for 10 seconds, followed by a 30 second sonication.</p> <p>7.6.4.4 Tease fibers apart and soak in a small amount of distilled water.</p> <p>7.6.4.5 Soak a small portion of a stain in distilled water and vortex.</p> <p>7.6.4.6 Soak a small portion of a stain in distilled water on a microscope slide, stain side down (may be followed by mastication).</p> <p>7.6.4.7 Cut the stain into small pieces, place the pieces on a microscope slide, and soak in a small amount of distilled water (may be followed by mastication).</p> <p>7.6.4.8 For the OneStep ABACard® p30 Test extraction method, refer to 7.11.6.1 through 7.11.6.7.</p> <p>NOTES: Always soak the material first; prolong the soaking for difficult stains. Use the sonicator on low (high setting will disintegrate spermatozoa). To concentrate an extract, after soaking a small portion of a stain or swab, centrifuge and make a smear of the sediment. DNA extracts can also be used to search for spermatozoa.</p>	

7 SEMEN IDENTIFICATION	Page 7 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>7.7 KERNECHTROT-PICROINDIGOCARMINE STAIN (CHRISTMAS TREE STAIN) (Reference 7, p. 141, Appendix B)</p> <p>7.7.1 Kernechtrot-Picroindigocarmine (KPIC) differential biological stain is used to assist in the identification of spermatozoa. The solutions for this procedure can either be purchased (SERI) or prepared in-house.</p> <p>7.7.2 Safety Considerations</p> <p>7.7.2.1 Aluminum sulfate - Caution! Harmful by inhalation, in contact with skin, and if swallowed! Emits toxic fumes under fire conditions!</p> <p>7.7.2.2 Nuclear fast red - Caution! Irritant! Emits toxic fumes under fire conditions!</p> <p>7.7.2.3 Saturated picric acid solution - Caution! Toxic! Explosive when dry! Emits toxic fumes under fire conditions!</p> <p>7.7.2.4 Indigocarmine dye - Caution! Harmful if swallowed! Emits toxic fumes under fire conditions!</p> <p>7.7.3 Equipment</p> <p>7.7.3.1 Flame or heat block (37° C)</p> <p>7.7.4 Materials</p> <p>7.7.4.1 Fixative (optional)</p> <p>7.7.5 Reagents</p> <p>7.7.5.1 Kernechtrot staining solution (KS)</p> <p>7.7.5.2 Picroindigocarmine staining solution (PICS)</p> <p>7.7.5.3 Distilled water</p> <p>7.7.5.4 95% ethanol or methanol</p> <p>7.7.6 Stock Solutions (In-house Preparation)</p> <p>7.7.6.1 Equipment</p> <p>7.7.6.1.1 Filtration apparatus</p> <p>7.7.6.1.2 500 ml glass beakers</p> <p>7.7.6.1.3 Balance</p> <p>7.7.6.1.4 Spatula</p>	

7 SEMEN IDENTIFICATION		Page 8 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES		Amendment Designator:
		Effective Date: 14-March-2006
7.7.6.1.5	Glass rod	
7.7.6.1.6	Plastic bottles	
7.7.6.2	Materials	
7.7.6.2.1	Filter paper	
7.7.6.2.2	Weigh boats or weigh paper	
7.7.6.3	Reagents	
7.7.6.3.1	Aluminum sulfate	
7.7.6.3.2	Nuclear Fast Red	
7.7.6.3.3	Distilled water	
7.7.6.3.4	Picroindigocarmine dye	
7.7.6.3.5	Saturated picric acid solution (Purchase saturated solution. DO NOT PURCHASE DRY PRODUCT! See Safety Considerations, 7.7.2.3.)	
7.7.6.4	Kernechtrot Solution (KS)	
	<ul style="list-style-type: none"> • In a beaker dissolve 5 g of aluminum sulfate in 100 ml of hot distilled water. • Immediately add 0.1 g of Nuclear Fast Red and stir with a glass rod. • Allow to cool and filter through filter paper. 	
7.7.6.4.1	Storage	
7.7.6.4.1.1	The Kernechtrot Solution is stable at room temperature for up to 6 months, but <u>may need to be refiltered after standing</u> .	
7.7.6.4.2	Labeling	
7.7.6.4.2.1	Label the bottle as KS with the expiration date and a lot number (the date of preparation followed by the initials of person preparing the solution). Example: KS Lot Number 100899JD was prepared by Jane Doe on October 8, 1999.	
7.7.6.5	Picroindigocarmine Solution (PICS)	
	<ul style="list-style-type: none"> • Dissolve 1 g of Indigocarmine dye in 300 ml of a commercially purchased saturated solution of picric acid. • Filter and store. 	

7 SEMEN IDENTIFICATION		Page 9 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES		Amendment Designator:
		Effective Date: 14-March-2006
<p>7.7.6.5.1 Storage</p> <p>7.7.6.5.1.1 The Picroindigocarmine Solution is stable at room temperature for up to 6 months, but <u>may need to be refiltered after standing.</u></p> <p>7.7.6.5.2 Labeling</p> <p>7.7.6.5.2.1 Label the bottle as PICS with an expiration date and a lot number (the date of preparation followed by the initials of person preparing the solution). Example: PICS Lot Number 100899JD was prepared by Jane Doe on October 8, 1999.</p> <p>7.7.7 SERI Christmas Tree Stain (R540) Kit</p> <p>7.7.7.1 Contents</p> <p>7.7.7.1.1 Solution A (Kernechtrot Solution - KS) - 30 ml</p> <p>7.7.7.1.2 Solution B (Picroindigocarmine Solution - PICS) - 30 ml</p> <p>7.7.7.1.3 Directions for use.</p> <p>7.7.7.2 Store under refrigeration in bottles provided.</p> <p>7.7.7.3 Shelf life: 6 months</p> <p>7.7.8 KPICS/Christmas Tree Staining Procedure</p> <p>7.7.8.1 Prepare a thin smear of an extract of a suspected semen stain and allow to dry, or examine a smear from the Physical Evidence Recovery Kit (PERK). Fix the smear with a quick flame or fixative, or by placing it on a 37° C heat block overnight.</p> <p>7.7.8.2 Add a sufficient amount (2-5 drops) of KS (red reagent) to cover the stained portion of the microscope slide.</p> <p>7.7.8.3 Let the slide stand at room temperature for at least 15 minutes.</p> <p>7.7.8.4 Wash KS off of the slide with a gentle stream of distilled water and drain the slide.</p> <p>7.7.8.5 Add a sufficient amount (2-5 drops) of PICS (green reagent) to cover the stained portion of the slide.</p> <p>7.7.8.6 Allow PICS to stain the smear for 5-15 seconds.</p> <p>7.7.8.7 Wash PICS off of the slide with 95% ethanol or methanol.</p> <p>7.7.8.8 Dry the slide at room temperature.</p>		

7 SEMEN IDENTIFICATION	Page 10 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>7.8 MICROSCOPIC EXAMINATION OF STAINED SLIDES FOR SPERMATOOZOA</p> <p>7.8.1 Equipment</p> <p>7.8.1.1 Microscope (with approximately 200X – 400X total magnification, with or without phase capability)</p> <p>7.8.2 Materials</p> <p>7.8.2.1 Distilled water, xylene substitute, or other appropriate mounting medium</p> <p>7.8.2.2 Coverslips</p> <p>7.8.3 Procedure for the Microscopic Examination of Stained Slides for Spermatozoa</p> <p>7.8.3.1 Quickly scan at approximately 200X total magnification. Confirm at approximately 400X total magnification.</p> <p>7.8.3.1.1 With phase microscopy: Spermatozoa heads are neon-like pink/red with darker pink/purple acrosomal caps and green tails. Epithelial cells and most bacteria stain green with some of the nuclei pink/red; however, these are shaped differently than spermatozoa. Yeast cells take on the same color as spermatozoa, but are shaped differently.</p> <p>7.8.3.1.2 Without phase microscopy: Spermatozoa heads are neon-like pink/red with pale pink (almost colorless) acrosomal caps, blue-green neck/midpieces, and green tails. Epithelial cells appear bright blue with red to purple nuclei.</p> <p>7.8.3.2 Document the approximate number of spermatozoa and spermatozoa heads on the smear per hpf (approximately 400X total magnification), per lpf (approximately 200X total magnification), per length of slide, or per slide, as appropriate. If only 1 spermatozoon or spermatozoon head is observed, there must be documented confirmation of its presence by a second qualified examiner.</p> <p>7.8.3.3 Place all smears submitted in the PERK back into the PERK. Properly label and return all other spermatozoa positive smears with the evidence. Note: If a stain is consumed in the preparation of a smear, properly label and return the smear even when no spermatozoa are identified.</p> <p>7.8.4 Refer to section 7.11 for reporting results.</p> <p>7.9 MICROSCOPIC EXAMINATION OF UNSTAINED SLIDES FOR SPERMATOOZOA</p> <p>7.9.1 Unstained smears may be examined using phase contrast microscopy.</p> <p>7.9.2 Equipment</p> <p>7.9.2.1 Microscope (approximately 200X – 400X total magnification) with phase capability</p>	


7 SEMEN IDENTIFICATION	Page 11 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>7.9.3 Materials</p> <p>7.9.3.1 Microscope slides</p> <p>7.9.3.2 Coverslips</p> <p>7.9.3.3 Applicator sticks</p> <p>7.9.4 Reagents</p> <p>7.9.4.1 Distilled water</p> <p>7.9.5 Procedure for the Microscopic Examination of Unstained Slides for Spermatozoa</p> <p>7.9.5.1 Place a small amount of an extract of a suspected semen stain on a microscope slide and cover with a coverslip, or add a drop of distilled water to a smear from the PERK, mix the water and the material on the smear, and cover with a coverslip.</p> <p>7.9.5.2 Scan quickly with phase at approximately 200X total magnification. Confirm with phase at approximately 400X total magnification.</p> <p>7.9.5.3 When the coverslip is touched gently, the spermatozoa and/or spermatozoa heads will roll, exhibiting their characteristic 3-dimensional shape. Use the distinctive size and morphology to identify the spermatozoa/spermatozoa heads.</p> <p>7.9.5.4 Document the approximate number of spermatozoa and spermatozoa heads on the smear per hpf (approximately 400X total magnification), per lpf (approximately 200X total magnification), per length of slide, or per slide, as appropriate. If only 1 spermatozoon or spermatozoon head is observed, there must be documented confirmation of its presence by a second qualified examiner.</p> <p>7.9.5.5 Place all smears submitted in the PERK back into the PERK. Properly label and return all other spermatozoa positive smears with the evidence. Note: If a stain is consumed in the preparation of a smear, properly label and return the smear even when no spermatozoa are identified.</p> <p>7.9.6 Refer to section 7.11 for reporting results.</p> <p>7.10 ABACard® OneStep p30 DETECTION TEST</p> <p>7.10.1 Technical Notes</p> <p>7.10.1.1 This test for the prostate protein p30 was characterized by Hochmeister, et al. in a paper entitled "Evaluation of Prostate-Specific Antigen (PSA) Membrane Test Assays for the Forensic Identification of Seminal Fluid" (reference 9, Appendix B) as an immunochromatographic PSA membrane test. A stain can be extracted for microscopic sperm identification using the sample preparation procedure set forth in section 7.11.6. If no spermatozoa are identified, the analyst can then proceed to identify semen with the OneStep ABACard® Test.</p>	

7 SEMEN IDENTIFICATION	Page 12 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p data-bbox="354 296 1118 331">7.10.1.2 Principle of the ABACard® OneStep p30 Detection Test</p> <p data-bbox="467 363 1549 867">Sample is added to the sample well “S” and if p30 is present, it will react with the mobile monoclonal antihuman p30 antibody and a mobile antibody-antigen complex is thus formed. This mobile antibody-antigen complex migrates through the absorbent device toward the test area “T”. In the test area “T”, a polyclonal antihuman p30 antibody is immobilized. This immobilized antibody captures the above complex so that an antibody-antigen-antibody sandwich is formed. The conjugated pink dye particles concentrate in a narrow zone on the membrane. When the p30 concentration in the sample exceeds 4 ng/ml the pink dye particles will form a pink colored band in the test area “T” indicating a positive test result. As an internal positive control, p30 antibody-dye conjugates cannot bind to the antibody in the test area “T”, but are captured by an immobilized anti-immunoglobulin antibody present in the control area “C” forming a complex. The captured pink dye particles will thus form a band in the control area “C” indicating that the test has worked properly and proper procedures have been followed. The presence of two colored lines, one in the test area “T” and the other in the control area “C”, indicates a positive result, while a line only in the control area “C” would indicate a negative result.</p> <p data-bbox="256 898 527 934">7.10.2 Quality Control</p> <p data-bbox="354 966 1549 1066">7.10.2.1 Before using a new lot number of the ABACard® OneStep p30 Detection Test, its specificity must be tested and appropriately documented in the laboratory’s quality control records. It is also desirable to test dilutions of semen to determine the sensitivity of the test.</p> <p data-bbox="354 1098 1549 1199">7.10.2.2 The ABACards® (“test devices”) must be tested against human blood, vaginal fluid, saliva, feces, urine, a positive control (semen), and a negative control (distilled water) to ensure that the test is semen specific.</p> <p data-bbox="467 1230 1549 1302">7.10.2.2.1 Samples of human blood, vaginal fluid, saliva, feces, urine, and semen will be prepared in-house.</p> <p data-bbox="467 1333 1549 1501">7.10.2.2.2 Label the known samples with the name of the substance (i.e., human semen, etc.) and the lot number (the date of preparation followed by the initials of the person preparing the sample). Example: human semen Lot Number 100899JD was prepared by Jane Doe on October 8, 1999.</p> <p data-bbox="467 1533 1027 1568">7.10.2.2.3 Store known samples in the freezer.</p> <p data-bbox="354 1600 1024 1635">7.10.2.3 The quality control documentation will include:</p> <p data-bbox="467 1667 1424 1738">7.10.2.3.1 The lot number, receipt date, expiration date, and manufacturer of the ABACard® OneStep p30 Detection Test.</p> <p data-bbox="467 1770 834 1806">7.10.2.3.2 The date of testing.</p> <p data-bbox="467 1837 1122 1873">7.10.2.3.3 Initials of the person conducting the testing.</p> <p data-bbox="467 1904 862 1940">7.10.2.3.4 Results of the testing.</p>	

7 SEMEN IDENTIFICATION	Page 13 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>7.10.2.4 Once the appropriate tests have been performed on a lot number of the ABACard® OneStep p30 Detection Test, they need not be repeated for each case. If another shipment of the same lot number is received on a different date, the QC testing described above must be repeated.</p> <p>7.10.3 “High Dose Hook Effect”</p> <p>7.10.3.1 The “High Dose Hook Effect” is a false negative result that is obtained in the presence of high concentrations of p30 (usually undiluted semen). This effect results from large amounts of human p30 binding to the antibody to form an antigen-antibody complex and free p30 migrating toward the test area “T”. The antibody in the test area “T” is blocked by this free p30. Therefore, the mobile antigen-antibody complex cannot bind to the antibody. As a result no pink line will form in the test area “T”. To confirm the presence of “High Dose Hook Effect”, repeat the test using a 10-10,000 fold dilution of the sample.</p> <p>7.10.4 Stability, Storage and Shelf Life</p> <p>7.10.4.1 The OneStep ABACard® p30 Detection Test should be stored below 82° F (28° C).</p> <p>7.10.4.2 The test can be stored in the sealed pouch below 82° F (28° C) until the expiration date as printed on the sealed test pouch.</p> <p>7.10.4.3 DO NOT FREEZE.</p> <p>7.10.4.4 Do not use the test after the expiration date.</p> <p>7.10.5 p30 BY OneStep ABACard® (References 9, 10, 11, Appendix B) Test Kit</p> <p>7.10.5.1 Reagents and Materials Provided</p> <p>7.10.5.1.1 Test Device (25 pieces, each individually sealed in a test pouch) - one device needed per sample tested</p> <p>7.10.5.1.2 A dropper and a desiccant sealed inside each of the test pouches</p> <p>7.10.5.1.3 Test Instructions</p> <p>7.10.5.2 Equipment Required But Not Provided</p> <p>7.10.5.2.1 Microcentrifuge</p> <p>7.10.5.2.2 Timer</p> <p>7.10.5.2.3 Scissors</p> <p>7.10.5.2.4 Tweezers</p> <p>7.10.5.2.5 Microcentrifuge tube rack</p>	

7 SEMEN IDENTIFICATION	Page 14 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>7.10.5.2.6 Pipettors (1000 µl and/or 200 µl)</p> <p>7.10.5.2.7 Dissecting needle</p> <p>7.10.5.3 Materials Required But Not Provided</p> <p>7.10.5.3.1 Microcentrifuge tubes</p> <p>7.10.5.3.2 Microcentrifuge tube lids</p> <p>7.10.5.3.3 Pipette tips</p> <p>7.10.5.4 Reagents Required But Not Provided</p> <p>7.10.5.4.1 Known semen sample</p> <p>7.10.5.4.2 Reagent blank</p> <p>7.10.5.4.3 Distilled water</p> <p>7.10.5.5 Minimum Standards and Controls</p> <p>7.10.5.5.1 On the day of use a positive reagent control (known seminal fluid) and a negative reagent control (distilled water) must be tested to ensure that the reagents and test device are working properly. The results of this testing must be documented in the case file.</p> <p>7.10.5.5.2 If either control does not give the expected result, do not proceed with testing evidence samples until the problem has been resolved as demonstrated by testing another set of positive and negative reagent controls and achieving the expected results with both controls.</p> <p>7.10.5.5.3 A substrate control (when available) must also be tested, unless the stain is on a cotton swab, and the results of the testing documented in the case file. It is not necessary to test submitted control swabs.</p> <p>7.10.5.6 p30 BY OneStep ABACard® Procedure</p> <p>7.10.5.6.1 Cut a portion of the stain into small pieces (size based upon the substrate and the intensity of the acid phosphatase test) and place into a labeled microcentrifuge tube.</p> <p>7.10.5.6.2 Add 200 µl of distilled water (250 µl if a sperm search is also being conducted) and cap the tube.</p> <p>7.10.5.6.3 Allow the sample to extract at room temperature for a minimum of 2 hours. Extraction can be done overnight if desired.</p> <p>7.10.5.6.4 Punch holes in the lid of the tube.</p>	

7 SEMEN IDENTIFICATION	Page 15 of 21								
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:								
	Effective Date: 14-March-2006								
<p>7.10.5.6.5 Place the cuttings into the lid.</p> <p>7.10.5.6.6 Centrifuge for 5 minutes at approximately 10,000 rpm to recover the liquid.</p> <p>7.10.5.6.7 If a microscopic sperm search is to be conducted, remove approximately 220 μl of the extract and place into a new labeled microcentrifuge tube. This aliquot will be used for the test procedure and may be stored between 2-8°C or frozen if not used immediately. The remaining extract and pellet can be used for the sperm search.</p> <p>7.10.5.6.8 Allow the sample to warm to room temperature if it has been refrigerated or frozen.</p> <p>7.10.5.6.9 Remove the device and dropper from the sealed pouch.</p> <p>7.10.5.6.10 Add approximately 200 μl (or 8 drops with the dropper) of the sample to the sample well “S” on a labeled test device.</p> <p>7.10.5.6.11 Record result at 10 minutes. A positive result can be seen as early as 1 minute. For negative results, one must wait for the full 10 minutes. All control samples must give the expected results before the result on an unknown sample can be called, i.e., the substrate control is negative, the reagent blank is negative, and the known semen sample is positive. A diagrammatic representation of the results is located below.</p> <p>7.10.5.6.12 Interpretation</p> <table border="0"> <tr> <td>7.10.5.6.12.1</td><td>Positive Result = 2 pink lines, one in the test area “T” and one in the control area “C” p30 level is at or above 4 ng/ml</td></tr> <tr> <td>7.10.5.6.12.2</td><td>Negative Result = 1 pink line in the control area “C” No p30 is present above 4 ng/ml <u>OR</u> presence of “High Dose Hook Effect”.</td></tr> <tr> <td>7.10.5.6.12.3</td><td>Invalid Result = No pink line in the control area “C” The test is inconclusive. Repeat the test.</td></tr> <tr> <td>7.10.5.6.12.4</td><td>Refer to the diagrammatic representation of the results on the next page.</td></tr> </table>		7.10.5.6.12.1	Positive Result = 2 pink lines, one in the test area “T” and one in the control area “C” p30 level is at or above 4 ng/ml	7.10.5.6.12.2	Negative Result = 1 pink line in the control area “C” No p30 is present above 4 ng/ml <u>OR</u> presence of “High Dose Hook Effect”.	7.10.5.6.12.3	Invalid Result = No pink line in the control area “C” The test is inconclusive. Repeat the test.	7.10.5.6.12.4	Refer to the diagrammatic representation of the results on the next page.
7.10.5.6.12.1	Positive Result = 2 pink lines, one in the test area “T” and one in the control area “C” p30 level is at or above 4 ng/ml								
7.10.5.6.12.2	Negative Result = 1 pink line in the control area “C” No p30 is present above 4 ng/ml <u>OR</u> presence of “High Dose Hook Effect”.								
7.10.5.6.12.3	Invalid Result = No pink line in the control area “C” The test is inconclusive. Repeat the test.								
7.10.5.6.12.4	Refer to the diagrammatic representation of the results on the next page.								

7 SEMEN IDENTIFICATION	Page 16 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p data-bbox="323 331 1373 363" style="text-align: center;"><i>OneStep</i> ABACard® p30 TEST RESULTS DIAGRAMMATIC REPRESENTATION</p> <div data-bbox="472 541 963 709" style="text-align: center;">  </div> <p data-bbox="149 905 1528 972">Note: OneStep ABACard® p30 Test results diagrammatic representation is taken from Abacus Diagnostics, OneStep ABACard® p30 Test For Identification of Semen, Technical Information Sheet (Revised 10/98).</p> <p data-bbox="199 1005 553 1035">7.11 REPORTING RESULTS</p> <p data-bbox="256 1073 1161 1104">7.11.1 Report the results of semen testing using the statements which follow:</p> <p data-bbox="352 1140 678 1171">7.11.1.1 Positive findings</p> <p data-bbox="467 1207 1008 1239">7.11.1.1.1 “Spermatozoa were identified ...”</p> <p data-bbox="467 1274 1036 1306">7.11.1.1.2 “A spermatozoon was identified ...”</p> <p data-bbox="467 1341 1252 1373">7.11.1.1.3 “Seminal fluid, but no spermatozoa, was identified ...”</p> <p data-bbox="352 1409 675 1440">7.11.1.2 Negative findings</p> <p data-bbox="467 1476 1203 1507">7.11.1.2.1 “No spermatozoa or seminal fluid was detected...”</p> <p data-bbox="467 1543 1490 1610">7.11.1.2.2 “No seminal fluid was detected...” (This wording will also be used if only a visual exam, with or without ALS or AP, was conducted.)</p> <p data-bbox="352 1646 719 1677">7.11.1.3 Inconclusive findings</p> <p data-bbox="467 1713 1146 1745">7.11.1.3.1 “Tests for seminal fluid were inconclusive...”</p> <p data-bbox="467 1780 1479 1848">7.11.1.3.2 “Tests for seminal fluid were inconclusive and the stain was insufficient for further body fluid identification testing...”</p>	

7 SEMEN IDENTIFICATION	Page 17 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p style="text-align: center;">SEMEN IDENTIFICATION STUDY QUESTIONS</p> <ol style="list-style-type: none"> 1. What is semen? 2. What glands contribute to seminal fluid? 3. What is p30 and where is it found? 4. What is the significance of p30 and under what circumstances would you test for it? 5. What factors can lead to a diminished sperm count in the male ejaculate? 6. Describe the mechanism and the purpose of the chemicals for the AP test. What would a positive result look like and what would a positive result tell you? 7. Describe the mechanism and the purpose of the chemicals for the p30 test. 8. Compare and contrast the different methods for detecting semen stains. 9. How does an alternate light source assist in locating stains? What alternate light sources are used by DFS (include filters used and wavelengths)? 10. What is the name of the stain used to stain smears for spermatozoa examination? What is the purpose of each chemical? 11. Describe the appearance of stained spermatozoa using phase contrast and bright field. 12. Describe the morphology of a spermatozoon. 13. What factors may affect the persistence of sperm in a living rape victim? What, if any, differences would one expect to find with regard to the persistence of sperm in a victim of rape and murder? 14. On average what is the total volume of seminal fluid per normal ejaculate? What is considered a normal sperm count per ml of seminal fluid? 15. What, if any, is the significance of observing only sperm heads versus intact sperm on a slide? 16. Explain the difference between seminal acid phosphatase and vaginal acid phosphatase. 17. If you do not detect a positive AP result on a swab or in a stain, is it possible to identify sperm? Explain your answer. 18. You get a call from an investigator saying he has a girl who is pregnant due to a rape that occurred about 6 weeks ago. She wants to get an abortion now. What do you advise the investigator? 19. How do you preserve a used condom? 20. If you have some bedding with stains and the stains test AP NEG what would be your next step? 	

7 SEMEN IDENTIFICATION	Page 18 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006
<p>21. If you have some swabs that test AP POS and an extract of the swabs is NEG for the sperm search, what is your next step?</p> <p>22. How long would you expect there to be sperm in the female reproductive tract? How about in a stain on bedding?</p> <p>23. You get a call from an investigator who says he's working a case in which a victim was raped by her husband. Her previous intercourse with him was 3 days ago. What do you advise the investigator?</p>	

7 SEMEN IDENTIFICATION	Page 19 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006

CHECKLIST FOR SEMEN IDENTIFICATION

Name of Trainee: _____

- Examination of several stained and unstained smears for spermatozoa using phase contrast microscopy. Results compared.

Date: _____ Training Coordinator: _____

Comments: _____
- Completion of presumptive and confirmatory testing on 50 (minimum) known semen samples in the following categories:

Semen of varying ages (approximately 10)

Date: _____ Training Coordinator: _____

Comments: _____

Semen on varying substrates (approximately 10)

Date: _____ Training Coordinator: _____

Comments: _____

Mixture samples (approximately 10)

Date: _____ Training Coordinator: _____

Comments: _____

Diluted Semen (neat to 1:100) (approximately 10)

Date: _____ Training Coordinator: _____

Comments: _____

Semen stains subjected to various contaminants and environmental conditions (approximately 10)

Date: _____ Training Coordinator: _____

Comments: _____
- Examination and comparison of 20 different animal semen samples.

Date: _____ Training Coordinator: _____

Comments: _____

7 SEMEN IDENTIFICATION	Page 20 of 21
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:
	Effective Date: 14-March-2006

- Completion of presumptive and confirmatory testing on various known physiological fluids, including different semen dilutions, mixtures, and aspermic semen samples (20 minimum).

Date: _____ Training Coordinator: _____

Comments: _____
- Determined the sensitivity of p30 by testing various dilutions of semen using the p30 By OneStep ABA Card® (12 minimum).

Date: _____ Training Coordinator: _____

Comments: _____
- Accurately examined a series of unknown samples of varying dilutions of seminal fluid including samples with no spermatozoa, using presumptive and confirmatory tests, as appropriate (25 minimum).

Date: _____ Training Coordinator: _____

Comments: _____
- Trainee has developed a thorough understanding of the theory behind the acid phosphatase test and the p30 By OneStep ABA Card® test for semen, including the use of controls, sources of error, and the specificity and limitations of the tests.

Date: _____ Training Coordinator: _____

Comments: _____
- Trainee has become proficient in the use of the acid phosphatase test and the p30 By OneStep ABA Card® test for semen, as well as in the extraction, staining and microscopic examination of spermatozoa.

Date: _____ Training Coordinator: _____

Comments: _____
- Trainee has become proficient in the use of alternate light sources for locating semen stains.

Date: _____ Training Coordinator: _____

Comments: _____
- Trainee has become proficient in techniques used to prevent cross-contamination between samples.

Date: _____ Training Coordinator: _____

Comments: _____

7 SEMEN IDENTIFICATION	Page 21 of 21	
TRAINING MANUAL: CASE APPROACH AND IDENTIFICATION OF BIOLOGICAL SUBSTANCES	Amendment Designator:	
	Effective Date: 14-March-2006	

11. Trainee's notebook is organized and complete.

Date: _____ Training Coordinator: _____

Comments: _____

12. Trainee has participated in a mock trial and/or practical or oral examinations. Performance was satisfactory.

Date: _____ Training Coordinator: _____

Comments: _____

13. Trainee has read and understands all applicable literature.

Date: _____ Training Coordinator: _____

Comments: _____

◆END