

10 VERTICAL TYPING GEL	Page 1 of 3
TRAINING PROGRAM FOR THE ANALYSIS OFFORENSIC CASEWORK USING PCR-BASED STR FLUORESCENCE IMAGING ANALYSIS AT THE POWERPLEX® 16 BIO LOCI	Issue No. 2
	Effective Date: 1-August-2005
<p>10 VERTICAL TYPING GEL</p> <p>10.1 GOALS:</p> <p>10.1.1 To become familiar with the theories of electrophoresis as they apply to vertical gels used in DNA STR analysis.</p> <p>10.1.2 To learn the parameters used for electrophoresis of the PowerPlex® 16 BIO typing gel.</p> <p>10.1.3 To develop an understanding and working knowledge of the use of the typing gel, including the limitations and proper documentation.</p> <p>10.1.4 To become familiar with the controls run on the PowerPlex® 16 BIO typing gel.</p> <p>10.2 TASKS:</p> <p>10.2.1 Prepare reagents and gels necessary to perform vertical gel electrophoresis of the DNA samples. Refer to the <u>Commonwealth of Virginia Department of Forensic Science Forensic Biology Section Procedure Manual, Section III - Fluorescent Detection PCR-Based STR DNA Protocol: PowerPlex® 16 BIO System</u> for the procedure.</p> <p>10.2.2 Run vertical PowerPlex® 16 BIO typing gels.</p> <p>10.2.3 Read applicable literature and become familiar with glossary terms. Refer to Appendices A, B, and C.</p> <p>10.2.4 Continue on to Chapter 11, FLUORESCENCE DETECTION.</p> <p>10.3 TRAINING EVALUATION:</p> <p>10.3.1 Knowledge</p> <p>10.3.1.1 Review of notes and worksheets in training notebook by training coordinator.</p> <p>10.3.1.2 Mini-mock trials and/or question and answer sessions.</p> <p>10.3.2 Skills</p> <p>10.3.2.1 The trainee should demonstrate an unquestionably sound technique for running consistently interpretable vertical PowerPlex® 16 BIO gels using proper documentation. This will be monitored by review of the documentation in the training notebook and continual observation by the training coordinator.</p> <p>10.3.3 Completion of the trainee checklist by the training coordinator.</p>	

10 VERTICAL TYPING GEL	Page 2 of 3
TRAINING PROGRAM FOR THE ANALYSIS OF FORENSIC CASEWORK USING PCR-BASED STR FLUORESCENCE IMAGING ANALYSIS AT THE POWERPLEX® 16 BIO LOCI	Issue No. 2
	Effective Date: 1-August-2005
<p>STUDY QUESTIONS:</p> <ol style="list-style-type: none"> 1. Which chemicals are used to cause the Page Plus typing gel to polymerize? 2. The STR loading buffer contains which dyes? 3. Why is xylene cyanol not used in the loading dye? 4. Why is Page Plus used instead of agarose? 5. What controls are used on the PowerPlex® 16 BIO gel? Why are these used? 6. Does the preparation of the glass plates affect the results? Explain. 7. What are the components of the Page Plus gel? What is the gel concentration of the Page Plus gel that is used to analyze the PowerPlex® 16 BIO System? Why is this concentration used? 	

10 VERTICAL TYPING GEL	Page 3 of 3
TRAINING PROGRAM FOR THE ANALYSIS OFFORENSIC CASEWORK USING PCR-BASED STR FLUORESCENCE IMAGING ANALYSIS AT THE POWERPLEX® 16 BIO LOCI	Issue No. 2
	Effective Date: 1-August-2005

CHECKLIST FOR VERTICAL TYPING GEL

Name of Trainee: _____

- Trainee has prepared reagents and gels necessary to perform vertical gel electrophoresis.
Date:_____ Training Coordinator:_____

Comments:_____
- Trainee has successfully and accurately completed all appropriate paperwork associated with the vertical typing gel.
Date:_____ Training Coordinator:_____

Comments:_____
- Trainee has successfully run vertical typing gels for samples amplified using the PowerPlex® 16 BIO System.
Date:_____ Training Coordinator:_____

Comments:_____
- Trainee has developed an understanding of the theory and limitations of the vertical typing gel.
Date:_____ Training Coordinator:_____

Comments:_____
- Notebook is organized and complete.
Date:_____ Training Coordinator:_____

Comments:_____
- Trainee has read and understands all applicable literature.
Date:_____ Training Coordinator:_____

Comments:_____
- Trainee has participated in mini-mock trials and/or question and answer sessions.
Date:_____ Training Coordinator:_____

Comments:_____

◆END