

<p align="center">15 GABAPENTIN QUANTITATION AND CONFIRMATION BY LCMS</p>	<p align="center">Page 1 of 4</p>
<p align="center">Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL</p>	<p>Amendment Designator:</p>
	<p>Effective Date: 31-March-2004</p>
<p align="center">15 GABAPENTIN QUANTITATION AND CONFIRMATION BY LCMS</p> <p>15.1 Summary</p> <p>15.1.1 Gabapentin is extracted from biological samples with an acetonitrile precipitation and analyzed by high performance liquid chromatography-electrospray ionization mass spectrometry (LC-ESI-MS).</p> <p>15.2 Specimen Requirements</p> <p>15.2.1 One mL blood, biological fluid or tissue homogenate.</p> <p>15.3 Reagents and Standards</p> <p>15.3.1 Ammonium acetate</p> <p>15.3.2 Methanol</p> <p>15.3.3 Acetonitrile</p> <p>15.3.4 Gabapentin (Neurontin® from Parke-Davis)</p> <p>15.3.5 Phenacetin (Merck)</p> <p>15.4 Solutions, Internal Standards, Calibrators and Controls</p> <p>15.4.1 10 mM Ammonium acetate: Weigh 0.38 g ammonium acetate. Transfer to 500 mL volumetric flask and QS to volume with dH₂O</p> <p>15.4.2 Working standard solutions for gabapentin</p> <p>15.4.2.1 10 mg/mL gabapentin stock solution. Weigh 100 mg gabapentin, transfer to 10 mL volumetric flask and QS to volume with dH₂O.</p> <p>15.4.2.2 0.1 mg/mL gabapentin working solution: Pipet 100 µl of 10 mg/mL stock solution of gabapentin into 10 mL volumetric flask and QS to volume with dH₂O</p> <p>15.4.3 Quality Control (QC) standard solutions</p> <p>15.4.3.1 0.1 mg/mL gabapentin QC solution: Pipet 100 µl of separate 10 mg/mL stock solution of gabapentin into 10 mL volumetric flask and QS to volume with dH₂O</p> <p>15.4.4 Internal standard working solution</p> <p>15.4.4.1 1 mg/mL phenacetin stock solution. Weigh 10 mg phenacetin, transfer to 10 mL volumetric flask and QS to volume with methanol.</p> <p>15.4.4.2 0.1 mg/mL phenacetin: Pipet 1 mL of 1 mg/mL phenacetin stock solution into 10 mL volumetric flask and QS to volume with dH₂O</p> <p>15.4.5 Calibrators. To prepare the calibration curve, pipet the following volumes of working solution into appropriately labeled 16 x 125 mm screw cap tubes:</p>	

15 GABAPENTIN QUANTITATION AND CONFIRMATION BY LCMS

Page 2 of 4

Division of Forensic Science

TOXICOLOGY TECHNICAL PROCEDURES MANUAL

Amendment Designator:

Effective Date: 31-March-2004

- 15.4.5.1 Cal 1: 50 mg/L gabapentin: 500 µL of 0.1 mg/mL gabapentin working solution
- 15.4.5.2 Cal 2: 20 mg/L gabapentin: 200 µL of 0.1 mg/mL gabapentin working solution
- 15.4.5.3 Cal 3: 10 mg/L gabapentin: 100 µL of 0.1 mg/mL gabapentin working solution
- 15.4.5.4 Cal 4: 5 mg/L gabapentin: 50 µL of 0.1 mg/mL gabapentin working solution
- 15.4.5.5 Cal 5: 2 mg/L gabapentin: 20 µL of 0.1 mg/mL gabapentin working solution
- 15.4.5.6 Cal 6: 1 mg/L gabapentin: 10 µL of 0.1 mg/mL gabapentin working solution
- 15.4.5.7 For each calibrator, add 1 mL blank blood to each tube to achieve final concentration.

15.4.6 Gabapentin Control (QC)

- 15.4.6.1 15 mg/L gabapentin: Pipet 150 µL of 0.1 mg/mL gabapentin QC Solution into appropriately labeled 16 x 125 mm screw cap tube and add 1 mL blank blood.
- 15.4.6.2 Negative blood control: Blood bank blood (or equivalent) previously determined not to contain gabapentin.

15.5 Apparatus

- 15.5.1 Screw cap test tubes, 16 x 125 mm
- 15.5.2 Screw cap test tubes, glass, conical bottom
- 15.5.3 Centrifuge capable of 2,000-3,000 rpm
- 15.5.4 Nitrogen evaporator with heating block
- 15.5.5 Vortex mixer
- 15.5.6 GC autosampler vials with inserts
- 15.5.7 LC/MS: Agilent Model 1100 LC-MSD
- 15.5.7.1 LCMS Instrument Conditions. The following instrument conditions may be modified or adjusted to improve separation and sensitivity.
- 15.5.7.1.1 Elution conditions:
- 15.5.7.1.1.1 Column: Agilent Hypersil BDS 125 mm X 3 mm, 3 µm particle size
- 15.5.7.1.1.2 Column thermostat: 30° C
- 15.5.7.1.1.3 Solvent A: 10 mM ammonium acetate in dH₂O
- 15.5.7.1.1.4 Solvent B: methanol
- 15.5.7.1.1.5 Isocratic elution, stop time: 6.00 min

Time	Solv. B	Flow
0.00	40	0.5

15.5.7.1.2 Spray Chamber

15 GABAPENTIN QUANTITATION AND CONFIRMATION BY LCMS

Page 3 of 4

Division of Forensic Science

TOXICOLOGY TECHNICAL PROCEDURES MANUAL

Amendment Designator:

Effective Date: 31-March-2004

- 15.5.7.1.2.1 Ionization Mode: Electrospray
- 15.5.7.1.2.2 Gas Temperature: 350° C
- 15.5.7.1.2.3 Drying Gas (N₂): 12.0 l/min
- 15.5.7.1.2.4 Nebulizer pressure: 30 psig
- 15.5.7.1.2.5 Vcap (Positive): 4000 V

15.5.7.1.3 Selected Ion Monitoring

- 15.5.7.1.3.1 Polarity: Positive
- 15.5.7.1.3.2 SIM parameters (quantitation ion)

Time	Group Name	SIM Ion	Frag-Mentor	Gain EMV	SIM Resol.	Actual Dwell
0	gabapentin	95.0	170	1.0	Low	218
		137.0	170			218
		<u>154.0</u>	170			218
		172.0	170			218
3.0	Phenacetin	110.0	170	1.0	Low	292
		138.0	170			292
		<u>180.0</u>	170			292

15.6 Procedure

- 15.6.1 Label 16 x 125 mm screw cap tubes appropriately with blank, calibrators, controls and case sample IDs.
- 15.6.2 Prepare calibrators and controls.
- 15.6.3 Add 1 mL case specimens to the appropriately labeled tubes.
- 15.6.4 Add 50 µL 0.1 mg/mL phenacetin internal standard working solution to each tube.
- 15.6.5 Slowly, add dropwise 2 mL cold (freezer temperature) acetonitrile to each tube while vortexing. Continuous vortexing, not mere mixing, is essential.
- 15.6.6 Vortex an additional 30 seconds.
- 15.6.7 Place tubes in freezer for at least 30 minutes to facilitate separation.
- 15.6.8 Centrifuge at approximately 2500 rpm for 15 minutes.
- 15.6.9 Transfer top (acetonitrile) layer to clean conical bottom tubes taking care not to transfer any lower layers.
- 15.6.10 Evaporate to dryness at approximately 50° C under nitrogen.
- 15.6.11 Reconstitute samples in 100 µL methanol. Vortex briefly. Transfer to GC autosampler vials.
- 15.6.12 Inject 1 µL of each sample on LC/MS in the API-ES/SIM Mode

15.7 Calculation

<p align="center">15 GABAPENTIN QUANTITATION AND CONFIRMATION BY LCMS</p>	<p align="center">Page 4 of 4</p>
<p align="center">Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL</p>	<p>Amendment Designator:</p>
	<p>Effective Date: 31-March-2004</p>
<p>15.7.1 Drug concentrations are calculated by linear regression analysis using the ChemStation software.</p> <p>15.8 Quality Control</p> <p>15.8.1 See Toxicology Quality Guidelines</p> <p>15.9 References</p> <p>15.9.1 C MacDonald and P Zavitsanos. Determination of Gabapentin in Human Serum by LC/MS. Agilent Application Note.</p> <p>15.9.2 D Ifa, M Falci, M Moraes, F Bezerra and G deNucci. Gabapentin Quantification in Human Plasma by LC-ESI-MS. <i>J Mass Spec</i> 36: 188-194, 2001.</p> <p>15.9.3 C Wolf, J Saady and A Poklis. Determination of Gabapentin in Serum using Solid Phase Extraction and GC/MS. <i>J Anal Tox</i> 20:498-501, 1996.</p> <p>15.9.4 J Pearson and R Steiner, in-house development.</p>	