

17 LSD QUANTITATION AND CONFIRMATION BY GCMS	Page 1 of 5
<div>Division of Forensic Science</div> <div>TOXICOLOGY TECHNICAL PROCEDURES MANUAL</div>	Amendment Designator:
	Effective Date: 31-March-2004
<div>17 LSD QUANTITATION AND CONFIRMATION BY GCMS</div> <div>17.1 Summary</div> <div>17.1.1 LSD and LSD-d<sub>3</sub> (internal standard) are extracted from biological samples using solid phase extraction (SPE) and injected into a GCMS for confirmation and quantitation by selected ion monitoring.</div> <div>17.2 Specimen Requirements</div> <div>17.2.1 2 mL of whole blood, biological fluid or tissue homogenate</div> <div>17.3 Reagents and Standards</div> <div>17.3.1 LSD, 1 mg/mL</div> <div>17.3.2 Methadone, 1 mg/mL</div> <div>17.3.3 Mepivacaine, 1 mg/mL</div> <div>17.3.4 LSD-d<sub>3</sub> 0.1 mg/mL</div> <div>17.3.5 Methanol</div> <div>17.3.6 Dichloromethane</div> <div>17.3.7 Isopropanol</div> <div>17.3.8 Toluene</div> <div>17.3.9 Hexane</div> <div>17.3.10 Isoamyl alcohol</div> <div>17.3.11 Potassium phosphate buffer solution concentrate (1 M, pH 6.0, e.g. Fisher)</div> <div>17.3.12 Ammonium hydroxide, concentrated</div> <div>17.3.13 Acetic Acid</div> <div>17.4 Solutions, Internal Standard, Calibrators and Controls</div> <div>17.4.1 Toluene:Hexane:Isoamyl Alcohol (THIA) ( 78:20:2, v:v:v) Mix 78 mL toluene, 20 mL hexane and 2 mL isoamyl alcohol</div> <div>17.4.2 Phosphate buffer 0.1 M, pH 6.0: Dilute one volume potassium phosphate buffer solution concentrate with nine volumes of dH<sub>2</sub>O.</div> <div>17.4.3 1 M Acetic acid: Add 100-200 mL dH<sub>2</sub>O to a 1 L volumetric flask. Add 57.5mL glacial acetic acid and QS to volume with dH<sub>2</sub>O.</div> <div>17.4.4 Elution Solvent (prepare fresh daily): Dichloromethane/isopropanol/ammonium hydroxide (78:20:2). Mix 78 mL dichloromethane with 20 mL isopropanol. Mix well. In hood, add 2 mL ammonium hydroxide. Mix gently.</div>	

17 LSD QUANTITATION AND CONFIRMATION BY GCMS		Page 2 of 5
Division of Forensic Science  TOXICOLOGY TECHNICAL PROCEDURES MANUAL		Amendment Designator:
		Effective Date: 31-March-2004
17.4.5	LSD calibrator stock solution 2 µg/mL: Pipet 20 µL LSD standard (1 mg/mL) into a 10 mL volumetric flask. QS to volume with methanol.	
17.4.6	LSD calibrator working solution 0.2 µg/mL: Pipet 1.0 mL LSD Stock Solution (2 µg/mL) into a 10 mL volumetric flask. QS to volume with methanol.	
17.4.7	LSD-d <sub>3</sub> 2 µg/mL internal standard solution: Pipet 200 µL LSD-d <sub>3</sub> (1mg/mL) into a 10 mL volumetric flask. QS to volume with methanol.	
17.4.8	LSD control solution (0.2 µg/mL):	
17.4.8.1	2 µg/mL LSD QC solution: Pipet 20 µL LSD stock (1 mg/mL, manufacturer or lot number different than that used for calibrators) into a 10 mL volumetric flask and QS to volume with methanol.	
17.4.8.2	0.2 µg/mL LSD QC solution: Pipet 1.0 mL 2 µg/mL LSD QC solution into 10 mL volumetric flask and QS to volume with methanol.	
17.4.9	Mepivacaine solution 25 µg/mL: Pipet 50 µL mepivacaine standard (1 mg/mL) into a 2 mL volumetric flask. QS to volume with methanol.	
17.4.10	Methadone solution 25 µg/mL: Pipet 50 µL methadone standard (1mg/mL) into a 10 mL volumetric flask. QS to volume with methanol.	
17.4.11	Calibrators. To 16 x 125 mm screw cap test tubes, add 2 mL blank blood and then add the following amounts of LSD solutions	
17.4.11.1	0.100 mg/L = 100 µL of 2 µg/mL LSD stock solution	
17.4.11.2	0.050 mg/L = 50 µL of 2 µg/mL LSD stock solution	
17.4.11.3	0.010 mg/L = 100 µL of 0.2 µg/mL LSD stock solution	
17.4.11.4	0.005 mg/L = 50 µL of 0.2 µg/mL LSD stock solution	
17.4.11.5	0.002 mg/L = 20 µL of 0.2 µg/mL LSD stock solution	
17.4.11.6	0.001 mg/L = 10 µL of 0.2 µg/mL LSD stock solution	
17.4.12	Controls	
17.4.12.1	Negative control blood: blood bank blood or equivalent determined not to contain LSD	
17.4.12.2	In house LSD control is prepared from a different lot number or different manufacturer of LSD.	
17.5 Apparatus		
17.5.1	Agilent GC/MSD, Chemstation software, compatible computer and printer	
17.5.2	Test tubes, 16 x 125 mm round bottom, screw cap tubes, borosilicate glass with Teflon caps	
17.5.3	Test tubes, 16 x 114 mm (10 mL) glass centrifuge, conical bottom	
17.5.4	Centrifuge capable of 2,000 – 3,000 rpm	
17.5.5	Cleanscreen® Extraction Cartridges (ZSDAU020) from United Chemical Technologies (200 mg columns)	

17 LSD QUANTITATION AND CONFIRMATION BY GCMS	Page 3 of 5
<div>Division of Forensic Science</div> <div>TOXICOLOGY TECHNICAL PROCEDURES MANUAL</div>	Amendment Designator:
	Effective Date: 31-March-2004
<div>17.5.6 Solid phase extraction manifold</div> <div>17.5.7 Vortex mixer</div> <div>17.5.8 Heating block</div> <div>17.5.9 Evaporator/concentrator</div> <div>17.5.10 GC autosampler vials and inserts</div> <div>17.5.11 GC/MSD Parameters. Instrument conditions may be changed or modified to improve performance and sensitivity.</div> <div>17.5.11.1 Acquisition Mode: SIM (<u>quantitation ions</u>)</div> <div>17.5.11.1.1 LSD: <u>221</u>, 323, 181</div> <div>17.5.11.1.2 LSD d<sub>3</sub>: <u>224</u></div> <div>17.5.11.1.3 Mepivacaine: <u>98</u></div> <div>17.5.11.2 Column: HP5 MS 25 m x 0.25 mm x 0.25 µm</div> <div>17.5.11.3 Detector Temperature: 280° C</div> <div>17.5.11.4 Oven Program</div> <div> <ul style="list-style-type: none"> <li>Equilibration time: 0.50 minutes</li> <li>Initial temp: 130° C</li> <li>Initial time: 1 minutes</li> <li>Ramp: 17° C/min</li> <li>Final Temp: 280° C</li> <li>Final Time: 7 minutes</li> <li>Run Time: 17 minutes</li> </ul> </div> <div>17.5.11.4.1 Inlet</div> <div> <ul style="list-style-type: none"> <li>Mode: Splitless</li> <li>Temperature: 250° C</li> <li>Injection volume: 2.0 µL</li> <li>Purge Time: ON at 2.0 minute</li> </ul> </div> <div>17.6 Procedure</div> <div>17.6.1 Label clean 16 x 125 mm screw cap tubes accordingly, blank blood (no IS), negative, calibrators, control(s) and case sample IDs.</div> <div>17.6.2 Prepare calibrators and controls.</div> <div>17.6.3 Pipet 2 mL of case samples into appropriately labeled tubes.</div> <div>17.6.4 Add 4 mL dH<sub>2</sub>O to each tube. Vortex briefly.</div> <div>17.6.5 Add 50 µL LSD-d<sub>3</sub> internal standard to each tube except tube labeled blank (no IS).</div>	

<b>17 LSD QUANTITATION AND CONFIRMATION BY GCMS</b>	Page 4 of 5
<b>Division of Forensic Science</b>  <b>TOXICOLOGY TECHNICAL PROCEDURES MANUAL</b>	Amendment Designator:
	Effective Date: 31-March-2004
<p>17.6.6 Add 20 µL of 25 µg/mL mepivacaine to each tube except tube labeled blank (no IS). Mepivacaine is an alternate internal standard used when LSD-d<sub>3</sub> internal standard contains detectable levels of LSD (assessed by comparing LSD SIM ions in blank blood with and without IS).</p> <p>17.6.7 Add 20 µL of 25 µg/mL methadone to each vial. Methadone is added as a carrier drug for the SPE extraction.</p> <p>17.6.8 Add 3.0 mL of pH 6 phosphate buffer to each tube. Vortex for 30 seconds.</p> <p>17.6.9 Centrifuge at approx 2000 rpm for 10 minutes.</p> <p>17.6.10 Condition the solid phase extraction columns. Throughout the SPE procedure, it is important not to permit the SPE sorbent bed to dry, unless specified. If necessary, add additional solvent/buffer to re-wet.</p> <p>17.6.10.1 Add 3 mL hexane to each column and aspirate on vacuum manifold</p> <p>17.6.10.2 Add 3 mL methanol to each column and aspirate on vacuum manifold.</p> <p>17.6.10.3 Add 3 mL dH<sub>2</sub>O and aspirate.</p> <p>17.6.10.4 Add 1 mL of 0.1 M pH 6.0 phosphate buffer and aspirate</p> <p>17.6.11 Without delay, pour specimens into appropriate SPE columns (leaving blood pellet formed during centrifugation in the bottom of each tube). Elute specimens from cartridges with ~ 1-2 mL/ minute flow.</p> <p>17.6.12 Wash the solid phase extraction columns:</p> <p>17.6.12.1 Add 3 mL dH<sub>2</sub>O and aspirate at ≤ 3 inches of mercury.</p> <p>17.6.12.2 Repeat the dH<sub>2</sub>O wash a second time.</p> <p>17.6.12.3 Wash with 2.0 mL 1.0 M acetic acid and aspirate.</p> <p>17.6.12.4 Wash with 1 mL methanol. Do not dry.</p> <p>17.6.12.5 Wipe the SPE column tips with Kimwipes®. Place labeled 10 mL conical test tubes in the manifold test tube rack. Be sure SPE column tips are in the designated conical tube.</p> <p>17.6.13 Elute drugs by adding 3 mL of freshly prepared dichloromethane/isopropanol/ammonium hydroxide solution to each column. Collect eluate by gravity drain (no vacuum).</p> <p>17.6.14 Cap conical tubes containing eluate. Refrigerate overnight. If an aqueous layer forms overnight, aspirate aqueous layer.</p> <p>17.6.15 Evaporate eluates to dryness at approximately 40° C under nitrogen.</p> <p>17.6.16 Reconstitute samples with 40µL THIA.</p> <p>17.6.17 Transfer to GC autosampler vials. Inject 2.0 µL on GC/MS in the SIM mode.</p>	
<p><b>17.7 Calculation</b></p> <p>17.7.1 Calculate the concentrations by interpolation of a linear plot of the response curve based on ion abundance ratios (using the target ions listed under GCMS conditions) versus calibrator concentration.</p>	

17 LSD QUANTITATION AND CONFIRMATION BY GCMS	Page 5 of 5
<div> <div>Division of Forensic Science</div> <div>TOXICOLOGY TECHNICAL PROCEDURES MANUAL</div> </div>	Amendment Designator:
	Effective Date: 31-March-2004
<div>17.8 Quality Control and Reporting</div> <div>17.8.1 See Toxicology Quality Guidelines</div> <div>17.9 References</div> <div>17.9.1 United Chemical Technologies, Inc. Clean Screen® solid phase extraction procedure for LSD in serum, plasma or whole blood.</div>	