

29 COCAINE, COCAETHYLENE AND BENZOYLECGONINE QUANTITATION AND CONFIRMATION BY LIQUID/LIQUID EXTRACTION AND GCMS	Page 1 of 6
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p style="text-align: center;">29 COCAINE, COCAETHYLENE AND BENZOYLECGONINE QUANTITATION AND CONFIRMATION BY LIQUID/LIQUID EXTRACTION AND GCMS</p> <p>29.1 Summary</p> <p>29.1.1 Cocaine (COC) and cocaethylene (CE) are extracted from biological samples by making the samples basic with saturated borate buffer and extracting with toluene/hexane/isoamyl alcohol (THIA). The remaining aqueous layer is re-extracted with a more polar organic solvent (chloroform/ethanol) to recover benzoylecgonine (BE). The BE extracts are derivatized with n-propyl iodide in dimethylsulfoxide (DMSO) in the presence of trimethylsulfonium hydroxide (TMSH). Both the COC/CE and BE extracts are analyzed by GCMS for confirmation and quantitation by selected ion monitoring.</p> <p>29.2 Specimen Requirements</p> <p>29.2.1 3 mL of whole blood, biological fluids or tissue homogenates.</p> <p>29.3 Reagents And Standards</p> <p>29.3.1 Cocaine hydrochloride</p> <p>29.3.2 Cocaethylene</p> <p>29.3.3 Benzoylecgonine</p> <p>29.3.4 Benzoylecgonine-d₃, 100 µg/mL</p> <p>29.3.5 Methapyrilene</p> <p>29.3.6 Sodium tetraborate decahydrate</p> <p>29.3.7 Toluene</p> <p>29.3.8 Hexane</p> <p>29.3.9 Isoamyl alcohol</p> <p>29.3.10 Sodium hydrogen carbonate</p> <p>29.3.11 Potassium carbonate</p> <p>29.3.12 Sulfuric Acid</p> <p>29.3.13 Dimethyl sulfoxide (DMSO)</p> <p>29.3.14 Acetonitrile</p> <p>29.3.15 Chloroform (do NOT use Burdick-Jackson chloroform with amylene preservative)</p> <p>29.3.16 Ethanol</p> <p>29.3.17 Trimethylsulfonium iodide</p>	

29 COCAINE, COCAETHYLENE AND BENZOYLECGONINE QUANTITATION AND CONFIRMATION BY LIQUID/LIQUID EXTRACTION AND GCMS	Page 2 of 6
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p>29.3.18 Silver oxide</p> <p>29.3.19 Methanol</p> <p>29.3.20 1-Iodopropane</p> <p>29.3.21 Sodium sulfite</p> <p>29.4 Solutions, Internal Standards, Calibrators, Controls</p> <p>29.4.1 Saturated borate buffer solution. Add sodium tetraborate decahydrate to dH₂O until no more dissolves after shaking vigorously.</p> <p>29.4.2 Toluene:Hexane:Isoamyl Alcohol (THIA) extraction solvent (78:20:2), v:v:v: Mix toluene (780 mL), hexane (200 mL), and isoamyl alcohol (20 mL).</p> <p>29.4.3 Sodium Hydrogen Carbonate/Potassium Carbonate (dry 3:2 w/w) Mix 300 g NaHCO₃ with 200 g K₂CO₃</p> <p>29.4.4 0.5 N Sulfuric Acid. Pipet 7 mL concentrated sulfuric acid into a 500 mL volumetric flask and QS to volume with dH₂O.</p> <p>29.4.5 Chloroform/ethanol (4:1, v/v): Mix 400 mL chloroform with 100 mL ethanol.</p> <p>29.4.6 Trimethyl sulfonium hydroxide (TMSH): Add 1.38 g trimethyl sulfonium iodide and 1.8 g silver oxide to 5 mL methanol in teflon capped foil covered tube (reacts to light). Rotate 4 hours or more. Centrifuge at approximately 2000 rpm for 10 minutes. Decant supernatant into clean Teflon capped foil covered tube. Store in freezer.</p> <p>29.4.7 Working stock solutions</p> <p>29.4.7.1 Cocaine stock solution (1 mg/mL). Weigh 11.2 mg cocaine hydrochloride into a 10 mL volumetric flask and QS to volume with acetonitrile.</p> <p>29.4.7.2 Benzoylecgonine stock solution (1 mg/mL). Weigh 10 mg benzoylecgonine into a 10 mL volumetric flask and QS to volume with methanol.</p> <p>29.4.7.3 Cocaethylene stock solution (1 mg/mL). Weigh 10 mg cocaethylene into a 10 mL volumetric flask and QS to volume with acetonitrile.</p> <p>29.4.7.4 COC/CE/BE working stock solution (0.1 mg/mL). Pipet 0.5 mL each of 1 mg/mL cocaine, cocaethylene and benzoylecgonine stock solution into a 5 mL volumetric flask and QS to volume with methanol. Prepare fresh daily.</p> <p>29.4.8 Internal Standard Solutions</p> <p>29.4.8.1 Methapyrilene stock solution (1 mg/mL). Weigh 28.5 mg methapyrilene into a 25 mL volumetric flask and QS to volume with methanol.</p> <p>29.4.8.2 Methapyrilene working stock solution (0.1 mg/mL, internal standard for COC and CE). Pipet 10 mL of 1 mg/mL methapyrilene stock solution into a 100 mL volumetric flask and QS to volume with methanol. Solution may be stored refrigerated for one year.</p> <p>29.4.8.3 BE-d₃ working solution (0.1 mg/mL, internal standard for BE). Pipet 1 mL of 1 mg/mL BE-d₃ into a 10 mL volumetric flask and QS to volume with methanol.</p>	

29 COCAINE, COCAETHYLENE AND BENZOYLECGONINE QUANTITATION AND CONFIRMATION BY LIQUID/LIQUID EXTRACTION AND GCMS	Page 3 of 6
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p>29.4.9 The following is an example of an acceptable procedure for the preparation of calibrators. Other quantitative dilutions may be acceptable to achieve similar results. To appropriately labeled 16 x 125 mm screw cap tubes, add the following volumes of the 0.1 mg/mL COC/CE/BE working solution and 3 mL blank blood to obtain the final calibrator concentrations:</p> <p>29.4.9.1 Cal 1: 0.03 mg/L: 1 µL of 0.1 mg/mL working standard + 3 mL blank blood</p> <p>29.4.9.2 Cal 2: 0.05 mg/L: 1.5 µL of 0.1 mg/mL working standard + 3 mL blank blood</p> <p>29.4.9.3 Cal 3: 0.10 mg/L: 3 µL of 0.1 mg/mL working standard + 3 mL blank blood</p> <p>29.4.9.4 Cal 4: 0.25 mg/L: 7.5 µL of 0.1 mg/mL working standard + 3 mL blank blood</p> <p>29.4.9.5 Cal 5: 0.50 mg/L: 15 µL of 0.1 mg/mL working standard + 3 mL blank blood</p> <p>29.4.9.6 Cal 6: 1.0 mg/L: 30 µL of 0.1 mg/mL working standard + 3 mL blank blood</p> <p>29.4.9.7 Cal 7: 2.0 mg/L: 60 µL of 0.1 mg/mL working standard + 3 mL blank blood</p> <p>29.4.9.8 Cal 8: 4.0 mg/L: 120 µL of 0.1 mg/mL working standard + 3 mL blank blood</p> <p>29.4.10 Controls</p> <p>29.4.10.1 Negative blood control. Blood bank blood (or comparable) determined not to contain cocaine, cocaethylene or benzoylecgonine.</p> <p>29.4.10.2 QAS Toxicology Control: 0.1 mg/L cocaine and cocaethylene and 1.0 mg/L benzoylecgonine</p> <p>29.4.10.3 In house control is prepared from a different lot number or a different manufacturer of cocaine, cocaethylene and benzoylecgonine.</p> <p>29.5 Apparatus</p> <p>29.5.1 Agilent GC/MSD, Chemstation software, compatible computer & printer</p> <p>29.5.2 Test tubes, 16 x 125 mm round bottom, screw cap tubes, borosilicate glass with Teflon caps</p> <p>29.5.3 Test tubes, 16 x 114 mm (10 mL) glass tubes, conical bottom</p> <p>29.5.4 Centrifuge capable of 2000 – 3000 rpm</p> <p>29.5.5 Vortex mixer</p> <p>29.5.6 Heating block</p> <p>29.5.7 Test tube rotator</p> <p>29.5.8 Evaporator/concentrator</p> <p>29.5.9 GC autosampler vials and inserts</p>	

29 COCAINE, COCAETHYLENE AND BENZOYLECGONINE QUANTITATION AND CONFIRMATION BY LIQUID/LIQUID EXTRACTION AND GCMS	Page 4 of 6										
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:										
	Effective Date: 31-March-2004										
<p>29.5.10 GC/MSD parameters. Instrument conditions may be changed to permit improved performance.</p> <p>29.5.10.1 Acquisition Mode: SIM</p> <p>29.5.10.2 SIM ions:</p> <table border="0"> <tr> <td>29.5.10.2.1 COC:</td> <td>303, 272, 182</td> </tr> <tr> <td>29.5.10.2.2 CE:</td> <td>317, 272, 196</td> </tr> <tr> <td>29.5.10.2.3 BE:</td> <td>210, 272, 331</td> </tr> <tr> <td>29.5.10.2.4 BE-d₃:</td> <td>213, 275, 334</td> </tr> <tr> <td>29.5.10.2.5 Methapyrilene</td> <td>261, 191, 97</td> </tr> </table> <p>29.5.10.3 Column: HP 5MS 25 m x 0.25 mm x 0.25 µm</p> <p>29.5.10.4 Detector Temperature: 280° C</p> <p>29.5.10.5 Oven Program for COC/CE</p> <ul style="list-style-type: none"> • Equilibration time: 0.50 minutes • Initial temp: 125° C • Initial time: 1.5 minutes • Ramp: 20° C/min • Final Temp: 280° C • Final Time: 3.5 minutes <p>29.5.10.6 Oven Program for BE</p> <ul style="list-style-type: none"> • Equilibration time: 0.50 minutes • Initial temp: 125° C • Initial time: 2 minutes • Ramp1: 35° C/min • Final Temp 1 : 180° C • Ramp 2: 40° C/min • Final Temp 2: 280° C • Final Time: 2 minutes <p>29.5.10.7 Inlet</p> <ul style="list-style-type: none"> • Mode: Splitless • Temperature: 250° C • Injection volume: 2.0 µL • Purge Time: ON at 2.0 minute 		29.5.10.2.1 COC:	303, 272, 182	29.5.10.2.2 CE:	317, 272, 196	29.5.10.2.3 BE:	210, 272, 331	29.5.10.2.4 BE-d ₃ :	213, 275, 334	29.5.10.2.5 Methapyrilene	261, 191, 97
29.5.10.2.1 COC:	303, 272, 182										
29.5.10.2.2 CE:	317, 272, 196										
29.5.10.2.3 BE:	210, 272, 331										
29.5.10.2.4 BE-d ₃ :	213, 275, 334										
29.5.10.2.5 Methapyrilene	261, 191, 97										
<p>29.6 Procedure</p> <p>29.6.1 Label clean 16 x 125 mm screw cap tubes accordingly, negative, calibrators, control(s) and case sample IDs.</p>											

29 COCAINE, COCAETHYLENE AND BENZOYLECGONINE QUANTITATION AND CONFIRMATION BY LIQUID/LIQUID EXTRACTION AND GCMS	Page 5 of 6
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p>29.6.2 Pipet 3 mL of blank blood, calibrators, controls and case sample bloods, fluids or tissue homogenates in appropriately labeled tubes.</p> <p>29.6.3 Add 10 µL of 0.1 mg/mL methapyrilene internal standard to all tubes and vortex briefly.</p> <p>29.6.4 Add 50 µL of 0.1 mg/mL BE-d3 internal standard to all tubes and vortex briefly.</p> <p>29.6.5 Add 3 mL saturated borate buffer to each tube.</p> <p>29.6.6 Add 6 mL THIA (78:20:2) extraction solvent to each tube.</p> <p>29.6.7 Cap and rotate tubes for 20 minutes.</p> <p>29.6.8 Centrifuge at approximately 2000 rpm for 15 minutes.</p> <p>29.6.9 Transfer upper (organic) layer to appropriately labeled 13 x 100 mm screw cap tubes and set aside for COC/CE extraction.</p> <p>29.6.10 To the lower aqueous layers (BE), add 7 mL chloroform/ethanol extraction solvent. Add solvent to no more than 4 tubes at a time. Cap immediately and invert tubes for 1 minute to prevent the formation of large clumps/clots.</p> <p>29.6.11 Rotate tubes for 10 minutes.</p> <p>29.6.12 Centrifuge at approximately 2000 rpm for 15 minutes.</p> <p>29.6.13 Aspirate and discard upper (aqueous) layer. Break protein plug. Carefully remove lower organic layer (filter if necessary) and transfer to appropriately labeled 13 x 100 mm screw cap tubes.</p> <p>29.6.14 Evaporate samples to dryness at 60-65° C under nitrogen.</p> <p>29.6.15 Add 500 µL DMSO, 100 µL 1-iodopropane and 100 µL TMSH to each tube.</p> <p>29.6.16 Cap each tube, swirl gently and heat at 60-65° C for 5 minutes.</p> <p>29.6.17 Remove tubes from heat, swirl gently and let stand 45-60 minutes for cooling.</p> <p>29.6.18 Add 0.2 mL of 0.5 N H₂SO₄ to each tube and swirl gently.</p> <p>29.6.19 Add several mg sodium sulfite to each until samples turn white.</p> <p>29.6.20 Add 2 mL THIA extraction solvent to each tube.</p> <p>29.6.21 Cap tubes and invert approximately 20-30 times.</p> <p>29.6.22 Centrifuge at approximately 2000 rpm for 10 minutes.</p> <p>29.6.23 Aspirate and discard upper (organic) layer.</p> <p>29.6.24 Adjust remaining aqueous layer to a basic pH by slowly adding solid 3:2 NaHCO₃/K₂CO₃ buffer until effervescence ceases. Then add approximately 0.3 gm excess NaHCO₃/K₂CO₃ buffer to saturate the aqueous layer.</p> <p>29.6.25 Add 400 µL THIA extraction solvent to each tube. Cap and vortex tubes for 10-15 seconds.</p>	

29 COCAINE, COCAETHYLENE AND BENZOYLECGONINE QUANTITATION AND CONFIRMATION BY LIQUID/LIQUID EXTRACTION AND GCMS	Page 6 of 6
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p>29.6.26 Centrifuge at approximately 2000 rpm for 5 minutes.</p> <p>29.6.27 Carefully transfer (upper) organic layer to appropriately labeled GC autosampler vials for BE analysis.</p> <p>29.6.28 To organic layers from 29.6.9, begin isolating COC/CE by adding 2 mL 0.5 N sulfuric acid.</p> <p>29.6.29 Cap and rotate tubes for 20 minutes.</p> <p>29.6.30 Centrifuge at approximately 2000 rpm for 15 minutes.</p> <p>29.6.31 Aspirate and discard top (organic) layer.</p> <p>29.6.32 Adjust remaining aqueous layer to a basic pH by slowly adding solid 3:2 NaHCO₃/K₂CO₃ buffer until effervescence ceases. Then add approximately 0.3 gm excess NaHCO₃/K₂CO₃ buffer to saturate the aqueous layer.</p> <p>29.6.33 Add 400 µL THIA extraction solvent to each tube. Cap and vortex tubes for 10-15 seconds.</p> <p>29.6.34 Centrifuge at approximately 2000 rpm for 5 minutes.</p> <p>29.6.35 Carefully transfer upper (organic) layer to appropriately labeled GC autosampler vials for COC/CE analysis.</p> <p>29.7 Calculations</p> <p>29.7.1 Calculate the concentrations by interpolation of a linear plot of the response curve based on peak height (or area) ratios (using the target ions listed under GCMS conditions) versus calibrator concentration.</p> <p>29.8 Quality Control</p> <p>29.8.1 See Toxicology Quality Guidelines</p> <p>29.9 References</p> <p>29.9.1 J Valentour, V Aggarwal, M McGee and S Goza. Cocaine and benzoylecgonine determinations in postmortem samples by GC. J Anal Tox 2: 134-137, 1978.</p> <p>29.9.2 V Spiehler and D Reed. Brain concentrations of cocaine and benzoylecgonine in fetal cases. J For Sci 30: 1003-1011, 1985</p>	