

<p align="center">38 VOLATILE TRACE CONGENERS IN ALCOHOLIC BEVERAGES BY GC/FID</p>	<p align="center">Page 1 of 4</p>
<p align="center">Division of Forensic Science</p> <p align="center">TOXICOLOGY TECHNICAL PROCEDURES MANUAL</p>	<p>Amendment Designator:</p>
	<p>Effective Date: 31-March-2004</p>
<p align="center">38 VOLATILE TRACE CONGENERS IN ALCOHOLIC BEVERAGES BY GC/FID</p> <p>38.1 Summary</p> <p>38.1.1 The identification and quantitation of the trace congeners is necessary to establish the identity and quality of distilled alcoholic spirits. Liquid samples are directly injected into a gas chromatograph equipped with a flame ionization detector. Using an external primary standard, the concentrations of the congeners are determined and used to help identify and characterize distilled alcoholic spirits.</p> <p>38.2 Specimen Requirements</p> <p>38.2.1 Approximately 3 mL liquid or 3 g solid material</p> <p>38.3 Reagents and Standards</p> <p>38.3.1 Ethanol, 95% (grain alcohol or equivalent)</p> <p>38.3.2 Acetaldehyde</p> <p>38.3.3 Methanol</p> <p>38.3.4 Acetone</p> <p>38.3.5 2-propanol</p> <p>38.3.6 Methyl acetate</p> <p>38.3.7 1-propanol</p> <p>38.3.8 Diacetyl</p> <p>38.3.9 Ethyl acetate</p> <p>38.3.10 2-butanol</p> <p>38.3.11 Isobutyl alcohol</p> <p>38.3.12 1-butanol</p> <p>38.3.13 Acetal</p> <p>38.3.14 Acetoin</p> <p>38.3.15 3-Methyl-1-butanol</p> <p>38.3.16 2-Methyl-1-butanol</p> <p>38.3.17 Ethyl lactate</p>	

**38 VOLATILE TRACE CONGENERS IN ALCOHOLIC
BEVERAGES BY GC/FID**

Page 2 of 4

**Division of Forensic Science
TOXICOLOGY TECHNICAL PROCEDURES MANUAL**

Amendment Designator:

Effective Date: 31-March-2004

38.4 Solutions, calibrators and controls

38.4.1 Stable Congener Standard. To a 1 L volumetric flask containing approximately 900 mL 95% ethanol, add the following volumes (or weights) of congener standard to obtain the final concentrations. Standard is stable for 2 years.

Congener	mL of standard	Final Concentration (g/100L)
1-butanol	0.3	24.30
1-propanol	0.2	16.08
2-butanol	0.2	16.20
2-methyl-1-butanol	0.5	40.80
2-propanol	0.2	15.60
3-methyl-1-butanol	1.5	121.95
acetoin	0.05 grams	5.0
acetone	0.2	15.80
isobutyl alcohol	0.5	40.10
methanol	0.5	39.60

38.4.2 Primary Congener Standard. Add approximately 9.5 mL of stable congener standard to a 10 mL volumetric flask. Add the following volumes of labile congeners and QS to volume with stable congener standard. Stable for 3-5 days.

Congener	µL of standard	Final Concentration (g/100L)
acetaldehyde	2	15.66
methyl acetate	2	18.66
ethyl acetate	7	63.07
diacetyl	2	19.62
acetyl	2	16.40
ethyl lactate	2	20.60

38.4.3 Quality control beverages

38.4.3.1 Vodka (negative control)

38.4.3.2 Bourbon whiskey quality control sample

38.4.3.3 Illegally distilled quality control sample

38.5 Apparatus

38.5.1 Agilent GC/FID, Chemstation software, compatible computer and printer

38.5.2 1.8 mL GC autosampler vials

38.5.3 GC/FID parameters. Conditions may be changed to permit improved performance.

38.5.3.1 Column: Restek RTX 1301, 60 m x 0.25 mm

38.5.3.2 Direct injection: 2 µL sample

38.5.3.3 Oven program

38 VOLATILE TRACE CONGENERS IN ALCOHOLIC BEVERAGES BY GC/FID		Page 3 of 4																																				
Division of Forensic Science		Amendment Designator:																																				
TOXICOLOGY TECHNICAL PROCEDURES MANUAL		Effective Date: 31-March-2004																																				
<table><tr><td>Equilibration Time</td><td>0.5 min</td></tr><tr><td>Initial Temp</td><td>50° C</td></tr><tr><td>Initial Time</td><td>10 min</td></tr><tr><td>Ramp 1</td><td>5° C/min</td></tr><tr><td>Final Temp 1</td><td>130° C</td></tr><tr><td>Ramp 2</td><td>40° C</td></tr><tr><td>Final Temp 2</td><td>210° C</td></tr><tr><td>Injector</td><td>150° C</td></tr><tr><td>Detector (FID)</td><td>200° C</td></tr><tr><td>Hydrogen flow</td><td>30 mL/min</td></tr><tr><td>Air flow</td><td>400 mL/min</td></tr><tr><td>Combined flow</td><td>25 mL/min</td></tr><tr><td>Make-up gas</td><td>helium</td></tr><tr><td>Inlet</td><td>Split</td></tr><tr><td>Split ratio</td><td>50:1</td></tr><tr><td>Split flow</td><td>143 mL/min</td></tr><tr><td>Total flow</td><td>148 mL/min</td></tr><tr><td>Pressure</td><td>41 psi</td></tr></table>			Equilibration Time	0.5 min	Initial Temp	50° C	Initial Time	10 min	Ramp 1	5° C/min	Final Temp 1	130° C	Ramp 2	40° C	Final Temp 2	210° C	Injector	150° C	Detector (FID)	200° C	Hydrogen flow	30 mL/min	Air flow	400 mL/min	Combined flow	25 mL/min	Make-up gas	helium	Inlet	Split	Split ratio	50:1	Split flow	143 mL/min	Total flow	148 mL/min	Pressure	41 psi
Equilibration Time	0.5 min																																					
Initial Temp	50° C																																					
Initial Time	10 min																																					
Ramp 1	5° C/min																																					
Final Temp 1	130° C																																					
Ramp 2	40° C																																					
Final Temp 2	210° C																																					
Injector	150° C																																					
Detector (FID)	200° C																																					
Hydrogen flow	30 mL/min																																					
Air flow	400 mL/min																																					
Combined flow	25 mL/min																																					
Make-up gas	helium																																					
Inlet	Split																																					
Split ratio	50:1																																					
Split flow	143 mL/min																																					
Total flow	148 mL/min																																					
Pressure	41 psi																																					
38.6 Procedure																																						
38.6.1 Transfer approximately 1.5 mL standard, controls and case samples to GC autosampler vials.																																						
38.6.2 Analyze case samples in duplicate.																																						
38.6.3 Run vodka samples (negative control) throughout the sequence to ensure no carryover of congeners.																																						
38.6.4 Inject 2 µL on GC/FID.																																						
38.7 Calculation																																						
38.7.1 Congeners are identified based on relative retention times compared to primary standard. Identification is performed by instrument software. Retention times for analytes must be within ± 2% of the retention time obtained from the primary standard.																																						
38.7.2 Concentration is calculated automatically by the software based on a single point calibration curve of the primary standard. Congener concentration is based on comparison of peak height (area) between case sample and primary standard.																																						
38.8 Quality Control and Reporting																																						
38.8.1 Positive controls. Acceptable tolerance for congener controls (Bourbon whiskey and illegally distilled quality control) is ± 20% of the target concentration (previously established mean for that control).																																						
38.8.2 Negative control. The negative control (vodka) is used to check for carryover. An acceptable negative control may not contain congeners other than ethanol.																																						
38.8.3 Duplicate tolerance. Calculate the average of the case sample duplicates. Calculate 5% of the average and a ± 5% range. Replicates must be within the ± 5% range. Reanalyze the sample if it is outside of tolerance.																																						

<p align="center">38 VOLATILE TRACE CONGENERS IN ALCOHOLIC BEVERAGES BY GC/FID</p>	<p align="center">Page 4 of 4</p>
<p align="center">Division of Forensic Science</p> <p align="center">TOXICOLOGY TECHNICAL PROCEDURES MANUAL</p>	<p>Amendment Designator:</p>
	<p>Effective Date: 31-March-2004</p>
<p>38.8.4 By comparison of congener distribution and congener ratios between the unknown sample and previously analyzed distilled spirits, a beverage class can often be determined.</p> <p>38.8.5 Congener values are never reported. The values are solely used to determine the class, brand or nature of the alcoholic beverage.</p> <p>38.9 References</p> <p>38.9.1 Amerine, M. <i>Laboratory Procedures for Enologists</i>. UC Davis, 1967.</p> <p>38.9.2 Barnett, J.H. and J.R. Einsman. <i>Occurrence and Distribution of Congeners in Distilled Alcohol Spirits</i>. J Assoc of Official Analytical Chemists Vol 60, 1977.</p> <p>38.9.3 Lange, N.A. <i>Lange's Handbook of Chemistry</i>. New York: McGraw-Hill, 1967.</p> <p>38.9.4 Lembeck, H. <i>Grossman's Guide to Wine, Beers and Spirits</i>. New York: Charles Scribner's Sons, 1983.</p> <p>38.9.5 Lichine, A. <i>Alexis Lichine's Encyclopedia of Wines and Spirits</i>. New York: Alfred Knopf, Inc., 1983.</p> <p>38.9.6 <i>Official Methods of Analysis of the Association of Official Analytical Chemists</i>. 15th Ed., 1990.</p>	