

<p align="center">42 DETERMINATION OF SPECIFIC GRAVITY AND SOLID EXTRACT OF ALCOHOLIC BEVERAGES</p>	<p align="center">Page 1 of 1</p>
<p align="center">Division of Forensic Science</p> <p align="center">TOXICOLOGY TECHNICAL PROCEDURES MANUAL</p>	<p>Amendment Designator:</p>
	<p>Effective Date: 31-March-2004</p>
<p align="center">42 DETERMINATION OF SPECIFIC GRAVITY AND SOLID EXTRACT OF ALCOHOLIC BEVERAGES</p> <p>42.1 Summary</p> <p>42.1.1 The specific gravity of alcoholic beverages is measured on a specific gravity balance. Using the specific gravity and ethanol concentration, the solid extract is calculated for the determination of the sugar content of the beverage.</p> <p>42.2 Reagents</p> <p>42.2.1 Quality control red wine</p> <p>42.3 Apparatus</p> <p>42.3.1 125 mL Erlenmeyer flasks</p> <p>42.3.2 Christian Becker Model SG-1 Specific Gravity Balance</p> <p>42.4 Procedure</p> <p>42.4.1 Fill a 125 mL Erlenmeyer flask with dH₂O. Submerge balance plummet in water and zero balance to 1.0000 (specific gravity of water).</p> <p>42.4.2 Fill 125 mL Erlenmeyer flasks with each case sample. Submerge balance plummet in the sample and read the specific gravity to nearest 0.0001. Record specific gravity on wine worksheet.</p> <p>42.4.3 Fill a 125 mL Erlenmeyer flask with quality control burgundy wine. Submerge balance plummet in the sample and read the specific gravity to nearest 0.0001. Record specific gravity on wine worksheet.</p> <p>42.5 Calculation</p> <p>42.5.1 Using the alcohol content and specific gravity of each sample, calculate solid extract from the following formula:</p> $\text{Solid extract (g/100 mL)} = \frac{\text{specific gravity of sample} - \text{specific gravity of ethanol solution}}{K}$ <p>where K = 0.00386 units of specific gravity per 1.0 g glucose per 100 mL</p> <p>42.5.2 Record solid extract to nearest 0.1 g/100 mL on wine worksheet.</p> <p>42.5.3 Specific gravity of quality control burgundy wine must agree within 10% of the target (previously established mean for that control).</p> <p>42.6 References</p> <p>42.6.1 AOAC 15th edition, 962.12, 1990.</p> <p>42.6.2 Laboratory Procedures for Enologists. Maynard Amerine, 1967, p 56.</p> <p>42.6.3 Handbook of Chemistry. Alfred Lange, 10th edition, 1967, p 1015.</p>	