

44 PROPOFOL CONFIRMATION AND QUANTITATION BY GCMS	Page 1 of 4
<div> Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL </div>	Amendment Designator:
	Effective Date: 3-August-2006
<div> <p style="text-align: center;">44 PROPOFOL CONFIRMATION AND QUANTITATION BY GCMS</p> <p>44.1 Summary</p> <p>44.1.1 Biological samples are made neutral with sodium phosphate buffer (pH 7.0) and extracted with a mixture of pentane and diethyl ether. The extract is derivatized with BSTFA in the presence of trimethylammonium hydroxide and analyzed by GCMS for confirmation and quantitation by selected ion monitoring.</p> <p>44.2 Specimen Requirements</p> <p>44.2.1 1 mL of fluid(s) or 1 g of tissue(s) or comparable amounts of fluid or tissue dilutions/homogenates</p> <p>44.3 Reagents and Standards</p> <p>44.3.1 Propofol</p> <p>44.3.2 Thymol</p> <p>44.3.3 0.1 M trimethylammonium hydroxide (TMAH)</p> <p>44.3.4 Diethyl ether</p> <p>44.3.5 Pentane</p> <p>44.3.6 Acetonitrile</p> <p>44.3.7 0.1 M disodium phosphate</p> <p>44.3.8 0.1 M monosodium phosphate</p> <p>44.3.9 Isopropyl alcohol</p> <p>44.4 Solutions, Internal Standard, Calibrators, Controls</p> <p>44.4.1 0.1 M Sodium phosphate buffer, pH 7.0: Mix 500 mL 0.1 M disodium phosphate and 250 mL 0.1 M monosodium phosphate and adjust pH as necessary.</p> <p>44.4.2 Diethyl ether: pentane (2:1, v:v): Mix 100 mL diethyl ether with 50 mL pentane.</p> <p>44.4.3 TMAH: Add 1.81 grams of trimethylammonium hydroxide pentahydrate to a 100 mL volumetric flask and qs to volum with isopropyl alcohol.</p> <p>44.4.4 Propofol Stock Solution (1.0 mg/mL): Add 10 mg of propofol to a 10 mL volumetric flask and qs to volume with acetonitrile</p> <p>44.4.5 Thymol Stock Solution (1.0 mg/mL): Add 10 mg of thymol to a 10 mL volumetric flask and qs to volume with acetonitrile.</p> <p>44.4.6 Working Standard Solution A (0.1 mg/mL): Add 0.1 ml of a 1 mg/mL propofol stock solution to a 1.0 mL volumetric flask and qs to volume with acetonitrile.</p> <p>44.4.7 Working Standard Solution B (0.01 mg/mL): Add 0.01 ml of a 1 mg/mL propofol stock solution to a 1.0 mL volumetric and qs to volume with acetonitrile.</p> </div>	

44 PROPOFOL CONFIRMATION AND QUANTITATION BY GCMS	Page 2 of 4
<div> <div>Division of Forensic Science</div> <div>TOXICOLOGY TECHNICAL PROCEDURES MANUAL</div> </div>	Amendment Designator:
	Effective Date: 3-August-2006
<div> <div>44.4.8 Working internal standard (0.02 mg/mL): Add 0.02 ml of a 1 mg/mL thymol stock solution to a 1.0 volumetric flask and qs to volume with acetonitrile.</div> <div>44.4.9 Blood calibrators, standards, and controls preparation:</div> <div> <div>44.4.9.1 To prepare the following calibration curve, pipet the following volumes of working standard solutions A or B into appropriately labeled 13 x 100 mm screw cap test tubes and then add 1 mL of whole blood.</div> <div> <div> <div>2.0 mg/L Calibrator</div> <div>20 µL of working standard solution A</div> </div> <div> <div>1.0 mg/L Calibrator</div> <div>100 µL of working standard solution B</div> </div> <div> <div>0.4 mg/L Calibrator</div> <div>40 µL of working standard solution B</div> </div> <div> <div>0.1 mg/L Calibrator</div> <div>10 µL of working standard solution B</div> </div> <div> <div>0.08 mg/L Calibrator</div> <div>8 µL of working standard solution B</div> </div> <div> <div>0.04 mg/L Calibrator</div> <div>4 µL of working standard solution B</div> </div> </div> <div>44.4.9.2 Controls</div> <div> <div>9.4.9.2.1 Negative control. Blood bank blood (or comparable) determined not to contain propofol or thymol.</div> <div>9.4.9.2.2 Positive control. In house control made from a different lot number, manufacturer, or aliquot.</div> </div> </div> <div>44.5 Apparatus</div> <div> <div>44.5.1 Agilent GC/MSD, Chemstation software, compatible computer & printer</div> <div>44.5.2 Test tubes, 13 x 100 mm round bottom, screw cap tubes, borosilicate glass with Teflon caps</div> <div>44.5.3 Centrifuge capable of 2,000 – 3,000 rpm</div> <div>44.5.4 Vortex mixer</div> <div>44.5.5 Evaporator/concentrator</div> <div>44.5.6 GC autosampler vials with inserts</div> <div>44.5.7 Test tube rotator</div> <div>44.5.8 GC/MSD parameters. Instrument conditions may be changed to permit improved performance.</div> <div>44.5.8.1 Oven program.</div> <div> <div>• Equilibration time: 0.5 minutes</div> <div>• Initial temp: 90° C</div> <div>• Initial time: 2.00 minutes</div> <div>• Ramp: 17° C/min</div> <div>• Temp: 150° C</div> <div>• Hold: 0.00</div> <div>• Ramp: 30° C/min</div> <div>• Final Temp: 290° C</div> <div>• Hold: 2.00 minutes</div> <div>• Run Time: 12.2 minutes</div> </div> </div> </div>	

44 PROPOFOL CONFIRMATION AND QUANTITATION BY GCMS	Page 3 of 4
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 3-August-2006
<p>44.5.8.2 Inlet.</p> <ul style="list-style-type: none"> • Mode: Pulsed Splitless • Temperature: 250° C • Constant pressure: 11.83 psi • Purge flow: 50.0 mL/min • Total flow: 54.1 mL/min • Injection volume: 2.0 µL <p>44.5.8.3 Column: HP-5 30 m x 0.25 mm x 0.25 µm</p> <p>44.5.8.4 Detector Temperature: 280° C</p> <p>44.5.8.5 Acquisition Mode: SIM</p> <p>44.5.8.6 SIM ions:</p> <p style="margin-left: 40px;">propofol: 235, 236, 250</p> <p style="margin-left: 40px;">thymol: 207, 222</p>	
44.6 Procedure	
<p>44.6.1 Label clean 13 x 100 mm screw cap tubes accordingly, negative, calibrators, control(s) and case sample IDs.</p> <p>44.6.2 Pipet 1 mL of blank blood, calibrators, controls and case sample bloods, fluids or tissue homogenates in appropriately labeled tubes.</p> <p>44.6.3 Add 20 µL of internal standard (thymol) working solution to each tube and vortex.</p> <p>44.6.4 Add 1 mL 0.1 M sodium phosphate buffer (pH 7.0) to each tube.</p> <p>44.6.5 Add 3 mL extraction solvent (diethyl ether:pentane) to each tube.</p> <p>44.6.6 Cap and rotate tubes for 30 minutes.</p> <p>44.6.7 Centrifuge at approx 2800 rpm for 15 minutes. Transfer organic (upper) layer to clean 5 mL conical bottom tubes.</p> <p>44.6.8 Add 20 µL TMAH to each tube and evaporate to dryness under a stream of nitrogen at room temperature.</p> <p>44.6.9 Add 100µL BSTFA to each tube, cap, vortex briefly, and heat at 85° C for 15 minutes.</p> <p>44.6.10 Transfer small aliquot to appropriately labeled GC vials and inject 2 µl on GC-MSD.</p>	
44.7 Calculation	
<p>44.7.1 Calculate the concentrations by interpolation of a linear plot of the response curve based on peak height (or area) ratios versus calibrator concentration.</p>	
44.8 Quality Control And Reporting	
<p>44.8.1 See Toxicology Quality Guidelines</p>	

44 PROPOFOL CONFIRMATION AND QUANTITATION BY GCMS	Page 4 of 4
<div> <div>Division of Forensic Science</div> <div>TOXICOLOGY TECHNICAL PROCEDURES MANUAL</div> </div>	Amendment Designator:
	Effective Date: 3-August-2006
<div> <div>44.9 References</div> <div> <div>44.9.1 Stetson, P.L., Domino, E.F. and Sneyd, J.R. Determination of plasma propofol levels using gas chromatography-mass spectrometry with selected-ion monitoring. <i>Journal of Chromatography</i>. 620 (1993) pp. 260-267.</div> <div>44.9.2 Dwight Flammia, PhD. and Henry Bateman in house development.</div> </div> <div>◆End</div> </div>	