

<p align="center">5 SMALL PARTICLE REAGENT</p>	<p align="center">Page 1 of 2</p>
<p align="center">Division of Forensic Science</p> <p align="center">LATENT FINGERPRINTS PROCEDURES MANUAL</p>	<p align="center">Amendment Designator:</p>
	<p align="center">Effective Date: 29-January-2004</p>
<div data-bbox="646 291 1055 321" data-label="Section-Header"> <h2 align="center">5 SMALL PARTICLE REAGENT</h2> </div> <div data-bbox="151 354 422 384" data-label="Section-Header"> <h3>5.1 INTRODUCTION</h3> </div> <div data-bbox="209 415 1528 569" data-label="Text"> <p>Small particle reagent (SPR) was devised and refined by the British Home Office as an effective procedure for processing wet surfaces. Surfaces, both porous and nonporous, which are wet at the time of latent print deposit or become wet after deposit, seldom retain sufficient water soluble material for conventional processing methods. Nonporous items which have been allowed to dry offer some potential if the deposit contains non-water soluble oily matter, but the drying process lessens the possibility of adequate adhesion for powders.</p> </div> <div data-bbox="209 598 1528 720" data-label="Text"> <p>Molybdenum disulfide is a lipid-sensitive reagent. Initial efforts to create a suspension of molybdenum disulfide in water used photoflo as a means of reducing surface tension. These met with limited success. Introduction of photoflo to the mixture requires a critical measurement as too much photoflo prevents complete adhesion of the molybdenum disulfide particles to the lipids. Organic solvents can not be used as these solvents may remove the lipid material.</p> </div> <div data-bbox="209 749 1528 903" data-label="Text"> <p>Refinements in the surfactant solution have not only improved the uniformity of suspension but have increased the application of SPR to other surfaces. SPR is very effective in the secondary treatment of cyanoacrylate ester developed impressions by adhering to faint impressions generally better than powders. Molybdenum disulfide is produced in various particle sizes. Smaller particle size is more effective: Lightning Powder Company provides the proper particle size. If distilled water is not available deionized water may be used.</p> </div> <div data-bbox="151 932 425 961" data-label="Section-Header"> <h3>5.2 PREPARATIONS</h3> </div> <div data-bbox="209 993 568 1022" data-label="Section-Header"> <h4>5.2.1 Surfactant Stock Solution</h4> </div> <div data-bbox="297 1054 1528 1115" data-label="List-Group"> <ol style="list-style-type: none"> 1. Dissolve 8 milliliters of Tergitol 7 in 500 milliliters of distilled water. This will make approximately 10 liters of working solution. </div> <div data-bbox="209 1146 670 1176" data-label="Section-Header"> <h4>5.2.2 SPR Suspension-Working Solution</h4> </div> <div data-bbox="297 1207 1528 1390" data-label="List-Group"> <ol style="list-style-type: none"> 1. Add 10 grams of molybdenum disulfide to 50 milliliters of the surfactant stock solution. Add the molybdenum disulfide slowly and stir continuously. 2. The mixture should be a creamy consistency free of any dry powder. 3. While stirring continuously, add the mixture to 900 milliliters of distilled water. </div> <div data-bbox="151 1421 734 1451" data-label="Section-Header"> <h3>5.3 MINIMUM STANDARDS AND CONTROLS</h3> </div> <div data-bbox="209 1482 1528 1543" data-label="Text"> <p>Molybdenum disulfide works by adhering to latent print residue. Due to the inherent ability to adhere and discolor these materials, there is no need for test impressions to be done prior to evidence application.</p> </div> <div data-bbox="151 1575 565 1604" data-label="Section-Header"> <h3>5.4 PROCEDURE OR ANALYSIS</h3> </div> <div data-bbox="209 1635 529 1665" data-label="Section-Header"> <h4>5.4.1 Immersion Technique</h4> </div> <div data-bbox="297 1696 1528 1938" data-label="List-Group"> <ol style="list-style-type: none"> 1. Shake the working solution well and place in a shallow tray such as a photographic tray. The tray should be filled until it will cover the item to be processed. 2. Stir the solution again and before each item is placed into the solution. 3. Place the item to be processed in the liquid to lie as flat as possible in the tray. 4. Allow the item to remain in the suspension and the molybdenum particles to settle on the item for 30 seconds. </div>	

<p align="center">5 SMALL PARTICLE REAGENT</p>	<p align="center">Page 2 of 2</p>
<p align="center">Division of Forensic Science</p> <p align="center">LATENT FINGERPRINTS PROCEDURES MANUAL</p>	<p align="center">Amendment Designator:</p>
	<p align="center">Effective Date: 29-January-2004</p>
<div data-bbox="297 264 1461 569"> <ol style="list-style-type: none"> 5. The item is then turned over and again allowed to set for 30 seconds. 6. This procedure is continued until all surfaces have been exposed to the solution. 7. The item is then placed into a tray of clear tap water. The tray can be rocked or a flow of tap water can be established in the tray. The excess SPR will readily be removed. 8. The item is allowed to dry. 9. All impressions should be photographed and can subsequently be lifted. </div> <div data-bbox="207 600 560 630"> <p>5.4.2 Wash Bottle Application</p> </div> <div data-bbox="297 661 1026 751"> <ol style="list-style-type: none"> 1. Spray a flow of SPR over the surface of the item. 2. Wash the surface with a light to moderate flow of clear tap water. </div> <div data-bbox="207 783 1549 873"> <p>Larger items may be processed using a wash bottle to spray a flow of SPR over the surface. For outdoor application of very large items, such as a wet automobile, a garden sprayer can be used. Generally light to moderate flows of rinse water will not dislodge the molybdenum disulfide particles.</p> </div> <div data-bbox="151 905 621 934"> <p>5.5 INTERPRETATION OF RESULTS</p> </div> <div data-bbox="207 966 1500 1087"> <p>SPR lifts easily from dried, processed, nonporous surfaces but all developed impressions should be photographed prior to lifting. Faint impressions may benefit from a reprocessing of the item. The intense black color generally facilitates photographic preservation. When SPR is used as a secondary technique after cyanoacrylate ester fuming, the results are sometimes superior to powders in both adhesion and clarity of detail.</p> </div> <div data-bbox="151 1119 388 1148"> <p>5.6 REFERENCES</p> </div> <div data-bbox="207 1180 1484 1392"> <ol style="list-style-type: none"> 1. Lee, Henry C.; Gaensslen, R. E., eds. <i>Advances in Fingerprint Technology</i>; Elsevier Science Publishers, NY, 1991. 2. Onstwedder, John III; Thomas E. Gamboe. "Small Particle Reagent: Developing Latent Prints on Water-Soaked Firearms and Effect on Firearms Analysis"; <i>Journal of Forensic Sciences</i>, 1989, 34, 2, 321-327. 3. Pounds, C.A.; R.J. Jones. "Physicochemical Techniques in the Development of Latent Fingerprints"; <i>Trends in Analytical Chemistry</i>, 1983, 2, 8, 180-183. </div> <div data-bbox="1490 1423 1549 1453"> <p align="right">◆End</p> </div>	