

<p align="center">8 COLOR TESTS</p>	<p align="center">Page 1 of 8</p>
<p align="center">Division of Forensic Science</p> <p align="center">CONTROLLED SUBSTANCES PROCEDURES MANUAL</p>	<p>Amendment Designator: A</p>
	<p>Effective Date: 15-June-2005</p>
<p align="center">8 COLOR TESTS</p> <p>8.1 Introduction:</p> <p>8.1.1 Color tests are used as a screening test at the beginning of an analysis. Most are performed on clean porcelain or disposable plastic spot plates; however some may be performed in disposable culture tubes (e.g., Scott's, Tannic Acid).</p> <p>8.1.2 The test reagent should be added to the plate or tube first, and then the questioned sample. This practice determines if the plate or tube was clean before the analysis.</p> <p>8.1.3 Several of the listed reagents have more than one recipe listed. Any of the listed, referenced recipes may be utilized in casework and should be reflected in the reagent logbook.</p> <p>8.1.4 Reagents, indicators and solutions listed in the USP-NF may be used for their published purposes.</p> <p>8.1.5 Positive color reactions are noted in each of the individual drug sections. These positive reactions may be recorded in the analytical notes by the use of a plus (+), a plus circled (⊕) or an abbreviation (e.g., pos). The color observed must be noted for drugs not routinely encountered. Negative reactions may be recorded in a similar fashion. (For Duquenois-Levine results see Section 6.5.1)</p> <p>8.1.6 Most color test reagents are comprised of strong acids and chemicals requiring careful handling. Appropriate safety precautions should be observed. Refer to MSDS's for storage and handling.</p> <p>8.1.7 Reagents should be made in quantities to minimize waste. The shelf life of color test reagents is two years.</p> <p>8.2 Color Tests and Reagents:</p> <p>The following lists the commonly used color test reagents and some examples of reactions with various drugs. The references for each test are in parenthesis.</p> <p>8.2.1 <u>Bates Test</u> (8.3.5) tests for cocaine base.</p> <p>8.2.1.1 Procedure: The Bates test is used as the second part of the Cobalt thiocyanate test. If the Cobalt thiocyanate test is negative, add Marquis reagent to spot well.</p> <p>8.2.1.2 Results: The formation of a very blue precipitate indicates cocaine base, other compounds give weaker blue or no reaction.</p> <p>8.2.2 <u>Benedict's Solution</u> (8.3.2) tests for reducing sugars and some antibiotics.</p> <p>8.2.2.1 Recipe: 1.73 g of copper sulfate in 10 mL of water. With the aid of heat, dissolve 17.3 g trisodium citrate and 10 g of anhydrous sodium carbonate in 80 mL of H₂O. Pour the two solutions together and let cool. Dilute to 100 mL with water.</p> <p>8.2.2.2 Procedure: Add 0.5 mL of the reagent to sample and heat.</p> <p>8.2.2.3 Results:</p> <ul style="list-style-type: none"> • Ascorbic acid, strong reducing agents, glucose, tetracycline – red • Streptomycin - orange/brown <p>8.2.3 <u>Chen's Test</u> (8.3.2) tests for phenethylamines.</p>	

8 COLOR TESTS	Page 2 of 8
Division of Forensic Science CONTROLLED SUBSTANCES PROCEDURES MANUAL	Amendment Designator: A
	Effective Date: 15-June-2005
<p>8.2.3.1 Recipe: 1 g copper sulfate and 1 mL glacial acetic acid in 100 mL H₂O.</p> <p>8.2.3.2 Procedure: Make an approximate 1% aqueous solution of the sample, add equal volumes of Chen's reagent and 2N NaOH.</p> <p>8.2.3.3 Results: ephedrine, PPA and pseudoephedrine – purple</p> <p>8.2.4 <u>Cobalt Thiocyanate</u> reacts with tertiary and quaternary amines to form a blue precipitate and used for general screening. May be used in conjunction with the Bates test (8.2.1) or the Stannous Chloride test (8.2.23).</p> <p>8.2.4.1 Recipes:</p> <ul style="list-style-type: none"> • 2 g cobalt thiocyanate in 100 mL H₂O or methanol (8.3.1) • 2 g cobalt thiocyanate in 100 mL H₂O and 100 mL of glycerine. (8.3.3) • 1.4 g CoCl₂ · 6H₂O and 0.9 g NH₄SCN in 100 mL H₂O. (8.3.7) <p>8.2.4.2 Procedure: Place reagent in well and add sample.</p> <p>8.2.4.3 Results:</p> <ul style="list-style-type: none"> • Cocaine HCl – blue precipitate forms, cocaine base may be initially negative or faintly blue, but blue intensifies upon the addition of dilute HCl. • PCP - blue • Amitriptyline / doxepin - blue • barbiturates with unsaturated side chain (i.e. butalbital) - faint blue <p>8.2.5 <u>Dille - Koppanyi Test</u> (8.3.9) reacts with barbiturates.</p> <p>8.2.5.1 Recipe:</p> <ul style="list-style-type: none"> • DK1: 0.1 g cobaltous acetate tetrahydrate in 100 mL methanol plus 0.2 mL glacial acetic acid • DK2: 5 mL isopropyl amine in 95 mL methanol <p>8.2.5.2 Procedure: This is a two part test. Place 2 drops of DK1 reagent in a well. Add sample. Add 1 drop of DK2 reagent. When doing multiple samples, they should be separated to avoid cross-contamination due to reagent spreading.</p> <p>8.2.5.3 Results:</p> <ul style="list-style-type: none"> • barbiturates - blue purple • theophylline, glutethimide and hydantoins - purple • ampicillin - brown <p>8.2.6 <u>Duquenois - Levine Test</u> (8.3.3, 8.3.4) reacts with marijuana and hash oil.</p> <p>8.2.6.1 Recipe: 4 g vanillin and 2.5 mL fresh acetaldehyde per 200 mL 95% ethanol</p> <p>8.2.6.2 Procedure: See Marijuana section.</p> <p>8.2.6.3 Results: marijuana/hash oil – blue/purple, pink/purple extracts into CHCl₃</p> <p>8.2.7 <u>Ehrlich's Reagent</u> (8.3.7) reacts with indole moiety and some amines.</p> <p>8.2.7.1 Recipe: 5 g p-dimethylaminobenzaldehyde to 50 mL of 95% ethanol and 50 mL of conc. HCl</p>	

8 COLOR TESTS	Page 3 of 8
Division of Forensic Science CONTROLLED SUBSTANCES PROCEDURES MANUAL	Amendment Designator: A
	Effective Date: 15-June-2005
<p>8.2.7.2 Procedure: Place reagent in well and add sample.</p> <p>8.2.7.3 Results:</p> <ul style="list-style-type: none"> • LSD, psilocyn - purple (beware of leaching of dyes in blotter paper or tablets) • benzocaine, procaine - yellow <p>8.2.8 <u>Fehlings Solution</u> (8.3.8) reacts with reducing compounds such as sugars.</p> <p>8.2.8.1 Recipe:</p> <ul style="list-style-type: none"> • Fehlings1 - 3.46 g copper sulfate per 50 mL H₂O • Fehlings2 - 86.5 g sodium potassium tartrate and 35 g of NaOH per 250 mL of H₂O <p>8.2.8.2 Procedure: Dissolve sample in water and mix. Add 5 drops of Fehlings1 and 5 drops of Fehlings2 and mix. Heat on steambath for approximately 5 minutes or until warm.</p> <p>8.2.8.3 Results: reducing sugars – yellow to red.</p> <p>8.2.9 <u>Ferric Chloride</u> (FeCl₃) tests for phenols and GHB.</p> <p>8.2.9.1 Recipes:</p> <ul style="list-style-type: none"> • 9% aqueous solution (8.3.14) • 5% aqueous solution (8.3.2) <p>8.2.9.2 Procedure: Place sample into a solution of water or methanol and add a drop of reagent.</p> <p>8.2.9.3 Results:</p> <ul style="list-style-type: none"> • salicylamide - dark purple • acetaminophen - blue • hydrolyzed aspirin – purple (to hydrolyze a sample, place in H₂O, add a little acid and heat) • GHB – red/brown <p>8.2.10 <u>Fiegel's / Nitroprusside (nitroferricyanide)</u> (8.3.3, 8.3.8) for secondary amines.</p> <p>8.2.10.1 Recipe: 1 g of sodium nitroprusside in 100 mL H₂O and 10 mL acetaldehyde</p> <p>8.2.10.2 Procedure: Dissolve sample in 2N Na₂CO₃ and add reagent.</p> <p>8.2.10.3 Storage: store in brown bottle and refrigerate.</p> <p>8.2.10.4 Results: secondary amines - deep blue color</p> <p>8.2.11 <u>Froehdes</u> (8.3.1, 8.3.2) reacts with narcotics and used for general screening.</p> <p>8.2.11.1 Recipe: 0.5 g ammonium molybdate per 100 mL H₂SO₄ (conc.)</p> <p>8.2.11.2 Procedure: Place reagent in well and add sample.</p> <p>8.2.11.3 Results:</p>	

8 COLOR TESTS	Page 4 of 8
Division of Forensic Science CONTROLLED SUBSTANCES PROCEDURES MANUAL	Amendment Designator: A
	Effective Date: 15-June-2005
<div data-bbox="428 268 799 491"> <ul style="list-style-type: none"> • heroin - purple → green • codeine - green → red/brown • morphine - deep purple → slate • aspirin - blue → purple • phenoxymethylpenicillin - blue • pentazocine – blue • acetaminophen – pale blue </div> <div data-bbox="219 520 1549 579"> <p>8.2.12 <u>GHB Color Test #3 (Smith Test)</u> (8.3.13) for GHB powders and solutions. This test will not react with GBL or 1,4-butanediol.</p> </div> <div data-bbox="315 613 509 642"> <p>8.2.12.1 Recipe:</p> </div> <div data-bbox="428 676 1549 957"> <ul style="list-style-type: none"> • Bromocresol Green – 0.03 g bromocresol green in 100 mL of 4:1 methanol:DI water adjusting the pH to 7.0 with 0.1 sodium hydroxide • Methyl Orange – 0.01 g methyl orange in 100 mL DI water adjusting the pH to 7.0 with 0.1 N sodium hydroxide • Modified Schweppes Reagent: Mix solutions A and B <ul style="list-style-type: none"> ○ Solution A – 2 g dextrose in 20 mL of DI water ○ Solution B – 2.4 g aniline hydrochloride in 20mL methanol • Mix Bromocresol Green and Methyl Orange together in a 1:1 ratio. Then, mix that indicator mixture with the modified Schweppes reagent in a 3:1 ratio. </div> <div data-bbox="315 987 1520 1045"> <p>8.2.12.2 Procedure: Add 0.5 mL of a liquid sample or a small amount of powder to a test tube. Add 2 drops of the mixed reagent and gently swirl.</p> </div> <div data-bbox="315 1079 514 1108"> <p>8.2.12.3 Results:</p> </div> <div data-bbox="428 1142 1549 1201"> <ul style="list-style-type: none"> • GHB – immediate green color • Negative results – pinkish orange (generally the same or slightly darker than the original test solution) </div> <div data-bbox="219 1234 1343 1264"> <p>8.2.13 <u>Marquis</u> (8.3.1, 8.3.2) reacts with opiates and phenethylamines and is used for general screening.</p> </div> <div data-bbox="315 1297 509 1327"> <p>8.2.13.1 Recipe:</p> </div> <div data-bbox="428 1360 1002 1419"> <ul style="list-style-type: none"> • 10 mL 37% formaldehyde in 90 mL H₂SO₄ (conc.) • 2 mL 37% formaldehyde in 75 mL H₂SO₄ (conc.) </div> <div data-bbox="315 1453 956 1482"> <p>8.2.13.2 Procedure: Place reagent in well and add sample.</p> </div> <div data-bbox="315 1516 743 1545"> <p>8.2.13.3 Storage: Keep tightly capped.</p> </div> <div data-bbox="315 1579 514 1608"> <p>8.2.13.4 Results:</p> </div> <div data-bbox="428 1642 979 1793"> <ul style="list-style-type: none"> • opiates - purple • amphetamine/methamphetamine - orange/brown • aspirin – pink → deep red on standing • phenoxymethylpenicillin - red • MDA/MDMA - black </div>	

8 COLOR TESTS	Page 5 of 8
Division of Forensic Science CONTROLLED SUBSTANCES PROCEDURES MANUAL	Amendment Designator: A
	Effective Date: 15-June-2005
<p>8.2.14 <u>Mayer's Reagent</u> (8.3.7) reacts with alkaloids</p> <p>8.2.14.1 Recipe: Dissolve 1 g of mercuric chloride in 100 mL H₂O. Add solid potassium iodide until the reddish precipitate first formed dissolves. Reagent should be clear and pale yellow in color.</p> <p>8.2.14.2 Procedure: Add 0.1 N HCl to a test tube. Add sample to acid and mix. Add Mayer's reagent to the acid solution.</p> <p>8.2.14.3 Results: alkaloids – a white to yellow precipitate is formed</p> <p>8.2.15 <u>Meckes</u> (8.3.1, 8.3.2) reacts with narcotics and used for general screening.</p> <p>8.2.15.1 Recipe: 1 g selenious acid per 100 mL H₂SO₄ (conc.)</p> <p>8.2.15.2 Procedure: Add reagent to well then sample.</p> <p>8.2.15.3 Results:</p> <ul style="list-style-type: none"> • heroin - green/blue • codeine - bright-green/blue green • PCP - light yellow • quinine - light yellow <p>8.2.16 <u>Methylene Blue</u> (8.3.10) reacts with vitamin C.</p> <p>8.2.16.1 Recipe: 12.5 mg of methylene blue dissolved in 25 mL of 95% ethanol.</p> <p>8.2.16.2 Procedure: Add reagent to well and add sample. It may be helpful to run a blank to compare the results.</p> <p>8.2.16.3 Results: Vitamin C - slowly decolorizes solution from dark blue to light blue.</p> <p>8.2.17 <u>Nitric Acid</u> (HNO₃) (8.3.1, 8.3.2) reacts with opiates and mescaline.</p> <p>8.2.17.1 Recipe: concentrated nitric acid</p> <p>8.2.17.2 Procedure: Place reagent in well and add sample.</p> <p>8.2.17.3 Results:</p> <ul style="list-style-type: none"> • heroin - yellow green • morphine - red • codeine - orange • mescaline - red • acetaminophen – fumes, orange brown <p>8.2.18 <u>Parri</u> (8.3.11) reacts with barbiturates.</p> <p>8.2.18.1 Recipe: cobaltous acetate (solid), barium oxide (solid), and methanol</p> <p>8.2.18.2 Procedure: Mix cobaltous acetate, BaO and powdered sample in equal parts in a spot plate well, add methanol.</p> <p>8.2.18.3 Results: barbiturates - blue</p>	

8 COLOR TESTS	Page 6 of 8
Division of Forensic Science CONTROLLED SUBSTANCES PROCEDURES MANUAL	Amendment Designator: A
	Effective Date: 15-June-2005
<p>8.2.19 <u>Scotts - Modified Cobalt Thiocyanate</u> (8.3.3) reacts with cocaine.</p> <p>8.2.19.1 Recipe: 2 g cobalt thiocyanate in 100 mL H₂O and 100 mL of glycerine.</p> <p>8.2.19.2 Procedure: Add reagent to well or tube and add sample. Dissolve the blue precipitate from the Co(SCN)₂ by the addition of HCl. Add CHCl₃</p> <p>8.2.19.3 Results: cocaine - blue color in the lower (CHCl₃) layer.</p> <p>8.2.20 <u>Silver Nitrate</u> (8.3.15) indicates the presence of chloride ions.</p> <p>8.2.20.1 Recipe: 5.0% w/v solution of silver nitrate in DI water.</p> <p>8.2.20.2 Caution: Poison; will cause staining.</p> <p>8.2.20.3 Storage: Store in the refrigerator in a dark environment.</p> <p>8.2.20.4 Procedure: Dissolve sample in water. Add silver nitrate solution. A white, curdy precipitate will form in the presence of chloride ions which will be insoluble in nitric acid. The precipitate will be soluble in 6N ammonium hydroxide.</p> <p>8.2.21 <u>Stannous Chloride modification for Co(SCN)₂ - HCl acidified</u> (8.3.9) differentiates between "caines".</p> <p>8.2.21.1 Recipe: 5 g SnCl₂ and 10 mL conc. HCl diluted to 100 mL with H₂O</p> <p>8.2.21.2 Procedure: The Stannous Chloride test is used as the second part of the Cobalt thiocyanate test. After performing the cobalt thiocyanate test, add a drop of stannous chloride reagent.</p> <p>8.2.21.3 Results:</p> <ul style="list-style-type: none"> • Cocaine salts - blue remains • Cocaine base - blue color forms (initially negative) • Other compounds which turned blue initially - blue fades <p>8.2.22 <u>Sulfuric Acid</u> (H₂SO₄) (8.3.1, 8.7.2)</p> <p>8.2.22.1 Recipe: concentrated sulfuric acid</p> <p>8.2.22.2 Procedure: Add reagent to well and add sample.</p> <p>8.2.22.3 Results:</p> <ul style="list-style-type: none"> • tetracycline - purple turning to yellow upon addition of water • 2,3-MDMA, 2,3-MDA – rose • 3,4-MDMA, 3,4-MDA – gray-brown <p>8.2.23 <u>Tannic Acid</u> (8.3.3) reacts with xanthines.</p> <p>8.2.23.1 Recipe: 1% aqueous solution of tannic acid</p> <p>8.2.23.2 Procedure: Add reagent to test tube then add powdered sample.</p> <p>8.2.23.3 Results: caffeine and theophylline - positive test will produce a precipitate which develops from "streamers" immediately visible in the solution.</p>	

8 COLOR TESTS	Page 7 of 8
Division of Forensic Science CONTROLLED SUBSTANCES PROCEDURES MANUAL	Amendment Designator: A
	Effective Date: 15-June-2005
<p>8.2.24 <u>TBPEE Solution</u> (8.3.8) differentiates between amines.</p> <p>8.2.24.1 Recipe:</p> <ul style="list-style-type: none"> • 0.01g Tetrabromophenolphthalein ethyl ester (TBPEE) in 100 mL CCl₄ • 10.6 g sodium carbonate in 100 mL H₂O (2N solution) <p>8.2.24.2 Caution: Carbon tetrachloride is carcinogenic. Use appropriate safety precautions.</p> <p>8.2.24.3 Procedure: Dissolve suspected amine in 2N Na₂CO₃ solution and add TBPEE solution. Note color change in the bottom TBPEE layer.</p> <p>8.2.24.4 Results:</p> <ul style="list-style-type: none"> • primary amines - violet • secondary amine - blue • tertiary amine – red <p>8.2.25 <u>Van Urk's</u> (8.3.12) reacts with the indole moiety and some amines.</p> <p>8.2.25.1 Recipe: 125 mg p-dimethylaminobenzaldehyde, 65 mL of concentrated H₂SO₄, and 2 drops of ferric chloride (USP T. S.) diluted to 100 mL with distilled water.</p> <p>8.2.25.2 Procedure: Add reagent to well then add sample.</p> <p>8.2.25.3 Results: LSD – blue/purple</p> <p>8.2.26 <u>Weber Test</u> (8.3.6) reacts with psilocyn.</p> <p>8.2.26.1 Recipe: 0.01 g of Fast Blue B or Fast Blue BB in 10 mL H₂O</p> <p>8.2.26.2 Procedure: Add 2 to 3 drops of reagent to a sample of mushrooms. Observe slow color change. Add 1 to 2 drops of conc. HCl, observe color change.</p> <p>8.2.26.3 Results: Psilocyn – Initially the solution turns red. The solution will turn from red to blue when the acid is added.</p>	
<p>8.3 References:</p> <p>8.3.1 Johns, S. H. et. al., "Spot Tests: A Color Chart Reference for Forensic Chemists," <i>Journal of Forensic Sciences</i>, Volume 24, No.3, July 1979, pp. 631-649.</p> <p>8.3.2 Moffat, A. C. , <i>Clarke's Isolation and Identification of Drugs</i>, The Pharmaceutical Press, London, 1986. Clarke, E. G. C., <i>Isolation and Identification of Drugs</i>, The Pharmaceutical Press, London, 1972, Vol. 1, 2.</p> <p>8.3.3 <i>Basic Training Program for Forensic Drug Chemists</i>, 2nd Edition , U. S. Dept. of Justice.</p> <p>8.3.4 Pitt, C. G. et. al., "The Specificity of the Duquenois Color Test for Marijuana and Hashish", <i>Journal of Forensic Science</i>, Vol. 17, No. 4, Oct. 1972, pp. 693-700.</p> <p>8.3.5 Oklahoma City Police Dept. Laboratory, Jane Bates et. al.</p> <p>8.3.6 Garrett, Steve A. et. al., "The Weber Tests," <i>N. E. A. F. S. Newsletter</i>, Vol. X , No. 2 June 1985.</p> <p>8.3.7 Gunn, John W. Jr., <i>Analysis of Drugs.</i>, United States Dept. of Justice.</p>	

8 COLOR TESTS	Page 8 of 8
<div>Division of Forensic Science</div> <div>CONTROLLED SUBSTANCES PROCEDURES MANUAL</div>	Amendment Designator: A
	Effective Date: 15-June-2005
<div> <div>8.3.8</div> <div>Feigel, Fritz., <i>Spot Tests in Organic Analysis</i>, 7th Edition, Elsevier Publishing Co., New York, 1966.</div> </div> <div> <div>8.3.9</div> <div>Butler, W. P., <i>Methods of Analysis</i>, Internal Revenue Service, Reprinted by BNDD, June 1967.</div> </div> <div> <div>8.3.10</div> <div><i>United States Pharmacopeia XIII</i></div> </div> <div> <div>8.3.11</div> <div>Connors, K. A., <i>A Textbook of Pharmaceutical Analysis</i>, 2nd Edition, John Wiley & Sons, New York, 1975</div> </div> <div> <div>8.3.12</div> <div>Beyer, E. and Dechert, D. D., "Identification of LSD and LSD Tartrate by Thin Layer Chromatography, "<i>Interbureau By-Lines</i> No. 1, BNDD, July 1967</div> </div> <div> <div>8.3.13</div> <div>Smith, Pam and Bozenko, Joseph, Jr. "New Presumptive Tests for GHB" <i>Microgram</i>, Vol. XXXV, No. 1, pp. 10-15 and DEA Communication, Pam Smith, DEA Special Testing and Research Laboratory, 2003</div> </div> <div> <div>8.3.14</div> <div><i>United States Pharmacopeia XXII</i></div> </div> <div> <div>8.3.15</div> <div><i>United States Pharmacopeia National Formulary</i>, United States Pharmacopeial Convention, Inc. Rockville, MD, 1995, p. 1722.</div> </div> <div> <div>◆ End</div> </div>	