

23 PSILOCYBIN AND PSILOCYN METHODOLOGY	Page 1 of 1
<div> Division of Forensic Science CONTROLLED SUBSTANCES PROCEDURES MANUAL </div>	Amendment Designator:
	Effective Date: 9-December-2003
<div> <p style="text-align: center;">23 PSILOCYBIN AND PSILOCYN METHODOLOGY</p> <p>23.1 Scheduling:</p> <p>23.1.1 Schedule I – psilocybin and psilocyn, which are found in mushrooms.</p> <p>23.1.2 The chemicals, rather than the botanical, are controlled.</p> <p>23.2 Extractions:</p> <p>23.2.1 Dry sample (in drying oven or microwave). Grind and soak in methanol for a period of 1 – 24 hours. Filter off plant material prior to analysis. This extraction method will allow for the analysis of either psilocyn or psilocybin.</p> <p>23.2.2 Acetic Acid Extraction Technique (recommended for mushrooms in chocolate or other matrices). This extraction will allow for the analysis of psilocyn.</p> <ul style="list-style-type: none"> • Dry approximately 3 grams material. (Easier to grind when dry.) • Grind up. • Let soak in 6% acetic acid for 30 minutes - 1 hour. • Filter off insoluble material. • Extract acid portion with three aliquots of CHCl₃.* (Discard CHCl₃.) • Basify acid portion with concentrated NH₄OH to pH 8 - 10. • Extract basic solution with three aliquots of CHCl₃.* • Combine aliquots of CHCl₃. • Evaporate CHCl₃ with air (<u>low heat</u>). • Resultant residue will yield psilocyn. <p>* Do not mix vigorously in the separatory funnel as an emulsion will probably form.</p> <p>23.3 Color Test Results:</p> <p>23.3.1 Ehrlich's- purple (positive for psilocyn and psilocybin).</p> <p>23.3.2 Weber Test- psilocyn - Fast Blue B or Fast Blue BB gives red color; addition of conc. HCl gives blue color.</p> <p>23.4 TLC:</p> <p>23.4.1 Baths: TLC1, TLC2, TLC3, TLC4 and TLC5 are recommended.</p> <p>23.4.2 Detection Sprays:</p> <p>23.4.2.1 p-DMAB and HCl - reddish violet area for psilocybin, blue for psilocyn.</p> <p>23.4.2.2 Weber Test (Fast Blue B or Fast Blue BB) can also be used as a TLC spray to detect psilocyn.</p> <p>23.5 GC:</p> <p>23.5.1 GC and GC/MS will give only psilocyn due to the dephosphorylation of the psilocybin caused by the GC injection port temperatures.</p> <p>23.5.2 After methanol extraction and drying, silylation with BSTFA prior to GC or GC/MS will allow differentiation of psilocybin from psilocyn. If psilocybin is to be confirmed, derivatization is required.</p> <p>23.6 FTIR: KBr pellet on extract from acetic acid extraction yields psilocyn.</p> <p style="text-align: right;">♦ End</p> </div>	