

24 CATHINONE METHODOLOGY	Page 1 of 1
Division of Forensic Science  CONTROLLED SUBSTANCES PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 9-December-2003
<p style="text-align: center;"><b>24 CATHINONE METHODOLOGY</b></p> <p><b>24.1 Scheduling:</b></p> <p>24.1.1 Schedule I – cathinone, which is found in <i>Catha Edulis</i> (Khat)</p> <p>24.1.2 The chemical, rather than the botanical, is controlled.</p> <p><b>24.2 Sample Handling:</b></p> <p>24.2.1 Suspected samples of Khat should be analyzed as soon as possible and the samples should be refrigerated, if possible.</p> <p>24.2.2 Cathinone degrades into cathine, a schedule IV substance.</p> <p><b>24.3 Extraction:</b></p> <p>24.3.1 Suggested <i>Catha Edulis</i> (Khat) Extraction Technique</p> <ul style="list-style-type: none"> <li>• Cut up approximately 5 grams of leaves and stems.</li> <li>• Homogenize (blender or bio-homogenizer) for 5 minutes in 100 mL of 0.1N HCl.</li> <li>• Filter off insoluble material.</li> <li>• Basify acid solution with 1N NaOH to pH 11-12.</li> <li>• Extract basic solution with aliquots of CHCl<sub>3</sub>.</li> <li>• Evaporate combined CHCl<sub>3</sub> aliquots with air (very low heat).</li> <li>• Analyze residue ASAP or refrigerate to avoid degradation.</li> <li>• Residue will give Cathinone and Cathine.</li> </ul> <p><b>24.4 TLC:</b></p> <p>24.4.1 Baths: TLC1, TLC2, TLC3, TLC4 and TLC5 are recommended.</p> <p>24.4.2 Detection Spray: Ninhydrin/heat gives red-brown color.</p> <p><b>24.5 References:</b></p> <p>24.5.1 Morselli <i>et al.</i>, “Gas-Chromatography/Mass Spectrometry Determination of the Active Principles of (<i>Catha Edulis</i>) African Vegetable,” <i>Microgram</i>, Vol. XXV, No. 11, November 1992, pp. 290-294.</p> <p>24.5.2 Lee, M. M., “The Identification of Cathinone in Khat (<i>Catha edulis</i>): A Time Study,” <i>Journal of Forensic Sciences</i>, JFSCA, Vol. 40, No. 1, January 1995, pp. 116-121.</p> <p style="text-align: right;">♦ End</p>	