

14 THC AND THC-CARBOXYLIC ACID CONFIRMATION IN BLOOD BY DPX EXTRACTION AND GCMS ANALYSIS	Page 1 of 5
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p>14 THC AND THC-CARBOXYLIC ACID CONFIRMATION IN BLOOD BY DPX EXTRACTION AND GCMS ANALYSIS</p> <p>14.1 Summary</p> <p>14.1.1 Delta-9-tetrahydrocannabinol (THC), principal metabolite 11-nor-9-carboxy-delta-9-tetrahydrocannabinol (THCA) and their deuterated internal standards are extracted and analyzed in biological samples using an acetonitrile precipitation. After centrifugation, the supernatant acetonitrile layer is made acidic and extracted with Disposable Pipette Extraction (DPX) tips. The DPX tips contain solid phase extraction (SPE) powder. The THC and THCA are adsorbed on to the SPE powder and extracted with Hexane/Ethyl Acetate. The extracted solvent is dried under nitrogen and derivatized with N,O-bis-(Trimethylsilyl)trifluoroacetamide (BSTFA) to form the trimethylsilyl derivatives of THC, THCA, deuterated THC, and deuterated THCA. The derivatized samples are injected into the GC/MS for confirmation by selected ion monitoring (SIM).</p> <p>14.2 Specimen Requirements</p> <p>14.2.1 1.0 mL of blood, biological fluid or tissue homogenate.</p> <p>14.3 Reagents And Standards</p> <p>14.3.1 Delta 9-THC, 1 mg/mL</p> <p>14.3.2 9-Carboxy-11-nor-delta 9-THC, 1 mg/mL</p> <p>14.3.3 Delta 9-THC-d₃, 100 µg/mL</p> <p>14.3.4 9-Carboxy-11-nor-delta 9-THC-d₃, 1 mg/mL</p> <p>14.3.5 N,O-bis(Trimethylsilyl)trifluoroacetamide with 1% Trimethylchlorosilane (BSTFA +1%TMCS)</p> <p>14.3.6 Acetonitrile</p> <p>14.3.7 Concentrated hydrochloric acid</p> <p>14.3.8 Methanol</p> <p>14.3.9 Hexane</p> <p>14.3.10 Ethyl acetate</p> <p>14.4 Solutions, Internal Standards, Calibrators, and Controls</p> <p>14.4.1 0.1 N Hydrochloric Acid: Pipet 8.25 mL of concentrated hydrochloric acid into a 1 L volumetric flask and QS to volume with dH₂O.</p> <p>14.4.2 Hexane/Ethyl Acetate (1:1): Mix 100 mL hexane with 100 mL ethyl acetate (v/v). Prepare fresh daily.</p> <p>14.4.3 Calibrators</p> <p>14.4.3.1 Working solution A (10 µg/mL): Pipet 100 µl each of THC and THCA stock solutions (1 mg/mL) into a 10 mL volumetric flask and QS to volume with methanol. Store in freezer.</p>	

14 THC AND THC-CARBOXYLIC ACID CONFIRMATION IN BLOOD BY DPX EXTRACTION AND GCMS ANALYSIS	Page 2 of 5
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<div> <div>14.4.3.2</div> <div>Working solution B (1.0 µg/mL): Pipet 1 mL of working solution A into a 10 mL volumetric flask and QS to volume with methanol. Store in freezer.</div> </div> <div> <div>14.4.3.3</div> <div>Working solution C (0.1 µg/mL): Pipet 1 mL working solution B into a 10 mL volumetric flask and QS to volume with methanol. Store in freezer.</div> </div> <div> <div>14.4.4</div> <div>Controls</div> <div> <div>14.4.4.1</div> <div>THC QC working solution A (10 µg/mL): Pipet 100 µL of THC stock solution (1 mg/mL) (<i>manufacturer or lot number different than that used for calibrators</i>) into a 10 mL volumetric flask and QS to volume with methanol. Store in freezer.</div> </div> <div> <div>14.4.4.2</div> <div>QC solution (0.2 µg/mL THC and 2.0 µg/mL THCA): Pipet 200 µl of THC QC working solution A into a 10 mL volumetric flask. Add 20 µl of THCA stock solution (1 mg/mL). QS to volume with methanol. Store in freezer.</div> </div> <div> <div>14.4.4.3</div> <div>Control (0.004 mg/L THC and 0.040 mg/L THCA): Measure 20 µL of QC solution into an appropriately labeled 16 x 125 mm screw cap tube containing 1 mL blank blood.</div> </div> <div> <div>14.4.4.4</div> <div>Negative blood control: Blood bank blood or equivalent previously determined not to contain THC or THCA.</div> </div> </div> <div> <div>14.4.5</div> <div>Internal Standard</div> <div> <div>14.4.5.1</div> <div>Internal standard working solution A (10 µg/mL THC-d₃ and THCA-d₃): Pipet 1 mL of the THC-d₃ stock solution (100 µg/mL) and 100 µl of the THCA-d₃ stock solution (1 mg/mL) into a 10 mL volumetric flask and QS to volume with methanol. Store in freezer.</div> </div> <div> <div>14.4.5.2</div> <div>Internal standard spiking solution (0.4 µg/mL THC-d₃ /THCA-d₃): Pipet 1 mL of the internal standard working solution A (10 µg/mL THC-d₃ /THCA-d₃) into a 25 mL volumetric flask and QS to volume with dH₂O. Prepare fresh daily.</div> </div> </div> <div> <div>14.4.6</div> <div>Calibrators. Pipet the following volumes of working standards into 16 x 125 mm tubes to achieve the following calibrator concentrations.</div> <div> <div>14.4.6.1</div> <div>Cal 1: 0.200 mg/L : 20 µL of 10 µg/mL working solution A</div> </div> <div> <div>14.4.6.2</div> <div>Cal 2: 0.100 mg/L: 10 µL of 10 µg/mL working solution A</div> </div> <div> <div>14.4.6.3</div> <div>Cal 3: 0.050 mg/L: 50 µL of 1.0 µg/mL working solution B</div> </div> <div> <div>14.4.6.4</div> <div>Cal 4: 0.010 mg/L: 10 µL of 1.0 µg/mL working solution B</div> </div> <div> <div>14.4.6.5</div> <div>Cal 5: 0.005 mg/L: 50 µL of 0.1 µg/mL working solution C</div> </div> <div> <div>14.4.6.6</div> <div>Cal 5: 0.002 mg/L: 20 µL of 0.1 µg/mL working solution C</div> </div> <div> <div>14.4.6.7</div> <div>Cal 5: 0.001 mg/L: 10 µL of 0.1 µg/mL working solution C</div> </div> <div> <div>14.4.6.8</div> <div>For each calibrator, evaporate standards to dryness under nitrogen and add 1 mL blank blood to each tube.</div> </div> </div>	

14 THC AND THC-CARBOXYLIC ACID CONFIRMATION IN BLOOD BY DPX EXTRACTION AND GCMS ANALYSIS

Page 3 of 5

Division of Forensic Science

TOXICOLOGY TECHNICAL PROCEDURES MANUAL

Amendment Designator:

Effective Date: 31-March-2004

14.5 Apparatus

- 14.5.1 Agilent GC/MSD, Chemstation software, compatible computer and printer
- 14.5.2 Screw cap test tubes, 16 x 125 mm
- 14.5.3 Screw cap test tubes, 13 x 100 mm
- 14.5.4 100 x 13 mm disposable test tubes
- 14.5.5 Centrifuge capable of 2,000-3,000 rpm
- 14.5.6 Nitrogen evaporator with heating block
- 14.5.7 Vortex mixer
- 14.5.8 12 x 75 mm disposable test tubes
- 14.5.9 GC auto sampler vials with inserts
- 14.5.10 DPX-CSP-05 (Disposable Pipette Extraction tips from EST Analytical)
- 14.5.11 GC/MSD Conditions. Instrument conditions may be changed to permit improved performance.
 - 14.5.11.1 Column: HP 5MS 25 m x 0.25 mm x 0.25 μ m
 - 14.5.11.2 Detector Temperature: 280° C
 - 14.5.11.3 SIM parameters

THC:	<u>386</u> ,371,303
THC-d ₃ :	<u>374</u>
THCA:	<u>371</u> ,473,488
THCA-d ₃ :	<u>374</u>
 - 14.5.11.4 Oven Program
 - Equilibration time: 0.50 minutes
 - Initial temp: 110° C
 - Initial time: 1 minutes
 - Ramp: 10° C/min
 - Final Temp: 290° C
 - Final Time: 9 minutes
 - Run Time: 28 minutes
 - 14.5.11.5 Inlet
 - Mode: Splitless
 - Temperature: 270° C
 - Injection volume: 1.0 μ L
 - Purge Time: ON at 1.0 minute

14 THC AND THC-CARBOXYLIC ACID CONFIRMATION IN BLOOD BY DPX EXTRACTION AND GCMS ANALYSIS	Page 4 of 5
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p>14.6 Procedure</p> <p>14.6.1 Label clean 16 x 125 mm screw cap tubes appropriately, blank, calibrators, controls and case sample IDs.</p> <p>14.6.2 Prepare calibrators and controls.</p> <p>14.6.3 Add 1 mL case specimen to the appropriately labeled tubes.</p> <p>14.6.4 Add 100 µL internal standard spiking solution to each tube. Vortex briefly.</p> <p>14.6.5 Acetonitrile extraction</p> <p>14.6.5.1 Add 2 mL acetonitrile to samples while vortexing.</p> <p>14.6.5.2 Vortex an additional 20 seconds.</p> <p>14.6.5.3 Centrifuge at 3000 rpm for 15 minutes.</p> <p>14.6.5.4 Pour off acetonitrile layer into 13 x 100 mL test tube</p> <p>14.6.5.5 Add 1 mL of 0.1 N hydrochloric acid to each sample extract</p> <p>14.6.6 Extraction with DPX tips</p> <p>14.6.6.1 Pre-wash the DPX tip with 500 µL hexane/ethyl acetate 1:1</p> <p>14.6.6.2 Draw sample extract into tip (allowing air to mix solution and powder)</p> <p>14.6.6.3 Let stand 20 seconds</p> <p>14.6.6.4 Expel into original test tube</p> <p>14.6.6.5 Redraw sample extract into tip</p> <p>14.6.6.6 Let stand 20 seconds</p> <p>14.6.6.7 Expel into original test tube</p> <p>14.6.6.8 Draw 500 µL hexane/ethyl acetate 1:1 into tip</p> <p>14.6.6.9 Let stand 20 seconds</p> <p>14.6.6.10 Expel this first extract into empty 12 x 75 mm test tube</p> <p>14.6.6.11 Draw another 500 µL hexane/ethyl acetate 1:1 into tip</p> <p>14.6.6.12 Let stand 20 seconds</p> <p>14.6.6.13 Expel this second extract into the 12 x 75 mm test tube containing the first extract</p> <p>14.6.7 Derivatization</p> <p>14.6.7.1 Remove water layer from bottom of combined extracts</p>	

14 THC AND THC-CARBOXYLIC ACID CONFIRMATION IN BLOOD BY DPX EXTRACTION AND GCMS ANALYSIS	Page 5 of 5
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<div data-bbox="308 344 932 621"> <p>14.6.7.2 Evaporate under nitrogen at 90° C</p> <p>14.6.7.3 Add 40 µL ethyl acetate to test tubes and vortex</p> <p>14.6.7.4 Transfer ethyl acetate to GC autosampler vials</p> <p>14.6.7.5 Add 20 µL BSTFA +1% TCMS</p> <p>14.6.7.6 Cap and inject 2 µL on GC/MSD in SIM mode</p> </div> <div data-bbox="175 651 375 680"> 14.7 Calculation </div> <div data-bbox="219 711 1534 865"> <p>14.7.1 Calculate the concentrations by interpolation of a linear plot of the response curve based on peak heights (or areas) ratios (using the target ions underlined above) versus calibrator concentration.</p> <p>14.7.2 Qualifier ion ratio range. The qualifier ion ratio range is calculated by determining the mean ion ratio from all calibrators used in the calibrations curve. Each drug has two qualifier ions.</p> </div> <div data-bbox="175 894 423 924"> 14.8 Quality Control </div> <div data-bbox="219 955 690 984"> <p>14.8.1 See Toxicology Quality Guidelines</p> </div> <div data-bbox="175 1014 367 1043"> 14.9 References </div> <div data-bbox="219 1075 1544 1262"> <p>14.9.1 Rapid and Sensitive Analysis of THC and COOH-THC in Whole Blood. Brandi L. Clelland and William E. Brewer, Ph.D., DFTCB, SOFT presentation October 2001</p> <p>14.9.2 Personal communiqué with William E. Brewer, Ph.D., January to August 2003</p> <p>14.9.3 Dwight Flammia, Ph.D., Randall Edwards, and Terri Woods, in-house development.</p> </div>	