

18 PANCURONIUM AND ROCURONIUM ANALYSIS BY LCMS	Page 1 of 5
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p style="text-align: center;">18 PANCURONIUM AND ROCURONIUM ANALYSIS BY LCMS</p> <p>18.1 Summary</p> <p>18.1.1 Quaternary nitrogen muscle relaxants (pancuronium and rocuronium) are extracted from biological samples using acetonitrile precipitation, solid phase extraction (SPE) and analyzed by high performance liquid chromatography-electrospray ionization mass spectrometry (LC-ESI-MS). Pancuronium, rocuronium and the internal standard, verapamil, are extracted and analyzed simultaneously. LC-ESI-MS analysis is achieved with a 10-90% acetonitrile gradient containing 0.1% trifluoroacetic acid.</p> <p>18.2 Specimen Requirements</p> <p>18.2.1 2 mL blood, biological fluid or tissue homogenate.</p> <p>18.3 Reagents and Standards</p> <p>18.3.1 Ammonium carbonate</p> <p>18.3.2 Ammonium acetate</p> <p>18.3.3 Methanol</p> <p>18.3.4 Acetonitrile</p> <p>18.3.5 Hexane</p> <p>18.3.6 Pancuronium bromide (Pavulon®, Baxter, 1 mg/mL)</p> <p>18.3.7 Pancuronium bromide (Sigma)</p> <p>18.3.8 Rocuronium bromide (Zemuron®, Organon, 10 mg/mL)</p> <p>18.3.9 Verapamil hydrochloride (e.g. Alltech)</p> <p>18.3.10 Trifluoroacetic acid</p> <p>18.3.11 Glacial Acetic Acid</p> <p>18.3.12 Potassium hydroxide</p> <p>18.4 Solutions, Internal Standards, Calibrators and Controls</p> <p>18.4.1 1.0 M Acetic Acid: Pipet 57.5 mL glacial acetic acid into a 1L volumetric flask. QS to volume with dH₂O.</p> <p>18.4.2 Ammonium Acetate Buffer (pH 5.0, 50mM): Weigh 3.85 g ammonium acetate. Transfer to 1 L volumetric flask and add approximately 900 mL dH₂O. Adjust pH to 5.0 with 1.0 M acetic acid. QS to volume with dH₂O.</p> <p>18.4.3 5.0 M Potassium hydroxide : Weigh 28 g potassium hydroxide. Transfer to 100 mL volumetric flask and QS to volume with dH₂O.</p> <p>18.4.4 Ammonium Carbonate Buffer (pH 9.3, 0.01M): Weigh 0.47 g ammonium carbonate. Transfer to 500 mL volumetric flask and add approximately 450 mL dH₂O. Adjust pH to 9.3 with 5 M potassium hydroxide. QS to volume with dH₂O.</p>	

18 PANCURONIUM AND ROCURONIUM ANALYSIS BY LCMS		Page 2 of 5
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL		Amendment Designator:
		Effective Date: 31-March-2004
18.4.5	0.1M acetic acid in methanol: Pipet 10 mL 1.0 M acetic acid into a 100 mL volumetric flask. QS to volume with methanol. Prepare fresh daily.	
18.4.6	Working standard solutions for pancuronium and rocuronium	
18.4.6.1	100 µg/mL pancuronium/rocuronium working solution: Pipet 1 mL of 1 mg/mL stock solution of pancuronium and 100 µL of 10 mg/mL stock solution of rocuronium into a 10 mL volumetric flask. QS to volume with acetonitrile. Prepare fresh daily.	
18.4.6.2	10 µg/mL pancuronium/rocuronium working solution: Pipet 1 mL of 0.1 mg/mL pancuronium/rocuronium working solution into a 10 mL volumetric flask. QS to volume with acetonitrile. Prepare fresh daily.	
18.4.7	Quality Control (QC) solution	
18.4.7.1	100 µg/mL pancuronium/rocuronium QC solution: Pipet 100 µL of 1 mg/mL stock solution of pancuronium and 10 µL of 10 mg/mL stock solution of rocuronium (both from different manufacturer or lot number than standards). Add 890 µL acetonitrile. Prepare fresh daily.	
18.4.8	Internal standard working solution	
18.4.8.1	10 µg/mL verapamil: Pipet 100 µL of 1 mg/mL verapamil stock solution into 10 mL volumetric flask and QS to volume with acetonitrile. Store in freezer.	
18.4.9	Calibrators. To prepare the following calibration curve, pipet the following volumes into appropriately labeled 16 x 125 mm screw cap tubes	
18.4.9.1	Cal 1: 10 mg/L pancuronium/rocuronium:	20 µL each of 1 mg/mL pancuronium and rocuronium
18.4.9.2	Cal 2: 5 mg/L pancuronium/rocuronium:	10 µL each of 1 mg/mL pancuronium and rocuronium
18.4.9.3	Cal 3: 2 mg/L pancuronium/rocuronium:	400 µL of 0.01 mg/mL pancuronium/rocuronium solution
18.4.9.4	Cal 4: 1 mg/L pancuronium/rocuronium:	200 µL of 0.01 mg/mL pancuronium/rocuronium solution
18.4.9.5	Cal 5: 0.5mg/L pancuronium/rocuronium:	100 µL of 0.01 mg/mL pancuronium/rocuronium solution
18.4.9.6	Cal 6: 0.1mg/L pancuronium/rocuronium:	20 µL of 0.01 mg/mL pancuronium/rocuronium solution
18.4.9.7	Add 2 mL blank blood to each tube.	
18.4.10	Pancuronium and rocuronium control (QC)	
18.4.10.1	1mg/L pancuronium/rocuronium QC: Pipet 200 µL of 0.01 mg/mL pancuronium/rocuronium QC solution into appropriately labeled 16 x 125 mm screw cap tube and add 2 mL blank blood.	
18.4.10.2	Negative control: blood bank blood (or equivalent) previously determined not to contain rocuronium, pancuronium or verapamil.	
18.5	Apparatus	
18.5.1	Test tubes, 16 x 125 mm round bottom, screw cap with Teflon caps	
18.5.2	Test tubes, 16 x 114 mm glass centrifuge, conical bottom	
18.5.3	Centrifuge capable of 2000-3000 rpm	
18.5.4	Nitrogen evaporator with heating block	
18.5.5	Vortex mixer	

18 PANCURONIUM AND ROCURONIUM ANALYSIS BY LCMS		Page 3 of 5												
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL		Amendment Designator:												
		Effective Date: 31-March-2004												
18.5.6 GC autosampler vials with inserts														
18.5.7 Solid Phase Extraction manifold														
18.5.8 Strata C18-E SPE columns (6 ml, bed volume 500 mg), Phenomenex														
18.5.9 LC/MS: Agilent Model 1100 LC-MSD														
18.5.10 LCMS Instrument Conditions. The following instrument conditions may be modified to adjust or improve separation and sensitivity.														
18.5.10.1 Elution conditions:														
18.5.10.1.1 Column: Agilent Hypersil BDS 125 mm X 3 mm, 3 µm particle size														
18.5.10.1.2 Column thermostat: 35° C														
18.5.10.1.3 Solvent A: Water with 0.1%Trifluoroacetic acid														
18.5.10.1.4 Solvent B: Acetonitrile														
18.5.10.1.5 Gradient elution, stop time: 13.00 min														
<table><tr><td>Time</td><td>Solv. B</td><td>Flow</td></tr><tr><td>0.00</td><td>10.0</td><td>0.65</td></tr><tr><td>8.00</td><td>90.0</td><td>0.65</td></tr><tr><td>9.00</td><td>10</td><td>0.65</td></tr></table>			Time	Solv. B	Flow	0.00	10.0	0.65	8.00	90.0	0.65	9.00	10	0.65
Time	Solv. B	Flow												
0.00	10.0	0.65												
8.00	90.0	0.65												
9.00	10	0.65												
18.5.10.2 Spray Chamber														
<ul style="list-style-type: none">• Ionization Mode: Electrospray• Gas Temperature: 350° C• Drying Gas (N₂): 12.0 L/min• Nebulizer pressure: 35 psig• Vcap (Positive): 4000 V														
18.5.11 Selected Ion Monitoring														
18.5.11.1 Polarity: Positive														
18.5.11.2 SIM parameters (<u>quantitation ions</u>)														
Rocuronium ions: 358, 413, 487, <u>529</u>														
Pancuronium ions: 412, 472, 571, <u>685</u>														
Verapamil IS ions: 165, 303, <u>455</u>														
18.6 Procedure														
18.6.1 Label 16 x 125 mm screw cap tubes appropriately (blank, calibrators, controls and case sample IDs).														
18.6.2 Prepare calibrators and controls.														

18 PANCURONIUM AND ROCURONIUM ANALYSIS BY LCMS	Page 4 of 5
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p>18.6.3 Pipet 2 mL of each case specimen into appropriately labeled tubes.</p> <p>18.6.4 Add 100 µL 0.01 mg/mL verapamil working internal standard solution to each tube. Vortex briefly.</p> <p>18.6.5 Slowly, add dropwise 2 mL cold (freezer temperature) acetonitrile to each tube while vortexing. Continuous vortexing, not mere mixing, is essential.</p> <p>18.6.6 Vortex an additional 30 seconds.</p> <p>18.6.7 Cap tubes.</p> <p>18.6.8 Place tubes in freezer for at least 30 minutes to facilitate separation.</p> <p>18.6.9 Centrifuge at approximately 2500 rpm for 15 minutes.</p> <p>18.6.10 Transfer top (acetonitrile) layer to clean 16x125 mL tubes taking care not to transfer any lower layers.</p> <p>18.6.11 Add 4 mL 0.01 M ammonium carbonate buffer to each tube. Vortex briefly.</p> <p>18.6.12 Prepare C18-E SPE columns</p> <p style="padding-left: 40px;">18.6.12.1 Add 2 mL methanol to each column. Aspirate slowly under vacuum (approx 1 mL/min).</p> <p style="padding-left: 40px;">18.6.12.2 Add 4 mL ammonium carbonate buffer to each column. Aspirate slowly under vacuum.</p> <p>18.6.13 Load buffered sample supernatants to columns. Aspirate slowly under vacuum.</p> <p>18.6.14 Wash columns with 4 mL ammonium carbonate. Aspirate slowly under vacuum.</p> <p>18.6.15 Repeat wash with 4 mL ammonium carbonate. Aspirate slowly under vacuum.</p> <p>18.6.16 Add 500 µL hexane to each column. Aspirate. Dry the columns at > 10 inches of Hg for at least 10 minutes.</p> <p>18.6.17 Elute drugs by adding 4 mL of freshly prepared 0.01M acetic acid in methanol. Collect eluants under gravity (no vacuum) into conical bottom screw cap tubes.</p> <p>18.6.18 Evaporate eluants to dryness at approximately 50° C under nitrogen.</p> <p>18.6.19 Reconstitute samples in 1 mL acetonitrile. Vortex briefly to ensure recovery of drugs from glass tube.</p> <p>18.6.20 Evaporate samples again to dryness at approximately 50° C under nitrogen.</p> <p>18.6.21 Reconstitute samples in 100 µL acetonitrile. Vortex briefly. Transfer to GC microvials.</p> <p>18.6.22 Inject 5 µL of each sample on LC/MS in the API-ES/SIM Mode</p> <p>18.7 Calculation</p> <p>18.7.1 Drug concentrations are calculated by linear regression analysis using the ChemStation software.</p> <p>18.8 Quality Control and Reporting</p> <p>18.8.1 See Toxicology Quality Guidelines</p> <p>18.9 References</p>	

18 PANCURONIUM AND ROCURONIUM ANALYSIS BY LCMS	Page 5 of 5
<div>Division of Forensic Science</div> <div>TOXICOLOGY TECHNICAL PROCEDURES MANUAL</div>	Amendment Designator:
	Effective Date: 31-March-2004
<p>18.9.1 CHM Kerskes, KJ Lusthof, PGM Zweipfenning and JP Franke. The Detection and Identification of Quaternary Nitrogen Muscle Relaxants in Biological Fluids and Tissues by Ion-Trap LC-ESI-MS. <i>J Anal Tox</i> 26:29-34, 2002.</p> <p>18.9.2 C Ferenc, C Enjalbal, P Sanchez, F Bressolle, M Audran, J Martinez and JL Aubagnac. Quantitative Determination of Rocuronium in Human Plasma by LC-ESI-MS. <i>J Chrom A</i> 910: 61-67, 2001.</p> <p>18.9.3 M Zecevic, LJ Zivanovic and A Stojkovic. Validation of HPLC Method for the Determination of Pancuronium in Pavulon Injections. <i>J Chrom A</i> 949: 61-64, 2002.</p> <p>18.9.4 L Gao, I Ramzan and B Baker. GCMS Assay for Rocuronium with Potential for Quantifying its Metabolite, 17-desacetylrocuronium, in Human Plasma. <i>J Chrom B</i> 757: 207-214, 2001.</p> <p>18.9.5 J Pearson and R Steiner, in-house development.</p>	