

20 CARISOPRODOL AND MEPROBAMATE QUANTITATION BY GC/FID	Page 1 of 5
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p style="text-align: center;">20 CARISOPRODOL AND MEPROBAMATE QUANTITATION BY GC/FID</p> <p>20.1 Summary</p> <p>20.1.1 Biological samples are buffered with phosphate buffer (pH 7) and extracted with a mixture of hexane and ethyl acetate. The extract is washed with hexane and reconstituted with toluene/hexane/isoamyl alcohol. An aliquot is injected into a GC equipped with an FID detector for quantitation of carisoprodol and meprobamate. The aliquot can be subsequently injected into a GCMS for confirmation, if necessary.</p> <p>20.2 Specimen Requirements</p> <p>20.2.1 200 µl biological fluid or comparable amount of tissue dilutions/homogenates</p> <p>20.3 Reagents and Standards</p> <p>20.3.1 Carisoprodol</p> <p>20.3.2 Meprobamate</p> <p>20.3.3 Cyclopal (cyclopentalbarbital)</p> <p>20.3.4 Trimethyl sulfonium iodide</p> <p>20.3.5 Disodium phosphate (Na_2HPO_4)</p> <p>20.3.6 Monosodium phosphate (NaH_2PO_4)</p> <p>20.3.7 Silver oxide</p> <p>20.3.8 Hexane</p> <p>20.3.9 Isoamyl alcohol</p> <p>20.3.10 Methanol</p> <p>20.3.11 Toluene</p> <p>20.3.12 Ethyl acetate</p> <p>20.3.13 Acetonitrile</p> <p>20.3.14 N-chlorobutane</p> <p>20.4 Solutions, Internal Standard, Calibrators and Controls</p> <p>20.4.1 0.1 M disodium phosphate: Weigh 1.70 g of disodium phosphate and transfer to 1L volumetric flask. QS to volume with dH_2O.</p> <p>20.4.2 0.1 M monosodium phosphate: Weigh 12.14 g monosodium phosphate and transfer to a 1 L volumetric flask. QS to volume with dH_2O.</p>	

20 CARISOPRODOL AND MEPROBAMATE QUANTITATION		Page 2 of 5
BY GC/FID		
Division of Forensic Science		Amendment Designator:
TOXICOLOGY TECHNICAL PROCEDURES MANUAL		Effective Date: 31-March-2004
20.4.3	0.1 M sodium phosphate buffer (pH 7.0): Mix 500 mL 0.1 M disodium phosphate with approximately 250 mL 0.1 M monosodium phosphate. Adjust pH to 7.0 ± 0.1 with 0.1 M monosodium phosphate (lowers pH) or 0.1 M disodium phosphate (raises pH).	
20.4.4	Toluene:Hexane:Isoamyl Alcohol (THIA) (78:20:2, v:v:v) Mix 78 mL toluene, 20 mL hexane and 2 mL isoamyl alcohol	
20.4.5	Hexane/ethyl acetate (50:50, v:v) extraction solvent: Mix 50 mL hexane with 50 mL ethyl acetate	
20.4.6	Trimethyl sulfonium hydroxide derivatizing reagent. Add 6.12 g trimethylsulfonium iodide, 7.39 g silver oxide, and 15 mL methanol to a 25 mL teflon capped test tube covered with aluminum foil (light sensitive reaction). Rotate for four or more hours, centrifuge, and decant the supernatant to an aluminum foil covered test tube. Store in freezer.	
20.4.7	Methanol/ dH ₂ O (50:50, v:v) Mix 50 mL methanol with 50 mL dH ₂ O	
20.4.8	Drug stock solutions:	
20.4.8.1	If 1 mg/mL commercially prepared stock solutions are not available, prepare 1 mg/mL solutions from powders. Weigh 10 mg of the free drug, transfer to a 10 mL volumetric flask and QS to volume with methanol.	
20.4.8.2	Internal standard solution (methylated cyclopal (cyclopentalbarbital)): To 2 mL of 1 mg/mL cyclopal (in methanol), add 1 mL TMSH. Cap and heat at 60° C for 2 hours. Let sit at room temperature overnight. Evaporate under nitrogen. Reconstitute with 10 mL methanol/ dH ₂ O (50:50, v:v).	
20.4.9	Blood calibrators, standards, and controls preparation:	
20.4.9.1	To prepare the calibration curve, pipet the following volumes of 1 mg/mL carisoprodol and meprobamate stock solutions into appropriately labeled 13 x 100 mm screw cap test tubes	
	<ul style="list-style-type: none"> ▪ 100 mg/L Calibrator ▪ 50 mg/L Calibrator ▪ 20 mg/L Calibrator ▪ 10 mg/L Calibrator ▪ 5 mg/L Calibrator ▪ 2 mg/L Calibrator 	<ul style="list-style-type: none"> 300 µL each of carisoprodol and meprobamate 150 µL each of carisoprodol and meprobamate 60 µL each of carisoprodol and meprobamate 30 µL each of carisoprodol and meprobamate 15 µL each of carisoprodol and meprobamate 6 µL each of carisoprodol and meprobamate
20.4.9.2	Evaporate standards to dryness under nitrogen. Add 3 mL blank blood to each tube. Calibrators may be stored in refrigerator for up to 1 year after preparation.	
20.4.9.3	Controls	
	20.4.9.3.1	Negative control. Blood bank blood (or comparable) determined not to contain carisoprodol or meprobamate
	20.4.9.3.2	Positive control. In house control containing each analyte of interest from a different lot number or manufacturer than standards, or prepared by a chemist different than the one performing the extraction.
20.5 Apparatus		
20.5.1	Agilent GC/MSD, Chemstation software (for confirmation, if necessary)	

20 CARISOPRODOL AND MEPROBAMATE QUANTITATION BY GC/FID	Page 3 of 5
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p>20.5.2 Agilent GC equipped with Flame Ionization Detector, Chemstation software, compatible computer & printer</p> <p>20.5.3 Test tubes, 13 x 100 mm round bottom, screw cap tubes, borosilicate glass with Teflon caps</p> <p>20.5.4 Test tubes, 16 x 114 mm (10 mL) glass centrifuge, conicals</p> <p>20.5.5 Centrifuge capable of 2,000 – 3,000 rpm</p> <p>20.5.6 Vortex mixer</p> <p>20.5.7 Evaporator/concentrator</p> <p>20.5.8 GC autosampler vials and inserts</p> <p>20.5.9 Test tube rotator</p> <p>20.5.10 GC/FID parameters. Conditions may be changed to permit improved performance.</p> <p>20.5.10.1 Oven program.</p> <ul style="list-style-type: none"> • Equilibration time: 0.50 minutes • Initial temp: 110° C • Initial time: 1.0 minutes • Ramp: 20° C/min • Final Temp: 260° C • Final Time: 1.5 minutes • Run Time: 15 minutes <p>20.5.10.2 Inlet.</p> <ul style="list-style-type: none"> • Mode: Splitless • Temperature: 250° C • Constant pressure: 25 psi • Purge flow: 1.9 mL/min • Total flow: 6.1 mL/min • Injection volume: 1.0 µL <p>20.5.10.3 Detector.</p> <ul style="list-style-type: none"> • Temperature: 290° C • Hydrogen flow: 50 mL/min • Air flow: 450 mL/min • Mode: Constant makeup flow • Makeup flow: 45 mL/min <p>20.5.10.4 Column: HP-5 30 m x 0.25 mm x 0.25 µm.</p> <p>20.5.11 GC/MSD parameters. Conditions may be changed to permit improved performance.</p> <p>20.5.11.1 Acquisition Mode: Scan (50 – 550 amu)</p> <p>20.5.11.2 Column: HP 5MS 25 m x 0.25 mm x 0.25 µm</p>	

20 CARISOPRODOL AND MEPROBAMATE QUANTITATION BY GC/FID	Page 4 of 5
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p>20.5.11.3 Detector Temperature: 280° C</p> <p>20.5.11.4 Oven Program</p> <ul style="list-style-type: none"> • Equilibration time: 0.50 minutes • Initial temp: 110° C • Initial time: 1 minutes • Ramp: 10° C/min • Final Temp: 290° C • Final Time: 9 minutes • Run Time: 28 minutes <p>20.5.11.4.1 Inlet</p> <ul style="list-style-type: none"> • Mode: Splitless • Temperature: 270° C • Injection volume: 1.0 µL • Purge Time: ON at 1.0 minute <p>20.6 Procedure</p> <p>20.6.1 Label clean 13 x 100 mm screw cap tubes accordingly, negative, calibrators, control(s) and case sample IDs.</p> <p>20.6.2 Prepare calibrators and controls</p> <p>20.6.3 Pipet 200µL of each calibrator, control, negative and case samples into appropriately labeled tubes.</p> <p>20.6.4 Add 30 µL methylated cyclopal internal standard to each tube.</p> <p>20.6.5 Add 0.5 mL sodium phosphate buffer (pH 7) to each tube.</p> <p>20.6.6 Add 3 mL extract solvent (hexane/ethyl acetate) to each tube.</p> <p>20.6.7 Cap and rotate tubes for 15 minutes.</p> <p>20.6.8 Centrifuge at approximately 2500 rpm for 15 minutes. Transfer organic upper layer to clean 10 mL conical bottom centrifuge tubes. Discard lower layers.</p> <p>20.6.9 Evaporate samples to dryness under nitrogen at 50-60° C.</p> <p>20.6.10 Reconstitute samples with 0.2 mL acetonitrile. Vortex briefly.</p> <p>20.6.11 Add 1 mL hexane to each tube. Vortex each sample for 30 seconds.</p> <p>20.6.12 Centrifuge at approximately 2500 rpm for 5 minutes.</p> <p>20.6.13 Aspirate (and discard) upper (hexane) layer.</p> <p>20.6.14 Evaporate lower (acetonitrile) layer under nitrogen at 50-60° C.</p> <p>20.6.15 Reconstitute samples with 75 µL of toluene/hexane/isamyl alcohol solvent or n-chlorobutane and vortex briefly</p> <p>20.6.16 Transfer samples to appropriately labeled GC vials and inject 1-2 µl on GC-FID.</p>	

<p>20 CARISOPRODOL AND MEPROBAMATE QUANTITATION BY GC/FID</p>	<p>Page 5 of 5</p>
<p>Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL</p>	<p>Amendment Designator:</p>
	<p>Effective Date: 31-March-2004</p>
<p>20.6.17 Save remainder of reconstituted samples for confirmation by GC-MSD (if not already confirmed).</p> <p>20.7 Calculation</p> <p>20.7.1 Calculate the concentrations by interpolation of a linear plot of the response curve based on peak height (or area) ratios versus calibrator concentration.</p> <p>20.8 Quality Control and Reporting</p> <p>20.8.1 See Toxicology Quality Guidelines</p> <p>20.9 References</p> <p>20.9.1 in house development, T England</p>	