

26 ACETAMINOPHEN QUANTITATION AND CONFIRMATION BY LCMS	Page 1 of 3
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL	Amendment Designator:
	Effective Date: 31-March-2004
<p style="text-align: center;">26 ACETAMINOPHEN QUANTITATION AND CONFIRMATION BY LCMS</p> <p>26.1 Summary</p> <p>26.1.1 Acetaminophen is extracted from biological samples by making the samples weakly acidic/neutral with saturated sodium carbonate buffer and extracting with hexane/ethyl acetate. An aliquot of the extract is analyzed by high performance liquid chromatography-electrospray ionization mass spectrometry (LC-ESI-MS).</p> <p>26.2 Specimen Requirements</p> <p>26.2.1 0.5 mL blood, fluid or tissue homogenate.</p> <p>26.3 Reagents and Standards</p> <p>26.3.1 Acetaminophen, 1 mg/mL</p> <p>26.3.2 Phenacetin, 1 mg/mL</p> <p>26.3.3 Sodium carbonate (NaHCO₃)</p> <p>26.3.4 Hexane</p> <p>26.3.5 Ethyl acetate</p> <p>26.3.6 Methanol</p> <p>26.3.7 Acetic Acid</p> <p>26.4 Solutions, Internal Standard, Calibrators and Controls</p> <p>26.4.1 Saturated sodium carbonate solution. Add sodium carbonate to dH₂O until no more dissolves after shaking vigorously.</p> <p>26.4.2 Hexane/Ethyl acetate (50:50, v:v). Mix 100 mL hexane with 100 mL ethyl acetate.</p> <p>26.4.3 Drug stock solutions:</p> <p>26.4.3.1 If 1 mg/mL commercially prepared stock solutions are not available, prepare 1 mg/mL solutions from powders. Weigh 10 mg of the free drug, transfer to a 10 mL volumetric flask and QS to volume with methanol. Note: If using the salt form, determine the amount of the salt needed to equal 10 mg of the free drug, and weigh this amount. Stock solutions are stored capped in a refrigerator and are stable for 2 years.</p> <p>26.4.4 Working acetaminophen standard solution (0.2 mg/mL). Pipet 1 mL of 1 mg/mL acetaminophen stock solution into a 5 mL volumetric flask and QS to volume with methanol.</p> <p>26.4.5 Working internal standard solution (50 µg/mL phenacetin): Pipet 500 µL of the 1 mg/mL stock solution of phenacetin into a 10 mL volumetric flask and QS to volume with methanol.</p> <p>26.4.6 To prepare the calibration curve, pipet the following volumes of the 0.2 mg/mL acetaminophen stock solution into appropriately labeled 16 x 125 mm screw cap test tubes. Add 0.5 mL blank blood to obtain the final concentrations listed below.</p>	

26 ACETAMINOPHEN QUANTITATION AND CONFIRMATION BY LCMS		Page 2 of 3																		
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL		Amendment Designator:																		
		Effective Date: 31-March-2004																		
<table> <tr> <td></td><td>Amount of 0.2 mg/mL acetaminophen (µL)</td><td>Final concentration of acetaminophen (mg/L)</td></tr> <tr> <td></td><td>250</td><td>100</td></tr> <tr> <td></td><td>150</td><td>60</td></tr> <tr> <td></td><td>100</td><td>40</td></tr> <tr> <td></td><td>50</td><td>20</td></tr> <tr> <td></td><td>25</td><td>10</td></tr> </table>				Amount of 0.2 mg/mL acetaminophen (µL)	Final concentration of acetaminophen (mg/L)		250	100		150	60		100	40		50	20		25	10
	Amount of 0.2 mg/mL acetaminophen (µL)	Final concentration of acetaminophen (mg/L)																		
	250	100																		
	150	60																		
	100	40																		
	50	20																		
	25	10																		
26.4.7 Controls																				
26.4.7.1 Acetaminophen Control. Control may be from an external source or prepared in house using drugs from different manufacturers, lot numbers or prepared by a chemist different than the individual performing the extraction.																				
26.4.7.2 Negative control. Blood bank blood or equivalent determined not to contain acetaminophen or phenacetin.																				
26.5 Apparatus																				
26.5.1 Test tubes, 16 x 125 mm, round bottom, borosilicate glass with Teflon caps																				
26.5.2 Test tubes, 16 x 114 mm, glass centrifuge, conical bottom																				
26.5.3 Centrifuge capable of 2000-3000 rpm																				
26.5.4 Nitrogen evaporator with heating block																				
26.5.5 Vortex mixer																				
26.5.6 GC autosampler vials with inserts																				
26.5.7 LC/MS: Agilent Model 1100 LC-MSD																				
26.5.7.1 LCMS Instrument Conditions. The following instrument conditions may be modified to adjust or improve separation and sensitivity.																				
26.5.7.1.1 Elution Condicitons																				
26.5.7.1.1.1 Column: Agilent Hypersil BDS 125 mm X 3 mm, 3 µM particle size																				
26.5.7.1.1.2 Column thermostat: 30° C																				
26.5.7.1.1.3 Solvent A: 45% water with 1% acetic acid																				
26.5.7.1.1.4 Solvent B: 55% methanol																				
26.5.7.1.1.5 Isocratic elution, stop time: 4.0 min																				
26.5.7.1.2 Spray Chamber																				
26.5.7.1.2.1 Ionization Mode: Electrospray																				
26.5.7.1.2.2 Gas Temperature: 350° C																				

26 ACETAMINOPHEN QUANTITATION AND CONFIRMATION BY LCMS						Page 3 of 3																																																								
Division of Forensic Science TOXICOLOGY TECHNICAL PROCEDURES MANUAL						Amendment Designator:																																																								
						Effective Date: 31-March-2004																																																								
26.5.7.1.2.3 Drying Gas (N ₂): 12.0 L/min 26.5.7.1.2.4 Nebulizer pressure: 30 psig 26.5.7.1.2.5 Vcap (Positive): 4000 V 26.5.7.1.3 Selected Ion Monitoring (<u>quantitation ions</u>) 26.5.7.1.3.1 Polarity: Positive 26.5.7.1.3.2 Injection volume: 2 µL																																																														
<table border="1"> <thead> <tr> <th>Time (min)</th><th>Group Name</th><th>SIM Ion</th><th>Frag-Mentor</th><th>Gain EMV</th><th>SIM Resol.</th><th>Actual Dwell</th></tr> </thead> <tbody> <tr> <td>0.00</td><td>Acetaminophen</td><td>93</td><td>175</td><td>1.0</td><td>Low</td><td>218</td></tr> <tr> <td></td><td></td><td><u>110</u></td><td>175</td><td></td><td>218</td><td></td></tr> <tr> <td></td><td></td><td>134</td><td>175</td><td></td><td>218</td><td></td></tr> <tr> <td></td><td></td><td>152</td><td>175</td><td></td><td>218</td><td></td></tr> <tr> <td>2.30</td><td>Phenacetin</td><td>110</td><td>170</td><td>1.0</td><td>Low</td><td>292</td></tr> <tr> <td></td><td></td><td>138</td><td>170</td><td></td><td>292</td><td></td></tr> <tr> <td></td><td></td><td><u>180</u></td><td>170</td><td></td><td>292</td><td></td></tr> </tbody> </table>						Time (min)	Group Name	SIM Ion	Frag-Mentor	Gain EMV	SIM Resol.	Actual Dwell	0.00	Acetaminophen	93	175	1.0	Low	218			<u>110</u>	175		218				134	175		218				152	175		218		2.30	Phenacetin	110	170	1.0	Low	292			138	170		292				<u>180</u>	170		292		
Time (min)	Group Name	SIM Ion	Frag-Mentor	Gain EMV	SIM Resol.	Actual Dwell																																																								
0.00	Acetaminophen	93	175	1.0	Low	218																																																								
		<u>110</u>	175		218																																																									
		134	175		218																																																									
		152	175		218																																																									
2.30	Phenacetin	110	170	1.0	Low	292																																																								
		138	170		292																																																									
		<u>180</u>	170		292																																																									
26.6 Procedure																																																														
26.6.1 Label clean 16 x 125 mm screw cap tubes appropriately with calibrators, controls and case sample IDs.																																																														
26.6.2 Prepare calibrators and controls.																																																														
26.6.3 Add 0.5 mL case specimens to the appropriately labeled tubes.																																																														
26.6.4 Add 50 µL of the 50 µg/mL phenacetin internal standard working solution to each tube for a final concentration of 5 mg/L.																																																														
26.6.5 Add 1 mL saturated sodium carbonate buffer and 3 mL extraction solvent (50:50 hexane/ethyl acetate) to each tube.																																																														
26.6.6 Cap and rotate tubes for 30 minutes.																																																														
26.6.7 Centrifuge at approx 2500 rpm for 15 minutes. Transfer organic (upper) layer to appropriately labeled conical bottom test tubes.																																																														
26.6.8 Evaporate samples to dryness at approximately 60° C under nitrogen.																																																														
26.6.9 Reconstitute samples in 100 µL methanol. Transfer to GC autosampler vials for analysis by LCMS.																																																														
26.7 Calculation																																																														
26.7.1 Drug concentrations are calculated by linear regression analysis using the ChemStation software.																																																														
26.8 Quality Control and Reporting																																																														
26.8.1 See Toxicology Quality Guidelines																																																														
26.9 References																																																														
26.9.1 M Kennedy, D Sullivan and R Steiner, in house development																																																														