

28 VALPROIC ACID CONFIRMATION BY GCMS	Page 1 of 3
<div> <div>Division of Forensic Science</div> <div>TOXICOLOGY TECHNICAL PROCEDURES MANUAL</div> </div>	Amendment Designator:
	Effective Date: 31-March-2004
<div> <div>28 VALPROIC ACID CONFIRMATION BY GCMS</div> <div> <div>28.1 Summary</div> <div> <div>28.1.1 Biological samples are slightly acidified with monosodium phosphate buffer (pH 5.5) and extracted with chloroform. An aliquot is injected into a GCMS for confirmation.</div> </div> <div>28.2 Specimen Requirements</div> <div> <div>28.2.1 1 mL biological fluid or comparable amount of tissue dilutions/homogenates</div> </div> <div>28.3 Reagents and Standards</div> <div> <div>28.3.1 Phensuximide</div> <div>28.3.2 Valproic acid</div> <div>28.3.3 Monosodium phosphate (NaH_2PO_4)</div> <div>28.3.4 Chloroform</div> <div>28.3.5 Methanol</div> </div> <div>28.4 Solutions, Internal Standard, Calibrators and Controls</div> <div> <div>28.4.1 1 M monosodium phosphate (pH 5.5): Weigh 13 g monosodium phosphate (NaH_2PO_4) and transfer to a 100 mL volumetric flask. QS to volume with dH_2O.</div> <div>28.4.2 20 mg/mL valproic acid stock solution. Weigh 200 mg valproic acid and transfer to 10 mL volumetric flask. QS to volume with methanol.</div> <div>28.4.3 0.2 mg/mL phensuximide internal standard solution. Weigh 10 mg phensuximide and transfer to 50 mL volumetric flask. QS to volume with methanol.</div> <div>28.4.4 Controls</div> <div> <div>28.4.4.1 Negative control. Blood bank blood (or comparable) determined not to contain valproic acid or phensuximide.</div> <div>28.4.4.2 Positive control. In house control containing valproic acid spiked at concentration similar to case specimens.</div> </div> </div> <div>28.5 Apparatus</div> <div> <div>28.5.1 Agilent GC/MSD, Chemstation software</div> <div>28.5.2 Test tubes, 16 x 114 mm (10 mL) glass centrifuge, conical bottom</div> <div>28.5.3 Centrifuge capable of 2000 – 3000 rpm</div> <div>28.5.4 Vortex mixer</div> <div>28.5.5 GC autosampler vials and inserts</div> </div> </div> </div>	

28 VALPROIC ACID CONFIRMATION BY GCMS	Page 2 of 3
<div>Division of Forensic Science</div> <div>TOXICOLOGY TECHNICAL PROCEDURES MANUAL</div>	Amendment Designator:
	Effective Date: 31-March-2004
<p>28.5.6 GC/MSD parameters. Conditions may be changed to permit improved performance.</p> <p>28.5.6.1 Acquisition Mode: Scan (50 – 550 amu)</p> <p>28.5.6.2 Column: HP 5MS 25 m x 0.25 mm x 0.25 µm</p> <p>28.5.6.3 Detector Temperature: 280° C</p> <p>28.5.6.4 Oven Program</p> <ul style="list-style-type: none"> • Equilibration time: 0.50 minutes • Initial temp: 110° C • Initial time: 1 minutes • Ramp: 10° C/min • Final Temp: 280° C • Final Time: 12 minutes • Run Time: 28 minutes <p>28.5.6.4.1 Inlet</p> <ul style="list-style-type: none"> • Mode: Splitless • Temperature: 260° C • Injection volume: 1.0 µL • Purge Time: ON at 1.0 minute <p>28.6 Procedure</p> <p>28.6.1 Label clean 16 x 114 mm screw cap tubes accordingly, negative and positive control and case sample IDs.</p> <p>28.6.2 Prepare negative and positive controls.</p> <p>28.6.3 Pipet 1 mL of each case sample into appropriately labeled tubes.</p> <p>28.6.4 Add 150 µL of 0.2 mg/mL phensuximide internal standard to each tube.</p> <p>28.6.5 Add 1 mL sodium monophosphate buffer (pH 5.5) to each tube.</p> <p>28.6.6 Add 1 mL chloroform to each tube.</p> <p>28.6.7 Vortex briefly.</p> <p>28.6.8 Centrifuge at approximately 2500 rpm for 15 minutes. Break protein plug and transfer bottom organic layer to appropriately labeled GC vials.</p> <p>28.6.9 Inject 1 µl of each sample onto GCMS.</p> <p>28.7 Calculation</p> <p>28.7.1 Quantitative valproic acid results are determined by FPIA.</p> <p>28.8 Quality Control and Reporting</p> <p>28.8.1 See Toxicology Quality Guidelines</p>	

28 VALPROIC ACID CONFIRMATION BY GCMS	Page 3 of 3
<div> <div>Division of Forensic Science</div> <div>TOXICOLOGY TECHNICAL PROCEDURES MANUAL</div> </div>	Amendment Designator:
	Effective Date: 31-March-2004
<div>28.9 References</div> <div> <div>28.9.1 I. Sunshine. Methodology for Analytical Toxicology. CRC Press, 1982.</div> </div>	