

<p style="text-align: center;">5 IMMUNOASSAY</p>	<p style="text-align: right;">Page 1 of 3</p>
<p style="text-align: center;">Division of Forensic Science</p> <p style="text-align: center;">TOXICOLOGY TRAINING MANUAL</p>	<p style="text-align: right;">Amendment No.:</p>
	<p style="text-align: right;">Effective Date: 26-March-2004</p>
<p style="text-align: center;">5 IMMUNOASSAY</p> <p>5.1 Objectives</p> <p>5.1.1 Understand and explain immunoassay.</p> <p>5.1.2 Understand the theory of commonly used immunoassay testing methods.</p> <p>5.1.3 Understand the theory and practice of fluorescence polarization immunoassay (FPIA), and specifics of operation of the Abbott AxSYM or TD_xTM instrument.</p> <p>5.1.4 Prepare tissue specimens for analysis by AxSYM or TD_x.</p> <p>5.1.5 Operate the Abbott AxSYM or TD_xTM instrument as per toxicology technical procedures manual.</p> <p>5.1.6 Interpret results by thoroughly explaining the instrument printout.</p> <p>5.1.7 Understand the quality control aspects of AxSYM or TD_x testing.</p> <p>5.2 Estimated Time: Four weeks (as ½ days)</p> <p>5.3 Methods of Instruction</p> <p>5.3.1 Lectures</p> <p style="padding-left: 40px;">5.3.1.1 Principles of immunoassay</p> <p style="padding-left: 40px;">5.3.1.2 Types of immunoassays</p> <p style="padding-left: 40px;">5.3.1.3 Components and operation of AxSYM or TD_x</p> <p style="padding-left: 40px;">5.3.1.4 Specimen preparation</p> <p style="padding-left: 40px;">5.3.1.5 Specimen analysis</p> <p style="padding-left: 40px;">5.3.1.6 Result interpretation</p> <p>5.3.2 Literature Review</p> <p style="padding-left: 40px;">5.3.2.1 Goldberger, B. & Jenkins A. 1992. <i>Testing of Abused Drugs in Urine by Immunological Techniques</i>. In Service Training, AACC.</p> <p style="padding-left: 40px;">5.3.2.2 Toxicology Technical Procedures Manual</p> <p style="padding-left: 40px;">5.3.2.3 AxSYM or TD_x Operator Manual, Abbott Laboratories Diagnostic Division.</p> <p style="padding-left: 40px;">5.3.2.4 AxSYM or TD_x System Continuing Education Seminar; Abbott Laboratories, Abbott Customer Commitment Network</p> <p style="padding-left: 40px;">5.3.2.5 Moffat, A.C., editor. <i>Clarke's Analysis of Drugs and Poisons</i>, 3rd edition. London: The Pharmaceutical Press, 2004 pp 301-312.</p> <p>5.3.3 Demonstration</p>	

<p align="center">5 IMMUNOASSAY</p>	<p align="center">Page 2 of 3</p>
<p align="center">Division of Forensic Science</p> <p align="center">TOXICOLOGY TRAINING MANUAL</p>	<p align="center">Amendment No.:</p>
	<p align="center">Effective Date: 26-March-2004</p>
<p>5.3.3.1 FPIA analysis will be observed from beginning to end and notes will be taken by the Trainee.</p> <p>5.3.4 Laboratory Exercises</p> <p>5.3.4.1 Analyze one batch of 15 blood specimens by FPIA screening for at least 2 drugs of abuse. At least 5 of the specimens will be above the cutoff concentration and at least one specimen below the cutoff.</p> <p>5.3.4.2 Perform analyses for at least 2 therapeutic drugs by using blank blood and one specimen known to be positive for acetaminophen, carbamazepine, phenytoin, salicylate or valproic acid.</p> <p>5.3.4.3 Calibrate one assay.</p> <p>5.3.4.4 Perform all daily, weekly and monthly routine maintenance.</p> <p>5.4 Evaluation</p> <p>5.4.1 Written Examination</p> <p>5.4.1.1 This will be administered as a “take home” exam.</p> <p>5.4.2 Laboratory Competency Testing</p> <p>5.4.2.1 Qualitative – a series of at least 10 previously analyzed blood specimens will be presented to the Trainee for a routine DUID panel according to the Toxicology Technical Procedures Manual. Qualitative results obtained by the Trainee must agree with previous results.</p> <p>5.4.2.2 Quantitative – 5 previously analyzed blood specimens will be presented to the Trainee for quantitation of both salicylate and acetaminophen. Results obtained by the Trainee must agree within 20% of the previous results per the Toxicology Technical Procedures Manual.</p> <p>5.4.3 Courtroom Exercise</p> <p>5.4.3.1 The Trainee must be capable of answering questions on this Module such as would be expected in a courtroom scenario.</p> <p>5.5 Examination Questions</p> <p>5.5.1 Explain the advantages and disadvantages of screening for the presence of drugs.</p> <p>5.5.2 Describe the following three different types of immunoassay: radioimmunoassay (RIA), enzyme immunoassay (EIA), and fluorescence polarization immunoassay (FPIA).</p> <p>5.5.3 Explain the following terms as they apply to FPIA: antigen, antibody, monoclonal/polyclonal antibody, fluorescence polarization, cross-reactivity, cutoff, true-positive, false-positive, sensitivity, false negative and specificity.</p> <p>5.5.4 Distinguish between homogeneous (e.g. enzyme multiplied immunoassay technique (EMIT)), and heterogeneous immunoassays (e.g. enzyme linked immunosorbent assay (ELISA)).</p> <p>5.5.5 Explain cross-reactivity stating advantages and disadvantages. Include the significance of immunoassay specificity for a specific drug vs. the specificity for a drug class.</p> <p>5.5.6 Name the chemical compound that is the primary target of the antibody in each of the FPIA assays.</p> <p>5.5.7 Explain the effect of the size and rotation rate of the tracer on polarization.</p>	

5 IMMUNOASSAY	Page 3 of 3
Division of Forensic Science TOXICOLOGY TRAINING MANUAL	Amendment No.:
	Effective Date: 26-March-2004
<p>5.5.8 Explain the relationship between net polarization and the concentration of the drug being determined.</p> <p>5.5.9 Describe the components of the AxSYM or TD_x kits and explain the purpose of each. Why are urine AxSYM or TD_x kits used for blood specimens?</p> <p>5.5.10 Explain the reasons for sample preparation (tissue dilution and homogenation , protein precipitation) prior to analysis.</p> <p>5.5.11 Describe the components of the AxSYM or TD_x instrument and explain the purpose(s) of each.</p> <p style="text-align: right;">◆ End</p>	