

10 QUANTITATION	Page 1 of 3
Division of Forensic Science TOXICOLOGY TRAINING MANUAL	Amendment No.:
	Effective Date: 26-March-2004
<p style="text-align: center;">10 QUANTITATION</p> <p>10.1 Objectives</p> <p>10.1.1 Generate accurate and precise quantitative results.</p> <p>10.1.2 Demonstrate techniques used for the quantitative determination of various drugs.</p> <p>10.1.3 Construct and apply calibration curves using internal or external standards.</p> <p>10.1.4 Understand and explain the criteria for acceptance of quantitative data.</p> <p>10.1.5 Demonstrate a working knowledge of reporting quantitative results in the manner used in the toxicology section.</p> <p>10.2 Estimated Time: One Month</p> <p>10.3 Methods of Instruction</p> <p>10.3.1 Lecture</p> <p style="padding-left: 40px;">10.3.1.1 Use of calibrators, controls and standards</p> <p style="padding-left: 40px;">10.3.1.2 Software calibration</p> <p style="padding-left: 40px;">10.3.1.3 Generating calibration curves</p> <p style="padding-left: 40px;">10.3.1.4 LOD, LOQ, ULOL, LLOL</p> <p style="padding-left: 40px;">10.3.1.5 Quality control</p> <p style="padding-left: 40px;">10.3.1.6 Statistical procedures</p> <p style="padding-left: 40px;">10.3.1.7 Reporting of quantitative results</p> <p>10.3.2 Literature Review</p> <p style="padding-left: 40px;">10.3.2.1 Toxicology Technical Procedures Manual</p> <p style="padding-left: 40px;">10.3.2.2 Agilent Technologies GCMS Instrument Manuals</p> <p style="padding-left: 40px;">10.3.2.3 HP6890 Series GC Operating Manual</p> <p style="padding-left: 40px;">10.3.2.4 Taylor, J.K. <i>Quality Assurance of Chemical Measurements</i>. 1987, Lewis Publishers, US.</p> <p style="padding-left: 40px;">10.3.2.5 Willard, H.H., Merritt, L.L.H., Dean, J. and F.A. Settle. <i>Instrumental Methods of Analysis</i>, 7th Ed., 1988, Wadsworth Pub Co, pp 540-578.</p> <p>10.3.3 Demonstration</p> <p style="padding-left: 40px;">10.3.3.1 Techniques used to quantitate drugs will be observed from beginning to end and notes will be taken by the Trainee.</p> <p>10.3.4 Laboratory Exercises</p>	

<p align="center">10 QUANTITATION</p>	<p align="center">Page 2 of 3</p>
<p align="center">Division of Forensic Science</p> <p align="center">TOXICOLOGY TRAINING MANUAL</p>	<p align="center">Amendment No.:</p>
	<p align="center">Effective Date: 26-March-2004</p>
<p>10.3.4.1 Determine the concentrations of two basic drugs in one batch of 15 blood specimens using gas chromatography. At least 10 of the specimens will be positive for at least one drug and at least one specimen will be negative.</p> <p>10.3.4.2 Determine the concentrations of benzodiazepines in one batch of 15 blood specimens using the GC/ECD. At least 10 of the specimens will be positive for at least one drug and at least one specimen will be negative.</p> <p>10.3.4.3 Determine the concentrations of cocaine, cocaethylene and benzoylecgonine in one batch of 15 blood specimens using GCMS. At least 10 of the specimens will be positive for at least one drug and at least one specimen will be negative.</p>	
<p>10.4 Evaluation</p>	
<p>10.4.1 Written Exam</p>	
<p>10.4.1.1 This will be administered as a “take home” exam.</p>	
<p>10.4.2 Laboratory Competency Testing</p>	
<p>10.4.2.1 Quantitation using GC/NPD – 5 previously analyzed cases will be presented to the trainee for quantitation.</p>	
<p>10.4.2.2 Quantitation using GC/ECD – 5 previously analyzed blood specimens will be presented to the trainee for the quantitation of benzodiazepines.</p>	
<p>10.4.2.3 Quantitation using GC/MS – 5 previously analyzed blood specimens will be presented to the trainee for the quantitation of cocaine, cocaethylene and benzoylecgonine.</p>	
<p>10.4.3 Courtroom Exercise</p>	
<p>10.4.3.1 The Trainee must be capable of answering questions on this Module such as would be expected in a courtroom scenario.</p>	
<p>10.5 Examination Questions</p>	
<p>10.5.1 Define and explain the following:</p>	
<p>10.5.1.1 Blank</p>	
<p>10.5.1.2 Internal standard</p>	
<p>10.5.1.3 External standard</p>	
<p>10.5.1.4 Control</p>	
<p>10.5.1.5 Calibrator</p>	
<p>10.5.2 Define and explain the following:</p>	
<p>10.5.2.1 Limit of detection (LOD)</p>	
<p>10.5.2.2 Limit of quantitation (LOQ)</p>	
<p>10.5.2.3 Upper limit of linearity (ULOL)</p>	

10 QUANTITATION	Page 3 of 3
Division of Forensic Science TOXICOLOGY TRAINING MANUAL	Amendment No.:
	Effective Date: 26-March-2004
<div data-bbox="297 291 766 325" data-label="Text"> <p>10.5.2.4 Lower limit of linearity (LLOL)</p> </div> <div data-bbox="211 352 906 386" data-label="Text"> <p>10.5.3 Describe the calibration curve used in the ethanol assay:</p> </div> <div data-bbox="297 413 609 447" data-label="Text"> <p>10.5.3.1 Curve generation</p> </div> <div data-bbox="297 474 630 508" data-label="Text"> <p>10.5.3.2 Acceptance criteria</p> </div> <div data-bbox="297 535 654 569" data-label="Text"> <p>10.5.3.3 Within run variability</p> </div> <div data-bbox="1461 598 1549 630" data-label="Text"> <p>◆ End</p> </div>	