

<p align="center">9 MBD</p>	<p align="center">Page 1 of 2</p>
<p align="center">Division of Forensic Science</p> <p align="center">LATENT FINGERPRINTS PROCEDURES MANUAL</p>	<p align="center">Amendment Designator:</p>
	<p align="center">Effective Date: 29-January-2004</p>
<div data-bbox="802 289 898 321" data-label="Section-Header"> <p align="center">9 MBD</p> </div> <div data-bbox="151 352 420 384" data-label="Section-Header"> <p>9.1 INTRODUCTION</p> </div> <div data-bbox="209 415 1503 506" data-label="Text"> <p>MBD is a supplemental processing procedure designed to enhance faint or indistinct impressions developed by super glue processing. MBD has an affinity for adhesion to polymerized latent impressions even at levels below visual observation. Excitation of MBD with the ALS produces extremely bright fluorescence at about 450 nm.</p> </div> <div data-bbox="151 537 420 569" data-label="Section-Header"> <p>9.2 PREPARATIONS</p> </div> <div data-bbox="209 600 722 632" data-label="Text"> <p>The examiner can mix one preparation of MBD.</p> </div> <div data-bbox="209 663 493 695" data-label="Section-Header"> <p>9.2.1 Methanol Formula</p> </div> <div data-bbox="297 726 911 758" data-label="List-Group"> <ol style="list-style-type: none"> 1. Dissolve 0.12 grams of MBD in 4.0 liters of methanol. </div> <div data-bbox="151 789 477 821" data-label="Section-Header"> <p>9.3 INSTRUMENTATION</p> </div> <div data-bbox="209 852 453 884" data-label="Section-Header"> <p>Alternate Light Source</p> </div> <div data-bbox="209 915 1481 968" data-label="Text"> <p>Alternate light sources can be used to illuminate the evidence and produce the desired fluorescence. The most common wavelength of light used is 450 nm.</p> </div> <div data-bbox="209 999 1528 1094" data-label="Text"> <p>Proper safety precautions including avoiding skin exposure and proper eye protection with appropriate optical densities should be utilized when operating ultraviolet light sources, or alternate light sources. Consult the appropriate users manuals for the safe use and appropriate eye protection for the specific piece of equipment being utilized.</p> </div> <div data-bbox="151 1125 696 1157" data-label="Section-Header"> <p>9.4 MINIMUM STANDARDS & CONTROLS</p> </div> <div data-bbox="209 1188 1528 1241" data-label="Text"> <p>Dye stains, such as MBD, work by discoloring latent impressions developed with cyanoacrylate ester. Due to their inherent ability to stain and discolor the ridge detail, there is no need for test impressions to be done prior to evidence application.</p> </div> <div data-bbox="151 1272 566 1304" data-label="Section-Header"> <p>9.5 PROCEDURE OR ANALYSIS</p> </div> <div data-bbox="209 1335 745 1367" data-label="Text"> <p>All applications should be done in a fume hood.</p> </div> <div data-bbox="209 1398 1183 1608" data-label="List-Group"> <ol style="list-style-type: none"> 1. Apply the solution to the item to be processed by immersion or squirt bottle. 2. Allow to dry. 3. Examine the item with the alternate light source at, 450 nm using the appropriate filters. 4. Have any impressions photographed. </div> <div data-bbox="151 1640 623 1671" data-label="Section-Header"> <p>9.6 INTERPRETATION OF RESULTS</p> </div> <div data-bbox="209 1703 1528 1797" data-label="Text"> <p>If the impressions are faint, repeated applications of the MBD solution may be attempted. Photographic preservation incorporating orange filters as used in the evidence examination or a Wratten #21 filter, and panchromatic films prove quite successful with even faint fluorescence.</p> </div>	

9 MBD	Page 2 of 2
<div>Division of Forensic Science</div> <div>LATENT FINGERPRINTS PROCEDURES MANUAL</div>	Amendment Designator:
	Effective Date: 29-January-2004
<div>9.7 REFERENCES</div> <div> <ol style="list-style-type: none"> 1. Lennard, Christopher J.; Pierre A. Margot. "Sequencing of Reagents for the Improved Visualization of Latent Fingerprints"; <i>Journal of Forensic Identification</i>, September/October 1988, 38, 5, 197-210. 2. Kent, Terry, ed. <i>Fingerprint Development Techniques</i>; Heanor Gate Publisher: Derbyshire, England, 1993. 3. Menzel, E. Roland. <i>Fingerprint Detection with Lasers</i>; Marcel Dekker: New York, 1980. 5. Menzel, E. Roland. "A Guide to Laser Latent Fingerprint Development Procedures"; <i>Identification News</i>, September 1983. 6. Menzel, E. Roland. "Detection Of Latent Fingerprints By Laser-excited Luminescence"; <i>Analytical Chemistry</i>, 1989, 61, 8, 557-561. <div>◆End</div> </div>	