

VIRGINIA DIVISION OF FORENSIC SCIENCE

QIAGEN BIOROBOT MONTHLY MAINTENANCE

Operator: _____

Date: _____

System liquid tubing check Flushed? YES / NO

If "no" is indicated, then no build-up was noted in the lines

Peristaltic pump tubing check

Record Inaccuracy (uL): _____

Inaccuracy corrected? YES / NO

Tubing replaced? YES / NO

If tubing replaced, re-run test

Record inaccuracy (uL): _____

Inaccuracy corrected? YES / NO

Dilutor check (mark if "OK" or "replaced")

Probe 1: _____

Probe 2: _____

Probe 3: _____

Probe 4: _____

Clean out old data DNA data file cleared: _____

Refrigerator cleared: _____

Other comments/work performed: