

<p align="center">4 BLOOD ALCOHOL</p>	<p align="center">Page 1 of 2</p>
<p align="center">Division of Forensic Science TOXICOLOGY TRAINING MANUAL</p>	<p align="center">Amendment No.:</p>
	<p align="center">Effective Date: 26-March-2004</p>
<p align="center">4 BLOOD ALCOHOL</p> <p>4.1 Objectives</p> <p>4.1.1 Understand the theory and application of headspace gas chromatography (GC).</p> <p>4.1.2 Comprehend the function and the specifics of operation of headspace GC.</p> <p>4.1.3 Prepare specimens for analysis by headspace GC.</p> <p>4.1.4 Operate the headspace GC.</p> <p>4.1.5 Calibrate the instrument and quantitate ethanol, methanol, acetone and 2-propanol.</p> <p>4.1.6 Interpret results by thoroughly examining and explaining the chromatograms.</p> <p>4.1.7 Understand the use of internal and external standards.</p> <p>4.1.8 Demonstrate proficiency by analyzing two full runs (20 samples each) of blood alcohol cases.</p> <p>4.1.9 Process results and record results of medical examiner, DUI/DUID and police casework.</p> <p>4.2 Estimated Time: One month</p> <p>4.3 Methods of Instruction</p> <p>4.3.1 Lectures</p> <p>4.3.1.1 Principles of headspace GC</p> <p>4.3.1.2 Operation of the headspace GC</p> <p>4.3.1.3 Specimen preparation (dilution, internal standard, external standard)</p> <p>4.3.1.4 Calibration and QC</p> <p>4.3.1.5 Result interpretation</p> <p>4.3.1.6 Paperwork processing in medical examiner, DUI/DUID and police casework</p> <p>4.3.2 Literature Review</p> <p>4.3.2.1 Garriott, J. C., <i>Medicolegal Aspects of Alcohol</i>. 4th^d Ed. 2003, Lawyers & Judges Pub. Co, Inc.</p> <p>4.3.2.2 Toxicology Technical Procedures Manual</p> <p>4.3.2.3 Code of Virginia (§18.2-266)</p> <p>4.3.2.4 Moffat, A.C., editor. <i>Clarke's Analysis of Drugs and Poisons</i>, 3rd edition. London: The Pharmaceutical Press, 2004 pp 53-67.</p> <p>4.3.3 Demonstration</p>	

<p align="center">4 BLOOD ALCOHOL</p>	<p align="center">Page 2 of 2</p>
<p align="center">Division of Forensic Science</p> <p align="center">TOXICOLOGY TRAINING MANUAL</p>	<p align="center">Amendment No.:</p>
	<p align="center">Effective Date: 26-March-2004</p>
<p>4.3.3.1 Blood alcohol analysis and operation of the headspace GC will be observed from beginning to end and notes will be taken by the Trainee.</p> <p>4.3.4 Laboratory Exercises</p> <p>4.3.4.1 Analyze one batch of 20 ME blood specimens for ethanol. At least 5 of the specimens will be positive for ethanol and at least one specimen will be negative.</p> <p>4.3.4.2 Analyze one batch of 20 DUI blood specimens for ethanol. At least 10 of the specimens will be positive for ethanol and at least one specimen will be negative.</p> <p>4.4 Evaluation:</p> <p>4.4.1 Written Examination</p> <p>4.4.1.1 This will be administered as a “take home” exam.</p> <p>4.4.2 Laboratory Competency Testing</p> <p>4.4.2.1 A series of at least 20 previously analyzed ME blood specimens will be presented to the Trainee for a routine blood alcohol analysis. Trainee’s results must fall within 10% of the previous results.</p> <p>4.4.2.2 A series of at least 20 previously analyzed DUI/DUID blood specimens will be presented to the Trainee for a routine blood alcohol analysis. Trainee’s results must fall within 10% of the previous results.</p> <p>4.5 Examination Questions</p> <p>4.5.1 Explain when calibration or recalibration of the headspace GC is necessary. How is recalibration accomplished?</p> <p>4.5.2 What is NIST? Why is it important?</p> <p>4.5.3 Discuss the relationship between the concentration of alcohol in blood with that in urine, serum, and vitreous humor.</p> <p>4.5.4 Explain the difference between serum and blood ethanol.</p> <p>4.5.5 Explain what causes the blood alcohol concentration in a specimen to either decrease or increase. What measures can be taken to prevent this?</p> <p>4.5.6 Explain the ethanol interconversion between mg/L, mg/dL, µg/mL and gm%. Present 5 examples of each.</p> <p>4.5.7 What is the purpose of running a mixed volatile control during the prerin? How can you calculate the acetone concentration, if present.</p> <p>4.5.8 What are the properties of a good internal standard?</p>	

◆ End