

<p align="center"><b>6 SPECTROPHOTOMETRY</b></p>	<p align="center">Page 1 of 2</p>
<p align="center"><b>Division of Forensic Science</b></p> <p align="center"><b>TOXICOLOGY TRAINING MANUAL</b></p>	<p align="center">Amendment No.:</p>
	<p align="center">Effective Date: 26-March-2004</p>
<p align="center"><b>6 SPECTROPHOTOMETRY</b></p> <p><b>6.1 Objectives</b></p> <p>6.1.1 Understand and explain the principles of ultraviolet (UV), visible (VIS), atomic absorption (AA) and fluorescence spectrophotometric measurements.</p> <p>6.1.2 Understand the practice of UV/VIS spectrophotometry and the specifics of operation of the spectrophotometers at DFS.</p> <p>6.1.3 Perform instrumental analysis of carboxyhemoglobin using a UV/VIS spectrophotometer.</p> <p>6.1.4 Interpret results by thoroughly examining and explaining the instrument printout.</p> <p>6.1.5 Understand the quality control aspects of spectrophotometric testing.</p> <p><b>6.2 Estimated Time:</b> Two weeks</p> <p><b>6.3 Methods of Instruction</b></p> <p>6.3.1 Lectures</p> <p>6.3.1.1 Principles of spectrophotometry and spectrofluorometry</p> <p>6.3.1.2 Components and operation of the UV/VIS spectrophotometer</p> <p>6.3.1.3 Specimen preparation</p> <p>6.3.1.4 Specimen analysis</p> <p>6.3.1.5 Result interpretation</p> <p>6.3.1.6 Palladium chloride diffusion confirmation test</p> <p>6.3.1.7 Salicylate confirmation by VIS spectrometry</p> <p>6.3.2 Literature Review</p> <p>6.3.2.1 DU Series 600 and 7000 Spectrophotometer: <i>Basic Instrument Training Manual and Advanced Application Training Manual</i>, Beckman</p> <p>6.3.2.2 Toxicology Technical Procedures Manual</p> <p>6.3.2.3 Moffat, A.C., editor. <i>Clarke's Analysis of Drugs and Poisons</i>, 3<sup>rd</sup> edition. London: The Pharmaceutical Press, 2004 pp 313-327.</p> <p>6.3.2.4 Willard, H. H., Merritt, L.L. Jr., Dean, J., Settle, F.A., <i>Instrumental Methods of Analysis</i>, 7<sup>th</sup> Ed. 1988, Wadsworth Pub. Co., pp. 159-223.</p> <p>6.3.3 Demonstration</p> <p>6.3.3.1 The use of UV/VIS spectrophotometry for the quantitative analyses of carbon monoxide and salicylates will be observed from beginning to end and notes will be taken by the Trainee.</p>	

<p align="center"><b>6 SPECTROPHOTOMETRY</b></p>	<p align="center">Page 2 of 2</p>
<p align="center"><b>Division of Forensic Science</b></p> <p align="center"><b>TOXICOLOGY TRAINING MANUAL</b></p>	<p align="center">Amendment No.:</p>
	<p align="center">Effective Date: 26-March-2004</p>
<p>6.3.4 Laboratory Exercises</p> <p>6.3.4.1 Analyze low, medium and high controls for the presence of carbon monoxide (CO).</p> <p>6.3.4.2 Screen one batch of 5 blood specimens for the presence of CO. At least 2 of the specimens will be positive and at least one specimen will be negative.</p> <p>6.3.4.2.1 Calculate the % saturation of each specimen.</p> <p>6.3.4.3 Confirm the presence of CO using the palladium chloride diffusion test.</p> <p><b>6.4 Evaluation</b></p> <p>6.4.1 Written Examination</p> <p>6.4.1.1 This will be administered as a “take home” exam.</p> <p>6.4.2 Laboratory Competency Testing</p> <p>6.4.2.1 A series of at least 5 previously analyzed blood specimens will be presented to the Trainee for CO analysis. The results obtained by the Trainee must agree within 20% of the previous results.</p> <p>6.4.3 Courtroom Exercise</p> <p>6.4.3.1 The Trainee must be capable of answering questions on this Module such as would be expected in a courtroom scenario.</p> <p><b>6.5 Examination Questions</b></p> <p>6.5.1 What are the wavelength ranges for visible and ultraviolet electromagnetic radiation?</p> <p>6.5.2 Explain what effects a change in solvent might have on the spectrum of a solute.</p> <p>6.5.3 Discuss why a change in the pH of a solution can be important when using UV for analysis.</p> <p>6.5.4 List and discuss some common sources of error in spectrophotometric measurements.</p> <p>6.5.5 Define the following terms: wavelength, absorbance, transmittance, excitation, emission, bandwidth and Beer's law.</p> <p>6.5.6 In the quantitative carboxyhemoglobin analysis, explain deoxyhemoglobin, oxyhemoglobin, methemoglobin and carboxyhemoglobin.</p> <p>6.5.7 How are the results reported on the certificate of analysis for CO?</p> <p>6.5.8 Explain the principle of the palladium chloride confirmation.</p> <p align="right">◆ End</p>	