

**PRIMER RESIDUE EVIDENCE HANDLING WORKSHEET FS#:** \_\_\_\_\_**Date:** \_\_\_\_\_**Examiner:** \_\_\_\_\_**Evidence Description:** One sealed \_\_\_\_\_ containing:

One sealed Primer Residue kit (item # \_\_\_\_\_) from \_\_\_\_\_

One sealed Primer Residue kit (item # \_\_\_\_\_) from \_\_\_\_\_

One sealed Primer Residue kit (item # \_\_\_\_\_) from \_\_\_\_\_

One sealed Primer Residue kit (item # \_\_\_\_\_) from \_\_\_\_\_

**Date:**

Opened sealed Primer Residue kit item/s \_\_\_\_\_, each individually. Visually examined samples. Label and initial sample vials. Make stub labels. Load SEM and run auto primer residue analysis as described in Primer Residue Worksheet SEM/EDS \_\_\_\_\_.

**Notes:** \_\_\_\_\_\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_**Date:**

Manually reviewed output from auto primer residue run

Item #	RH				LH			
	Indicative	HSPR	NPF	NNA	Indicative	HSPR	NPF	NNA

**Date:**

Opened sealed Primer Residue kit item/s \_\_\_\_\_, each individually. Visually examined samples. Label and initial sample vials. Make stub labels. Load SEM and run auto primer residue analysis as described in Primer Residue Worksheet SEM/EDS \_\_\_\_\_.

**Notes:** \_\_\_\_\_\_\_\_\_\_  
\_\_\_\_\_**Date:**

Manually reviewed output from auto primer residue run.

Item #	RH				LH			
	Indicative	HSPR	NPF	NNA	Indicative	HSPR	NPF	NNA

**Notes and Comments:** \_\_\_\_\_\_\_\_\_\_  
\_\_\_\_\_