

<p align="center"><b>9 SERIAL NUMBER RESTORATION</b></p>	<p align="center">Page 1 of 6</p>
<p align="center"><b>Division of Forensic Science</b></p> <p align="center"><b>FIREARM/TOOLMARK PROCEDURES MANUAL</b></p>	<p align="center">Amendment Designator:</p>
	<p align="center">Effective Date: 21-April-2003</p>
<p align="center"><b>9 SERIAL NUMBER RESTORATION</b></p> <p><b>9.1 Introduction</b></p> <p>Many items manufactured today have serial numbers for identification that are usually die stamped. This process compresses the material in the area immediately surrounding and a short distance below the penetration of the die. Serial numbers are removed and/or obliterated in a variety of ways and may be restored if the removal/obliteration is not taken past the previously mentioned compression zone, by using methods such as polishing, Magnaflux®, chemical or electrochemical restoration, or heat restoration.</p> <p><b>9.2 Safety Considerations</b></p> <p>Examinations performed in the Firearm and Toolmark Section are inherently hazardous. These procedures involve hazardous chemicals, firearms, ammunition, and power tools. All hazardous procedures must be performed in compliance with the DFS Safety Manual.</p> <p><b>9.3 Preparation</b></p> <p><b>NOTE: ALWAYS ADD ACID TO WATER. NEVER ADD WATER TO ACID.</b></p> <p>9.3.1 Fry's Reagent</p> <ul style="list-style-type: none"> <li>• To 90 grams of Cupric Chloride (CuCl<sub>2</sub>)</li> <li>• Add 100 milliliters of distilled water (H<sub>2</sub>O)</li> <li>• Add 120 milliliters of Hydrochloric Acid (HCl)</li> <li>• Store solution in an appropriate, sealed container that is marked with the date and initials of the preparer</li> <li>• Record in the Firearms Quality Record Book</li> </ul> <p>9.3.2 Turner's Reagent</p> <ul style="list-style-type: none"> <li>• To 2.5 grams of Cupric Chloride (CuCl<sub>2</sub>)</li> <li>• Add 40 milliliters of Hydrochloric Acid (HCl)</li> <li>• Add 25 milliliters of Ethyl Alcohol</li> <li>• Add 30 milliliters of distilled water (H<sub>2</sub>O)</li> <li>• Store solution in an appropriate, sealed container that is marked with the date and initials of the preparer</li> <li>• Record in the Firearms Quality Record Book</li> </ul> <p>9.3.3 Davis's Reagent</p> <ul style="list-style-type: none"> <li>• To 5 grams of Cupric Chloride (CuCl<sub>2</sub>)</li> <li>• Add 50 milliliters of distilled water (H<sub>2</sub>O)</li> <li>• Add 50 milliliters of Hydrochloric Acid (HCl)</li> <li>• Store solution in an appropriate, sealed container that is marked with the date and initials of the preparer</li> <li>• Record in the Firearms Quality Record Book</li> </ul> <p>9.3.4 25% Nitric Acid Solution</p> <ul style="list-style-type: none"> <li>• To 75 milliliters of distilled water (H<sub>2</sub>O)</li> <li>• Add 25 milliliters of Nitric Acid (HNO<sub>3</sub>)</li> <li>• Store solution in an appropriate, sealed container that is marked with the date and initials of the preparer</li> <li>• Record in the Firearms Quality Record Book</li> </ul>	

<p align="center"><b>9 SERIAL NUMBER RESTORATION</b></p>	<p align="center">Page 2 of 6</p>
<p align="center"><b>Division of Forensic Science</b></p> <p align="center"><b>FIREARM/TOOLMARK PROCEDURES MANUAL</b></p>	<p align="center">Amendment Designator:</p>
	<p align="center">Effective Date: 21-April-2003</p>
<p>9.3.5 Acidic Ferric Chloride Solutions</p> <ul style="list-style-type: none"> <li>• To 25 grams of Ferric Chloride (FeCl<sub>3</sub>)</li> <li>• Add 100 milliliters of distilled water (H<sub>2</sub>O)</li> <li>• Add 25 milliliters of Hydrochloric Acid (HCl)</li> <li>• Store solution in an appropriate, sealed container that is marked with the date and initials of the preparer</li> <li>• Record in the Firearms Quality Record Book</li> </ul> <p>9.3.6 Ferric Chloride Solution</p> <ul style="list-style-type: none"> <li>• To 25 grams of Ferric Chloride (FeCl<sub>3</sub>)</li> <li>• Add 100 milliliters of distilled water (H<sub>2</sub>O)</li> <li>• Store solution in an appropriate, sealed container that is marked with the date and initials of the preparer</li> <li>• Record in the Firearms Quality Record Book</li> </ul> <p>9.3.7 10% Sodium Hydroxide Solution</p> <ul style="list-style-type: none"> <li>• To 100 milliliters of distilled water (H<sub>2</sub>O)</li> <li>• <u>Slowly</u> add 10 grams of Sodium Hydroxide (NaOH)</li> <li>• Store solution in an appropriate, sealed container that is marked with the date and initials of the preparer</li> <li>• Record in the Firearms Quality Record Book</li> </ul> <p>9.3.8 Hydrofluoric Acid Solution</p> <p><b>WARNING! Concentrated Hydrofluoric Acid (HF) is a “particularly hazardous substance” and must be handled using appropriate PPE (laboratory coat, thick “rubber” gloves, and face shield). Calcium gluconate gel must be available in the work area. HF may not be handled when working alone.</b></p> <ul style="list-style-type: none"> <li>• To two (2) parts of Concentrated Hydrofluoric Acid (HF)</li> <li>• Add one (1) part of Nitric Acid (HNO<sub>3</sub>)</li> <li>• Add three (3) parts of Glycerol</li> <li>• Store solution in an appropriate, sealed container that is marked with the date and initials of the preparer</li> <li>• Record in the Firearms Quality Record Book</li> </ul> <p>9.3.9 Aqua Regia Solution</p> <ul style="list-style-type: none"> <li>• To 75 milliliters of Hydrochloric Acid (HCl)</li> <li>• Add 25 milliliters of Nitric Acid (HNO<sub>3</sub>)</li> <li>• Store solution in an appropriate, loosely sealed container that is marked with the date and initials of the preparer</li> <li>• Record in the Firearms Quality Record Book</li> </ul> <p>9.3.10 Cupric Chloride in Nitric Acid Solution</p> <ul style="list-style-type: none"> <li>• To five (5) grams of Cupric Chloride (CuCl<sub>2</sub>)</li> <li>• Add 100 milliliters of distilled water (H<sub>2</sub>O)</li> <li>• Add three (3) milliliters of Hydrochloric Acid (HCl)</li> <li>• Store solution in an appropriate, sealed container that is marked with the date and initials of the preparer</li> <li>• Record in the Firearms Quality Record Book</li> </ul>	

<p align="center"><b>9 SERIAL NUMBER RESTORATION</b></p>	<p align="center">Page 3 of 6</p>
<p align="center"><b>Division of Forensic Science</b></p> <p align="center"><b>FIREARM/TOOLMARK PROCEDURES MANUAL</b></p>	<p align="center">Amendment Designator:</p>
	<p align="center">Effective Date: 21-April-2003</p>
<p>9.3.11 Zinc Alloy Etching Solutions</p> <ul style="list-style-type: none"> <li>Solution 1 - To two (2) milliliters of distilled water (H<sub>2</sub>O), add 98 milliliters of Phosphoric Acid (H<sub>3</sub>PO<sub>4</sub>)</li> <li>Solution 2 - To 95 milliliters of distilled water (H<sub>2</sub>O), add five (5) milliliters of Nitric Acid (HNO<sub>3</sub>)</li> <li>Store solution in an appropriate, sealed container that is marked with the date and initials of the preparer</li> <li>Record in the Firearms Quality Record Book</li> </ul> <p>9.3.12 Griffin's Reagent</p> <ul style="list-style-type: none"> <li>To 30 grams of Cupric Chloride (CuCl<sub>2</sub>)</li> <li>Add 30 milliliters of distilled water (H<sub>2</sub>O)</li> <li>Add 30 milliliters of Hydrochloric Acid (HCl)</li> <li>Add 120 milliliters of Methanol (CH<sub>3</sub>OH)</li> <li>Store solution in an appropriate, sealed container that is marked with the date and initials of the preparer</li> <li>Record in the Firearms Quality Record Book</li> </ul> <p><b>9.4 Instrumentation</b></p> <ul style="list-style-type: none"> <li>Scale/Balance</li> <li>Low voltage DC power source</li> <li>UV light source (if 14AM Prepared Bath is being used)</li> <li>Yoke magnets</li> <li>Y-7 AC/DC Yoke electromagnet</li> </ul> <p><b>9.5 Minimum Analytical Standards and Controls</b></p> <p>NONE</p> <p><b>9.6 Procedure or Analysis</b></p> <p>The evidence will be marked in accordance with the Quality Manual. Initial inspection of the serial number area should include observations of coating, trace material, character remnants, and the method of obliteration. Initial observations should be recorded in the notes by documenting and/or photographing the serial number area.</p> <p>Serial numbers are removed and/or obliterated in a variety of ways and may be restored if the removal/obliteration is not taken past the previously mentioned compression zone, by using methods such as polishing, magnetic or Magnaflux®, chemical or electrochemical restoration, or heat restoration.</p> <p><b>9.6.1 Polishing Procedure</b></p> <p>The polishing procedure is a desirable method used to remove prior obliteration by polishing, grinding, and filing scratches that obscure the serial number. The polishing procedure can be effective independently, but is more often used in conjunction with various chemical or heat-restoration procedures.</p> <ul style="list-style-type: none"> <li>Polish the area of the obliteration using either a Dremel-type tool with a sanding/polishing disc or fine-grit sandpaper</li> <li>Depending on the extent of the obliteration, continue polishing until the surface is mirror-like, removing all scratches</li> <li>If the obliteration is severe, it may not be possible or desirable to remove all the scratches</li> <li>Observations should be recorded in examiner's notes by documenting and/or photographing the serial number area</li> </ul>	

<p align="center"><b>9 SERIAL NUMBER RESTORATION</b></p>	<p align="center">Page 4 of 6</p>
<p align="center"><b>Division of Forensic Science</b></p> <p align="center"><b>FIREARM/TOOLMARK PROCEDURES MANUAL</b></p>	<p align="center">Amendment Designator:</p>
	<p align="center">Effective Date: 21-April-2003</p>
<div data-bbox="248 296 743 327" data-label="Section-Header"> <p>9.6.2      Magnetic or Magnaflux® Procedure</p> </div> <div data-bbox="355 357 1546 510" data-label="Text"> <p>The magnetic procedure or Magnaflux® technique is used by metallurgists to detect surface or subsurface flaws in iron or steel. Magnetic particles, applied to a magnetized specimen, outline the obliterated characters in a successful restoration. This procedure, in conjunction with the polishing procedure, is an effective way to restore an obliterated serial number in magnetic metal. The Magnaflux® technique is nondestructive, and can be applied without affecting other restoration methods.</p> </div> <div data-bbox="391 541 1520 919" data-label="List-Group"> <ul style="list-style-type: none"> <li>• Determine the serial number medium's physical properties, i.e. magnetic or non-magnetic</li> <li>• Determine whether the specimen is suitable for testing with Magnaflux® by placing a magnet on the area of obliteration</li> <li>• The specimen is suitable if it can be magnetized</li> <li>• Clean the area of obliteration with the SKC-S Cleaner/Remover by spraying this onto the surface and wiping, allowing to dry before proceeding</li> <li>• Apply appropriately prepared 9CM or 7HF Bath to the area of obliteration with a disposable pipette</li> <li>• Place a magnet behind the area of obliteration, with the magnetic poles on either side of the area</li> <li>• This placement may be adjusted to reveal more or different areas of the obliteration</li> <li>• If 14AM (Fluorescent) prepared bath is being used, observe the characters under a black light</li> <li>• Observations should be recorded in examiner's notes by documenting and/or photographing the serial number area</li> </ul> </div> <div data-bbox="248 951 565 982" data-label="Section-Header"> <p>9.6.3      Chemical Procedure</p> </div> <div data-bbox="355 1012 1520 1075" data-label="Text"> <p>The chemical-restoration procedure is suitable for restoration of serial numbers in metal. This procedure, in conjunction with the polishing procedure, is an effective way to restore an obliterated serial number in metal.</p> </div> <div data-bbox="355 1104 1520 1167" data-label="Text"> <p>Selection of the appropriate chemical reagent, based on initial observations, may include magnetic media or non-magnet media.</p> </div> <div data-bbox="355 1197 829 1228" data-label="Text"> <p>Appropriate magnetic media reagent choices:</p> </div> <div data-bbox="391 1260 711 1419" data-label="List-Group"> <ul style="list-style-type: none"> <li>• Fry's Reagent</li> <li>• Turner's Reagent</li> <li>• Davis's Reagent</li> <li>• 25% Nitric Acid Solution</li> <li>• Aqua Regia</li> </ul> </div> <div data-bbox="345 1449 873 1480" data-label="Text"> <p>Appropriate non-magnetic media reagent choices:</p> </div> <div data-bbox="391 1512 794 1671" data-label="List-Group"> <ul style="list-style-type: none"> <li>• Ferric Chloride Solution</li> <li>• Acidic Ferric Chloride Solution</li> <li>• 25% Nitric Acid Solution</li> <li>• 10% Sodium Hydroxide Solution</li> <li>• Hydrofluoric Acid Solution</li> </ul> </div> <div data-bbox="342 1703 1546 1822" data-label="Text"> <p>As appropriate, apply the chemical solution to the area of obliteration utilizing cotton tip applicators or swabs that have been moistened with the appropriate chemical reagent or solution and note any numbers or characters that become visible. Observations should be recorded in examiner's notes by documenting and/or photographing the serial number area</p> </div>	

<p align="center"><b>9 SERIAL NUMBER RESTORATION</b></p>	<p align="center">Page 5 of 6</p>
<p align="center"><b>Division of Forensic Science</b></p> <p align="center"><b>FIREARM/TOOLMARK PROCEDURES MANUAL</b></p>	<p align="center">Amendment Designator:</p>
	<p align="center">Effective Date: 21-April-2003</p>
<p>9.6.4 Electrochemical Procedure</p> <p>The electrochemical procedure is a form of chemical restoration that is enhanced by the application of voltage that speeds the oxidation process of metal. This technique, in conjunction with the polishing procedure, is an effective way to restore an obliterated serial number in metal. Selection of the appropriate chemical reagent, based on initial observations, may include magnetic media or non-magnetic media.</p> <p>Appropriate magnetic media reagent choices:</p> <ul style="list-style-type: none"> <li>• Fry's Reagent</li> <li>• Turner's Reagent</li> <li>• Davis's Reagent</li> <li>• 25% Nitric Acid Solution</li> <li>• Aqua Regia</li> </ul> <p>Appropriate non-magnetic media reagent choices:</p> <ul style="list-style-type: none"> <li>• Ferric Chloride Solution</li> <li>• Acidic Ferric Chloride Solution</li> <li>• 25% Nitric Acid Solution</li> <li>• 10% Sodium Hydroxide Solution</li> <li>• Hydrofluoric Acid Solution</li> </ul> <p>The electrochemical procedure follows.</p> <ul style="list-style-type: none"> <li>• The electrochemical technique requires the attachment of the item to the positive terminal of a power supply via the use of metal alligator clips</li> <li>• Thoroughly soak the cotton tip of an applicator with the appropriate etching chemical solution and attach the moistened cotton tip to the negative terminal of the power supply via another metal alligator clip, being certain to do so on a moistened area at the base of the cotton tip</li> <li>• Turn on the power supply and increase the voltage gradually until the reaction appears</li> <li>• Wipe the area of obliteration with the moistened cotton tip, being careful to not touch the surface of the item with the metal alligator clips and note any numbers or characters that become visible</li> <li>• Observations should be recorded in the notes by documenting and/or photographing the serial number area</li> </ul> <p>9.6.5 Heat Procedure</p> <p>The Heat-Restoration procedure is suitable for restoration of serial numbers in plastic. The die stamping, or embossing process, is a form of "cold-working" plastic. This procedure, in conjunction with the polishing procedure, is an effective way to restore an obliterated serial number in plastic.</p> <ul style="list-style-type: none"> <li>• The heat technique requires the application heat to the area of obliteration utilizing a high intensity lamp or heat gun</li> <li>• Continue the application of heat until the plastic in the obliterated area starts to liquefy and note any numbers or characters that becomes visible</li> <li>• Observations should be recorded in the notes by documenting and/or photographing the serial number area</li> </ul> <p>9.6.6 Interpretation of Results</p> <p>Interpretation of results would include full restoration, partial restoration, or unsuccessful restoration. A full</p>	

<p align="center"><b>9 SERIAL NUMBER RESTORATION</b></p>	<p align="center">Page 6 of 6</p>
<p align="center"><b>Division of Forensic Science</b></p> <p align="center"><b>FIREARM/TOOLMARK PROCEDURES MANUAL</b></p>	<p align="center">Amendment Designator:</p>
	<p align="center">Effective Date: 21-April-2003</p>
<div data-bbox="342 296 1526 415"> <p>restoration would be a total recognition of all obliterated characters. A partial restoration would be recognition of all obliterated characters less than the total being sought. An unsuccessful restoration would be no recognition of any obliterated characters. Notes should include if the restoration procedure was full, partial, or unsuccessful.</p> </div> <div data-bbox="245 447 482 476"> <p>9.6.7 Verifications</p> </div> <div data-bbox="342 508 1503 600"> <p>A second examiner will verify all restoration results. Documentation in case notes shall be in the form of a notation, hand written by the verifying examiner, and shall delineate the obliterated characters that have been restored and/or not restored, along with the date verified and the verifier's initials.</p> </div> <div data-bbox="245 630 545 659"> <p>9.6.8 Reporting Formats</p> </div> <div data-bbox="342 690 854 720"> <p>Reporting of Serial Number Restoration Results</p> </div> <div data-bbox="342 751 513 781"> <p>Full restoration:</p> </div> <div data-bbox="438 812 1533 842"> <p>The obliterated serial number on the item __ pistol/revolver/rifle/shotgun was restored to read _____.</p> </div> <div data-bbox="342 873 540 903"> <p>Partial restoration:</p> </div> <div data-bbox="438 934 1533 999"> <p>The obliterated serial number on the item __ pistol/revolver/rifle/shotgun was partially restored to read        _ _ _ _ _</p> </div> <div data-bbox="342 1029 612 1058"> <p>Unsuccessful restoration:</p> </div> <div data-bbox="438 1089 1476 1146"> <p>Attempts to restore the obliterated serial number on the item __ pistol/revolver/rifle/shotgun were unsuccessful.</p> </div> <div data-bbox="151 1178 529 1207"> <p><b>9.7 Appropriate Appendices</b></p> </div> <div data-bbox="245 1239 518 1268"> <p>Appendix - Work Sheets</p> </div> <div data-bbox="245 1299 612 1329"> <p>Appendix - Calibration Standards</p> </div> <div data-bbox="151 1360 371 1390"> <p><b>9.8 References</b></p> </div> <div data-bbox="245 1421 1321 1451"> <p>Bureau of Alcohol, Tobacco and Firearms Laboratory. <u>Serial Number Restoration Handbook</u>. 1999.</p> </div> <div data-bbox="245 1482 1528 1543"> <p>Polk, Donald, E. and Giessen, Bill, C. "Metallurgical Aspects of Serial Number Recovery". <u>AFTE Journal</u>. Vol. 21, No. 2, p. 174.</p> </div> <div data-bbox="245 1575 1429 1604"> <p>Treptow, Richard, S. <u>Handbook of Methods for the Restoration of Obliterated Serial Numbers</u>. NASA. 1978.</p> </div> <div data-bbox="1479 1608 1549 1638"> <p align="right">♦ End</p> </div>	