

<p align="center"><b>14 GENTIAN VIOLET</b></p>	<p align="center">Page 1 of 1</p>
<p align="center"><b>Division of Forensic Science</b></p> <p align="center"><b>LATENT FINGERPRINTS PROCEDURES MANUAL</b></p>	<p align="center">Amendment Designator:</p>
	<p align="center">Effective Date: 29-January-2004</p>
<div data-bbox="711 291 990 323" data-label="Section-Header"> <h2 align="center">14 GENTIAN VIOLET</h2> </div> <div data-bbox="151 354 422 386" data-label="Section-Header"> <h3>14.1 INTRODUCTION</h3> </div> <div data-bbox="207 417 1544 596" data-label="Text"> <p>Gentian violet (crystal violet) is a sensitive stain which reacts with epithelial cells and other portions of latent print residue transferred upon surface contact. The presence of sebum appears to serve as an excellent transfer medium for sloughed epidermal cells and as a result, gentian violet is usually effective on surfaces which readily hold the deposited sebum, such as the adhesive side of tapes. The high sensitivity of gentian violet produces an immediate reaction upon skin contact, therefore, leak proof gloves are required for examinations. Accidental staining of hands is relatively harmless but usually cannot be de-stained. Disappearance of discoloration is a result of cell sloughing.</p> </div> <div data-bbox="151 617 418 648" data-label="Section-Header"> <h3>14.2 PREPARATIONS</h3> </div> <div data-bbox="207 680 880 711" data-label="Text"> <p>Gentian violet working solution- 0.1% concentration preferred.</p> </div> <div data-bbox="207 735 1520 795" data-label="Text"> <p>Higher concentrations are sometimes used, but increased amounts of gentian violet are difficult to dissolve and can create a increased background discoloration.</p> </div> <div data-bbox="207 821 946 911" data-label="List-Group"> <ol style="list-style-type: none"> <li>1. If distilled water is not available deionized water may be used.</li> <li>2. Dissolve 1.0 grams of gentian violet in one liter of distilled water.</li> </ol> </div> <div data-bbox="151 936 696 968" data-label="Section-Header"> <h3>14.3 MINIMUM STANDARDS &amp; CONTROLS</h3> </div> <div data-bbox="207 997 1520 1087" data-label="Text"> <p>Dye stains, such as Gentian Violet, work by discoloring latent impressions composed of epithelial cells and sebum. Due to their inherent ability to stain and discolor these materials, there is no need for test impressions to be done prior to evidence application.</p> </div> <div data-bbox="151 1113 565 1144" data-label="Section-Header"> <h3>14.4 PROCEDURE OR ANALYSIS</h3> </div> <div data-bbox="207 1173 1529 1467" data-label="List-Group"> <ol style="list-style-type: none"> <li>1. Immerse item to be processed in the working solution in a large tray.</li> <li>2. Allow the item to remain completely immersed for approximately 30 seconds while agitating.</li> <li>3. Remove the item from the working solution and rinse excess stain from the item by washing with a gentle flow of cold tap water.</li> <li>4. This process may be repeated until optimum contrast is reached between the impressions developed and the background.</li> <li>5. Photograph any developed impressions.</li> </ol> </div> <div data-bbox="151 1493 623 1524" data-label="Section-Header"> <h3>14.5 INTERPRETATION OF RESULTS</h3> </div> <div data-bbox="207 1551 1539 1642" data-label="Text"> <p>Developed latent impressions on light colored tapes should be photographed directly. Contrast may decrease as the substrate dries. Stained impressions which fade as the tape dries may be improved by immersing the tape in a tray of clear water and photographing the impressions while the tape is submerged.</p> </div> <div data-bbox="151 1669 389 1701" data-label="Section-Header"> <h3>14.6 REFERENCES</h3> </div> <div data-bbox="207 1728 1503 1906" data-label="List-Group"> <ol style="list-style-type: none"> <li>1. Arima, T. "Development of Latent Fingerprints on Sticky Surfaces by Dye Staining or Fluorescent Brightening"; <i>Identification News</i>, February 1981.</li> <li>2. Cowger, James F. <i>Friction Ridge Skin Comparison and Identification of Fingerprints</i>; Boca Raton: CRC Press, 1993.</li> <li>3. Kent, Terry, ed. <i>Fingerprint Development Techniques</i>. Heanor Gate Publisher: Derbyshire, England, 1993.</li> </ol> </div> <div data-bbox="1490 1906 1544 1938" data-label="Text"> <p align="right">◆End</p> </div>	