

GSR EVIDENCE HANDLING WORKSHEET**FS#:** _____**Date:** _____**Examiner:** _____**Evidence Description:** One sealed _____ containing:

One sealed GSR kit (item # _____) from _____

One sealed GSR kit (item # _____) from _____

One sealed GSR kit (item # _____) from _____

One sealed GSR kit (item # _____) from _____

Date:

Opened sealed GSR kit item/s _____, each individually. Visually examined samples. Label and initial sample vials. Make stub labels. Load SEM and run auto GSR analysis as described in GSR Worksheet SEM/EDS _____.

Notes: __________

_____**Date:**

Manually reviewed output from auto GSR run. Documented a _____ particle found on the _____ sample by photography and with a hard copy of its x-ray spectrum. Removed sample/s from SEM and return it/them to their GSR kit/s. Seal and initial the GSR kit/s. Return GSR kit/s to _____. Seal and initial the _____.

Date:

Opened sealed GSR kit item/s _____, each individually. Visually examined samples. Label and initial sample vials. Make stub labels. Load SEM and run auto GSR analysis as described in GSR Worksheet SEM/EDS _____.

Notes: __________

_____**Date:**

Manually reviewed output from auto GSR run. Documented a _____ particle found on the _____ sample by photography and with a hard copy of its x-ray spectrum. Removed sample/s from SEM and return it/them to their GSR kit/s. Seal and initial the GSR kit/s. Return GSR kit/s to _____. Seal and initial the _____.

Notes and Comments: __________

