

Case Number: _____

Scientist: _____

Date: _____

COPYRIGHT © 2011

Reagent	Amount	Lot#	Source
<u>Chelex Method</u>			
Ultra-Pure Water	1.0 ml	082410MR	DFS
Chelex ®- 100	200.0 µl	210006730	BioRad
<u>Organic Method</u>			
Extraction Buffer	400.0 µl	021411MR	DFS
Proteinase K	10.0 µl	021411MR	DFS
PCIAA (25-24-1)	400.0 µl	1007005	Ambion
η-Butanol	400.0 µl	082623	Fisher
Microcon®-30		R8HN90328	Millipore
TE Buffer	(2x) 400.0 µl	032811MR	DFS
TE Buffer	100.0 µl	032811MR	DFS

Number	Specimen Description	Amount Extracted

Starting Temp

Ending Temp

Heat Block: _____

COPY

Case Number: _____

Scientist: _____

Date: _____

COPYRIGHT © 2011

Reagent	Amount	Lot#	Source
Extraction Buffer	400.0 µl	021411MR	DFS
Proteinase K	10.0 µl	021411MR	DFS
PCIAA (25-24-1)	400.0 µl	1007005	Ambion
n-Butanol	400.0 µl	082623	Fisher
Microcon® -30		R8HN90328	Millipore
TE Buffer	(2x) 400.0 µl	032811MR	DFS
TE Buffer	100.0 µl	032811MR	DFS

VIRGINIA
DEPARTMENT
OF
FORENSIC SCIENCE

Number	Specimen Description	Amount Extracted

Starting Temp

Ending Temp

Heat Block: _____

UNCONTROLLED
COPY

Case Number: _____

Scientist: _____

Date: _____

COPYRIGHT © 2011

Reagent	Amount	Lot#	Source
5% Terg-a-zyme	1.0 ml	G1E1	Alconox
Ethanol	1.0 ml	RG0232	AAPER
Ultra-Pure Water	1.0 ml	082410MR	DFS
Extraction Buffer	130.0 / 57.0 µl	021411MR	DFS
Dithiothreitol (1M)	8.0 µl	021411MR	DFS
Proteinase K	5.0 µl	021411MR	DFS
PCIAA (25-24-1)	200.0 µl	1007005	Ambion
n-Butanol	200.0 µl	082623	Fisher
Microcon® -30		R8HN90328	Millipore
TE Buffer	300 / 100 µl	032811MR	DFS

DEPARTMENT OF FORENSIC SCIENCE

Number	Specimen Description	Amount Extracted

UNCONTROLLED

Starting Temp

Ending Temp

Heat Block: _____

COPY

Case Number: _____

Scientist: _____

Date: _____

COPYRIGHT © 2011

Reagent	Amount	Lot#	Source
Ultra-Pure Water	(2x) 25 ml	082410MR	DFS
Ethanol	(2x) Cover	RG0232	AAPER
Demineralization Buffer	3.0 ml	022311MR	DFS
Proteinase K	200.0 µl	021411MR	DFS
PCIAA (25-24-1)	(2x) 3.0 ml	1007005	Ambion
n-Butanol	3.0 ml	082623	Fisher
AMICON® ULTRA-4		R7MN97183	Millipore
TE Buffer	(2x) 2.0 ml	032811MR	DFS
TE Buffer	100.0 µl	032811MR	DFS

FORENSIC SCIENCE

Number	Specimen Description	Amount Extracted

UNCONTROLLED
COPY

ANALYST: _____

FS LAB#: _____

DATE: _____

Well #	SAMPLE	Primer	+/- Product	Concentration (ng/μL)	Photo
ORIGIN 1					
1	DNA MW Marker XIV (1:2)				
2					
3					
4					
5					
6					
7					
8					
9					
10					
ORIGIN 2					
11	DNA MW Marker XIV (1:2)				
12					
13					
14					
15					
16					
17					
18					
19					
20					

PRODUCT GEL PHOTO

REAGENT	LOT#	SOURCE
PRODUCT GEL		
AGAROSE – NuSieve® 3-1	0000229335	Lonza
TAE Gel Buffer (1X)	822058	Invitrogen
SYBR Green	880315	Invitrogen
Loading Buffer (5X)	011311MMB	DFS
TAE Tank Buffer (1X)	822058	Invitrogen
DNA MW Marker XIV	12264323	Roche

Comments:

ANALYST: _____
DATE: _____

FS LAB#: _____

Well #	SAMPLE	Primer	+/- Product	Concentration (ng/μL)	Photo
ORIGIN 1					
1	DNA MW Marker XIV (1:2)				<p>PRODUCT GEL PHOTO</p>
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
ORIGIN 2					
15	DNA MW Marker XIV (1:2)				<p>Comments:</p>
16					
17					
18					
19					
20					
21					
22					
23					
24					
25					
26					
27					
28					

REAGENT	LOT#	SOURCE
PRODUCT GEL		
AGAROSE – NuSieve® 3-1	0000229335	Lonza
TAE Gel Buffer (1X)	822058	Invitrogen
SYBR Green	880315	Invitrogen
Loading Buffer (5X)	011311MMB	DFS
TAE Tank Buffer (1X)	822058	Invitrogen
DNA MW Marker XIV	12264323	Roche

Case Number: _____

Scientist: _____

Purification/Plate Loading Date: _____ / _____

Plate Name: _____

COPYRIGHT © 2011

VIRGINIA
DEPARTMENT
OF
FORENSIC SCIENCE

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Purification Method

Reagent	Amount Vol(µl) x # = Total Vol(µl)	Lot #	Source
Running Buffer (1X) POP-6		1012159 1106064	AB AB

UNCONTROLLED
COPY

Case Number: _____

Scientist: _____

Date: _____

Denature Samples		Lot #	Source	Time: Immediate		
Denaturation Solution		N02865	Roche	15µl per sample		
PCR Product				15µl each well		
Hybridize Samples		Lot #	Source	Time: 15 minutes		
Sample				30µl total		
Wash Buffer		051011MR	Roche	3.0 ml each well		
Hybridize Rinse		Lot #	Source	Time: 5-10 seconds		
Wash Buffer		051011MR	Roche	3.0 ml each well		
Enzyme Conjugate		Lot #	Source	Time: 5 minutes	Amount	
Wash Buffer		051011MR	Roche	3.0 ml each well	#	vol
Streptavidin-POD		N06219	Roche			Vol
					0.0 x 3.0 ml	3ml
					0.0 x 8.0 µl	8µl

Enzyme Conjugate Rinse		Lot #	Source	Time: 5-10 seconds		
Wash Buffer		051011MR	Roche	3.0 ml each well		
Stringency Wash		Lot #	Source	Time: 12 minutes		
Wash Buffer		051011MR	Roche	3.0 ml each well		
Stringency Rinse		Lot #	Source	Time: 5-10 seconds		
Wash Buffer		051011MR	Roche	3.0 ml each well		
Color Development Rinse		Lot #	Source	Time: 5 minutes		
Citrate Buffer		050211MR	DFS	3.0 ml each well		
Color Development		Lot #	Source	Time: 15 minutes	Amount	
Citrate Buffer		050211MR	DFS	3.0 ml each well	#	vol
Hydrogen Peroxide (30%)		097162	Fisher			Vol
Chromogen-TMB		021711MR	DFS			
					0.0 x 3.0 ml	3ml
					0.0 x 0.4 µl	0.4µl
					0.0 x 0.15 ml	0.15ml

Color Development End				Repeat 2X
Water				Temperature: Room
				5.0 ml each well

