

		TECAN						Plate Washer			Reader		ALL	
		Daily			Weekly		Monthly	Daily	Monthly		Daily	Monthly	Annually	
Date	Initials	Clean surfaces w/ isopropanol	Wipe probe w/ isopropanol	Check fluid levels / Empty waste containers	Check tubing for leaks	Clean liquid container w/ mild soap solution	Acid/Base Wash	Clean surfaces w/ isopropanol	Disassemble Manifold. Sonicate in MEOH solution	Clean dispensing needles w/cleaning needles	Clean surfaces w/ isopropanol	Clean filter slide w/kimwipe	Manufacturer PM	comments