

PRIMER RESIDUE EVIDENCE HANDLING WORKSHEET FS#: _____

Date: _____ **Examiner:** _____

Evidence Description: One sealed _____ containing:
 One sealed Primer Residue kit (item # _____) from _____
 One sealed Primer Residue kit (item # _____) from _____
 One sealed Primer Residue kit (item # _____) from _____
 One sealed Primer Residue kit (item # _____) from _____

Date: _____
 Opened sealed Primer Residue kit item/s _____, each individually. Visually examined samples. Label and initial sample vials. Make stub labels. Load SEM and run auto primer residue analysis as described in Primer Residue Worksheet SEM/EDS _____.

Notes: _____

Date: _____
 Manually reviewed output from auto primer residue run:

Item #	RH				LH			
	Consistent	CHAR	NPF	NNA	Consistent	CHAR	NPF	NNA

Date: _____
 Opened sealed Primer Residue kit item/s _____, each individually. Visually examined samples. Label and initial sample vials. Make stub labels. Load SEM and run auto primer residue analysis as described in Primer Residue Worksheet SEM/EDS _____.

Notes: _____

Date: _____
 Manually reviewed output from auto primer residue run:

Item #	RH				LH			
	Consistent	CHAR	NPF	NNA	Consistent	CHAR	NPF	NNA

Notes and Comments: _____

Codes: CHAR: Characteristic of primer residue NPF: No primer residue particles found NNA: Not necessary for analysis