

ABC ALCOHOL CONTENT BATCH SUMMARY WORKSHEET

Specimens removed from TX ADM storage,
 aliquots removed for analysis and
 specimens returned to TX ADM storage by: _____ (analyst) _____ (date)

Vial Check: _____ Reviewer/Date: _____ QC File: _____

Control Chart Updated (Initials/Date) : _____

POS	SAMPLE ID#	ETOH	AVE/Range	POS	SAMPLE ID#	ETOH	AVE/Range
1	CHECK (50% STD)			26			
2	50% STD		47.50-52.50%	27	"		
3	NEG			28	NIST 7.5%		7.125-7.875%
4	25% STD		23.75-26.25%	29			
5	5% STD		4.75-5.25%	30	"		
6	0.5% STD		0.475-0.525%	31			
7				32	"		
8	"			33			
9				34	"		
10	"			35			
11				36	"		
12	"			37			
13				38	"		
14	"			39	NIST 31.25%		29.69-32.81%
15				40			
16	"			41	"		
17	WQA			42			
18				43	"		
19	"			44			
20				45	"		
21	"			46			
22				47	"		
23	"			48			
24				49	"		
25	"			50	5% QC		4.75-5.25%
Calibrators:				Controls:			
50%		5%		WQA		NIST 31.25%	
25%		0.50%		NIST 7.5%		5% QC	