

FLUORESCENCE WORKSHEET

EXAMINER: _____

DATE: _____

FS#: _____

OLYMPUS BX-40 PLM _____

MAGNIFICATION: _____

OLYMPUS BX-40 COMPARISON PLM: DFS # 998 _____ **or DFS # 1212** _____

MOUNTED IN: _____

Item(s) / Layer(s)	WU Wide Ultra-Violet (330-385 nm)	WBV Wide Blue-Violet (400-440 nm)	WB Wide-Blue (450-480 nm)	WG Wide-Green (510-550 nm)	Additional Comments
QC Sample: (+) pink fibers (-)	(+): orange (-): dark	(+): bright orange (-): dark	(+): bright orange (-): dark	(+: very bright orange-red (-): dark	Mounted in: Protexx; Magnification: QC check passed ? (Yes = v, No = X) (+) <input type="checkbox"/> (-) <input type="checkbox"/>